

HELPING PRACTICES GET **BEYOND 20/20 SINCE 1954!**







BERNELL OFFICIAL UK & EUROPEAN DISTRIBUTOR

ST£75 MIN ORDER (EX VAT & DELIVERY) NO VAT IS CHARGED ON EXPORT ORDERS

For all enquiries: sales@louisstone.co.uk

Ғ @louisstoneoptical 🔰 @LouisStoneOptic 🧿 @louisstoneoptic



Lightweight, premium frame designed to fit comfortably on both children and adults. Frames work with both our flat and corrected curve 48mm prisms. Prisms snap securely in place, yet can be easily rotated by hand to the desired location. No more thumb screws!

SEE PAGE 138 FOR DETAILS.

Louis Stone Optical Ltd are proud to have been distributing Bernell Vision Training VT products for over 15 years

Bernell are one of the primary suppliers of VT products worldwide and are currently celebrating over 25 years in business.

Louis Stone are receiving more demand for these products year on year due to the demand and interest amongst optical professionals for behavioural training items and tests to assist their patients.

Browse through the online catalogue which has recently been updated with many new products, and see hundreds of products including field and macular degeneration tests, red/ green goggles, prism related items including prism goggles, VT starter and patients kits, the list is considerable.

Ordering Bernell Vision Training Products from Louis Stone couldn't be easier

All enquiries can be made by email to sales@louisstone.co.uk.

Always include name of practice, full address, contact name and tel number. If you have a VAT and EORI NUMBER also insert. This information will greatly assist us when quoting by return.

Always list model number and page number with your enquiry; if you are unsure of item required refer to product index on last 3 online pages.

On receipt of order a proforma invoice will be sent for pre-payment by credit card or bank transfer.

ST£75 Minimum order (ex VAT & delivery)

No vat is charged on export orders, however you will be responsible for local taxes and fees in country of receipt.

Louis Stone prices these products as competitively as possible and Bernell are fully aware of our pricing. By purchasing from us you will avoid USA taxes /fees / courier, admin charges, exchange rate fluctuations and paperwork to fill in.

Products are not sold on a sale or return basis and can only be returned within 7 days of receipt if there is a workmanship or product fault only. Prior to a return you must contact in the first instance rma@louisstone.co.uk for return instructions.

Delivery can take up to 21 days, however updates will be given on your order's progress.

5

Email for all Bernell enquiries on sales@louisstone.co.uk



CATALOG CONTENTS

EXAM ROOM NEEDS -- PAGES 1-24

Pg. 1-2Pharmace	euticals Pg. 15-16	Occluders
Pg. 3 Post My	driatics Pg. 17-18	Eve Patches
Pg. 4 Bulbs & B	atteries Pg. 19	Fixation Devices
Pg. 5 Dry Eye So		
Pg. 6-9 Surgical & S	5	
Pg. 10 Miscell		
Pg. 11-14 Retinoscopy/Diagnostic	Lenses Pg. 24	Measurement Devices

ASSESSMENTS -- PAGES 25 - 50

Pg. 27 Pg. 28 Pg. 29-32	Vision Screening Field & Macular Degeneration Tests Acuity - Contrast Sensitivity/OKN Acuity - Distance	Pg. 43-44 Pg. 45-46	Depth Perception Color Vision Testing
0	Acuity - Distance Acuity - Near	1 g. 47-50	

OPTICAL -- PAGES 51 - 74

Pg. 51-54	Supplies & Equipment	Pg. 59-60	Readers & Sunglasses
Pg. 55-56	Optical Accessories	Pg. 61-73	Flex Frames
Pg. 57-58	Contact Lens Products	Pa. 74	Erin's World Frames

VISION THERAPY -- PAGES 75 - 134

Pg. 75-76	VT Starter Kits	Pg. 98-104	Anti-Suppression
Pg. 77-78	VT Patient Kits	Pg. 105-109	Vectograms, Vectographs, & Tranaglyphs™
Pg. 79-83	HTS Software	Pg. 110-112	Flippers
Pg. 84	Optics Trainer Software	Pg. 113-114	Accommodation
Pg. 85-86	RightEye Vision System	Pg. 115	Haidinger Brush Equipment
Pg. 87-88	Senaptec Sports Vision Equipment	Pg. 116-125	Binocular Training
Pg. 89	BTrackS™ Balance Tracking Systems	Pg. 126-128	Perceptual
Pg. 90-91	Balance/Reaction/Sports Vision	Pg. 129	Miscellaneous
Pg. 92		Pg. 130	Posture & Writing Aids
Pg. 93-94	Rotation Trainers	Pg. 131-132	
Pg. 95-96	Marsden & Saccadic Training	Pg. 133-134	Syntonics
Pg. 97	Brock Strings	-	-

NEURO REHABILITATION -- PAGES 135 - 139

Pg. 135-138 Spatial Localization & Prism Training	Pg. 139	Prism Training & Equipment
---	---------	----------------------------

Low Vision -- Pages 140 -145

Pg. 140 Practice Aids	Pg. 143Handheld Magnifiers
Pg. 141 Electronic Magnifiers	Pg. 144Stand Magnifiers
Pg. 142Page Magnifiers	Pg. 145 Spectacle Magnifiers


REFERENCE MATERIALS -- PAGES 146 - 151

Pg. 146-147Eye Charts &	Vodels Pg. 148-151
-------------------------	--------------------



Louis Stone
OPTICAL EST. 1950BERNELL OFFICIAL
UK & EUROPEAN D **UK & EUROPEAN DISTRIBUTOR**

STE75 MIN ORDER (EX VAT & DELIVERY) NO VAT IS CHARGED ON EXPORT ORDERS



..... Books

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY



BERNELL offers a wide variety of pharmaceuticals to save our customers both time & money!

As our customers are painfully aware, there has been much turmoil these past few years due to FDA changes involving ingredients of pharmaceuticals. We have seen removal from the market of formulations which had been available for years. This causes huge market price fluctuations as well as long backorders and sometimes permanent removal. We've built-in procedures to keep pricing and availability posted online as our yearly catalog can not keep up with weekly fluctuations. Bernell always tries to offer the lowest current prices in the market!

- Note: If a specific brand is required, please specify when placing your order.
- · Pharmaceuticals from this page will arrive separately from the rest of your order and are not available outside the US.

Important Note:

A shipping surcharge of \$6.95 will be added to pharmaceutical orders from this page. However, if orders contain over \$35 in pharmaceuticals from this page, then this fee will be waived!

💥 This symbol indicates that the item requires refrigeration, which may require express shipping per FDA regulations

PART #	DESCRIPTION	DRUG TYPE	SIZE	COMPARE TO
AK49D5DS	Fluorescein Sodium 0.25% with Benoxinate 0.4% Solution	Staining Anesthetic	5mL	Fluress
AK2D15DS	Proparacaine HCI 0.5% Ophthalmic Solution	Anesthetic	15mL	Alcaine, Ophthaine, or Ophthetic
AK102D15DS	Tetracaine HCI 0.5% Ophthalmic Solution	Anesthetic	15mL	
AK42D5DS	Gentamicin Sulfate 0.3% Ophthalmic Solution	Antibiotic	5mL	
AK43D10DS	Neomycin-Polymyxin B-Gramicidin Ophthalmic Solution	Antibiotic	10mL	
AK45O3DS	Erythromycin 0.5% Ophthalmic Ointment (Preservative Free)	Antibiotic	3.5gm	llotycin®
AK40D5DS	Neomycin-Polymixin B-Dexamethasone Ophthalmic Suspension	Antibiotic/Anti-Inflammatory	5mL	Maxitrol
AK40ODS	Neomycin-Polymixin B-Dexamethasone Ointment	Antibiotic/Anti-Inflammatory	3.5gm	Maxitrol
AK59D5DS	Prednisolone Acetate 1% Ophthalmic Suspension	Anti-Inflammatory/Steroidal	5mL	
AK60D5DS	Ketorolac Tromethamine 0.5% Opthalmic Solution	Anti-Inflammatory/ Non-Steroidal	5mL	Acular
AK44D5DS	Ciprofloxacin 0.3% Ophthalmic Solution	Antibiotic/Antibacterial	5mL	Ciloxan
AK3D15DS	Cyclopentolate Hydrochloride 1% Ophthalmic Solution	Cycloplegic	15mL	Cyclogyl
AK1D15DS	Tropicamide 0.5% Ophthalmic Solution	Cycloplegic	15mL	Mydriacyl
AK1AD15DS	Tropicamide 1% Ophthalmic Solution	Cycloplegic	15mL	Mydriacyl
AK55DS	Atropine Sulfate 1% Ophthalmic Solution	Cycloplegic	5mL	
AK57D15DS	Phenylephrine Hydrochloride 2.5% Ophthalmic Solution (Paragon)	Mydriatic	15mL	
AK57NRD15DS	Phenylephrine Hydrochloride 2.5% Ophthalmic Solution (AKORN)	Mydriatic	15mL	

Visit us online at www.bernell.com to view additional pharmaceuticals that we offer, as well as access the most up-to-date pricing and availability!

NEURO REHABILITATION

Oraganize your pharmaceuticals with our doctordesigned tray. It has compartments for dilation drops, Shrimer and Fluorescein Strips, Dacriose or CL solutions, ointments, gonio, Alger Brush, and RX or Amsler Grid Pad. The tray can sit on

your desk or be wall mounted. Made from high-

quality materials and built to last.

AP100





PHARMACEUTICALS



PART #	NAME	DESCRIPTION	SIZE	PRICE
AK12B100	BIO-GLO [™] FLUORESCEIN	Used to show defects or abrasions in thecorneal epithelium.	Box of 100	Call or Go Online for Pricing & Availability
AK12B300	STRIPS	Compares to Fluorets $^{\ensuremath{\mathbb{R}}}$. 1mg Fluorescein Sodium Ophthalmic Strip.	Box of 300	Call or Go Online for Pricing & Availability
AK8810	FUL-GLO® Compare to Fluorets [®] . 0.6mg Fluorescein Sodium individually wrapped, one sterile strip per pack.		Box of 300	Call or Go Online for Pricing & Availability
AK32100	FLUORESCEIN	0.6mg Fluorescein Sodium ophthalmic strips to stain the anterior segment of the eye when fitting contact lenses, disclosing corneal injury and in applana-	Box of 100	Call or Go Online for Pricing & Availability
AK32300	GLO-STRIPS [®]	tion tonometry. Narrow, rounded tip for exact placement. Tip is smooth, lint free. Protective barrier keeps fluid off finger during application. Sterile.	Box of 300	Call or Go Online for Pricing & Availability
AK64B100	Rose Bengal GloStrips®	Diagnostic agent in routine ocular exams or when superficial corneal or conjuctiva tissue change is suspected. Come as 100 individually wrapped strips with each containing 1.3mg of Rose Bengal. Compare to Rose Glo.	Box of 100	Call or Go Online for Pricing & Availability

PART #	NAME	DESCRIPTION	SIZE	PRICE
GDROP1		A revolutionary new way to use eye drops. GentleDrop is a silicone sleeve that rests on the nasal bridge and safely delivers eye drops. It is compatible with nearly all 5, 10, & 15 ml cylindrical bottles.	Universal	
ILEGD	AUTODROP [®]	Essential for patients who struggle to put drops in their eyes. This unit quickly fits most drop bottles & keeps the tip sterile! Unique anti-blink design helps to ensure drops are successfully delivered on the 1st attempt!	Universal	
TFTS100	TEAR FLOW TEST STRIPS 100 STRIPS	The Tear Flow Test Strips feature a printed metric scale for unmistakable test interpretation. For use as a diagnostic agent when reduced tear secretions are suspected. Each sterile, lint-free strip features a comfortable rounded tip. 100 individually wrapped sterile strips. Made in the U.S.A.	Box of 100	
AK82W	Ammonia Wipes Inhalant	To prevent or treat fainting, sometimes known as smelling salts, provide immediate relief of lightheadedness. Box of 10	Box of 10	
AK51D4BL	ADVANCED EYE RELIEF EYE WASH BY BAUSCH & LOMB	Gently cleanses and refreshes the eye to help relieve irritation, stinging and discomfort by removing loose foreign material, air pollutants (pollen, ragweed), dust or chlorinated water. This sterile solution is great to have on hand for daily or emergency eye cleansing. The convenient bottle tip permits use with or without the enclosed sterile eye cup.	4 fl oz	
PH9958	P H P APER	A rainbow of distinct colors indicates each pH unit over the entire pH scale. Single-Roll dispenser, 15ft roll. Every office should have pH paper when irrigating an eye in an emergency. It is important to know that the levels of the eye are neutral, especially not too high a base!	15ft roll	

Dozens of additional pharmaceutical products are listed on our website. Buy online at www.bernell.com to save time and money when ordering supplies!

QUICKVUE ADENOVIRAL CONJUNCTIVITIS TEST

DETERMINE VIRAL OR ALLERGY IN JUST 10 MINUTES!

Acute conjunctivitis - more commonly known as "pink eye" - affects nearly 6.9 million people in the United States every year. Approximately 1 in 4 patients with acute conjunctivitis have confirmed adenoviral conjunctivitis. QuickVue Adenoviral conjunctivitis Test is the first and only CLIA-waived, point-of-care test that detects all known serotypes of adenoviral conjunctivitis.

Clinical benefits

- Identify, isolate and manage contagious patients, preventing spread to other patients and staff.
- Reduce inappropriate antibiotic use.Offer an accurate, fast, disposable
- platform that is easy to use.
- Shelf Life: 30 Months
- Storage Conditions: 59°F to 86°F

Easy and Convenient

- Disposable test, requires just 4 easy steps.
 Just sample, dip and read!
- Just sample, dip and read!
 No additional equipment is
- needed to administer or interpret.
- Fast and Accurate
- Results in 10 minutes.
- 90% sensitivity and 96% specificity.

We wanted to sell individual tests, but the FDA does not allow since the test was not labeled that way. Relabeling a test for the FDA is not an easy or inexpensive process. We know of several offices who buy a case together and share the tests. That is not illegal by the FDA.


A Certificate of Waiver is required to perform this test in a CLIA-waived setting. To obtain a Certificate of Waiver, please contact your state health department

2



Order Now by Phone at (574) 259-2070 or Go To www.bernell.com





ROLLENSA Amber ROLLENS Gray Imprinting Available on Rollens - Contact us for details!

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

5 - 9

4 or less

SAVE WITH BERN			R BEAT ANY	ADVERTISED	PRICE! -	Bulbs &	BATT	ERIES	EXAM ROOM NEEDS
PROJECT	<mark>or Bulbs</mark>								MO
Nаме	DESCRIPTION	PART #	PRICE		-			N D	õ
(A) Eiko Projector Bulb	This is our most popular bulb! (130V, 50W)	LTCAX				A	B	C	XAM F
(B) Halogen Bulb for Projector	(130V, 50W) I	HALCAX		Ei Tree	KO. E Quality Lampi al + Scientific + Dubthalinir	Ê	in the second se	Î	
(C) B & L Projector Bul	b For use with B/L (120V 100W)	GECDD		wenn	ur • Eningelitin • Epolitemer		-	Ť	ITS
Slit Lam	p Bulbs								ASSESSMENTS
NAME	Descri	PTION	PART #		D	TT	E)	F	ŝ
(D) A/O Slit Lamp Bul	-		GE1468		1		- (1		SSI
(E) Slit Lamp Bulb	For Mentor, To	,	GE64138		B	- E		R	Š
(F) Haag-Streit Slit Lamp Bulb 366	For Haag-Streit 90 G10, GH12, S10, #1000, 800, 850. \	SH12, Burton	AM366	4		-	e	P	Ļ
	#1000, 000, 000, 000.	7011290: 0774.0A		-		T			S
Ophthali	Moscope (or Tra	NSILLUN	AINATOR)	Bulbs				OPTICAL
NAM	IE	DE	SCRIPTION	Part #	PRICE		3		
(G) Ophthalmoscope E	Bulb for Welch Allyn	For 3.5V Welc	h Allyn Halogen m	odel BD03000	0				≻
(H) Transilluminator Bu	ulb for Welch Allyn	For 3.5V Welc	h Allyn Halogen m	odel BD03100			-		٩
(I) Direct Ophthalmosc for Welch Allyn	ope Bulb		1720,11730,11735 : 3.5V/.720A Halog			HG			VISION THERAPY
Retinoso	OPE BUL	38							/ISIO
Nаме	Des		P	ART # PRICE		K	M	and then	
(J) Retinoscope Bulb	For 2.5V Welch Allyn Stre	eak model 17700 (n	on-halogen) BC	01800	A.S.			BEST	Z
(K) Retinoscope Bulb	For 2.5V Welch Allyn S	treak model 18100) (halogen) BC	03700	A		11 .		Õ
(L) Retinoscope Bulb	For Welch Allyn Spot -	23mm	В	T45A				1 ell	AT
(M) Retinoscope Bulb	For Bernell Spot Retine	oscope	BT4	4500BU					Ę
(N) Retinoscope Bulb	Machine/Model #: Stre Voltage/Watts: 3.5V H		, 18245 h: 20mm	IARB	J	₩ L		N	HABII
BATTERIES - FOR OPHTHALMOSCOPES, RETINOSCOPES, & TRANSILLUMINATORS									
A CARDINAL CJB7201	B	JINAL CIE				SUPREMERO	WER .	Reitan Beitan	NEURO REHABILITATIO

BRAND	DESCRIPTION	Part #	PRICE
(A) Cardinal Battery	Fits 2.5V Welch Allyn Models: 60700,70700/10/15/20/50/51/54/55, & 71500	CM720	
(B) Cardinal Battery	Fits 3.5V Welch Allyn Models: 71000; 71010; 71015/20/50/51/54/55 & 71670, 72200) (Fits All Proper Ophthalmoscopes)	CM722	
(C) Cardinal Battery	3.5V rechargeable battery for Welch-Allyn Handles	CM72300	
(D) Cardinal Battery	Fits 2.5V Optec 360 Non-Halogen Retinoscope	CM360	
(E) Cardinal Battery	Fits 2.5V Optec 360 Halogen Retinoscope	CM360H	
(F) Supreme Power Battery	3.5V 2.75" (70mm) Rechargeable Battery. Fits Ophthalmoscope 71000/10/15/20/50/51/54/55 and 71670	BA722	
(G) Supreme Power Battery	3.5V 2.70" (68mm) Welch Allyn 72300 Convertible Handle Ophthalmoscope 71050-C Audiometer Battery	BA723	

RDINAL

D

C18360

E

Ξ

CJB360H

LOW VISION

SUPREME POWER

G (

DRY EYE SOLUTIONS



120 Soft Gels - MRSP \$49.95

60 Soft Gels - MRSP \$49.95

Save 20% Off Retail

with Lemon Flavor

with FloraGLO

NNPO120

ProDHA Eye®

NNPD60L

Nordic Naturals[®] High Quality Omega-3

Help ensure that your patients use top quality omega-3 supplements to support their eye health! Over 10 million Americans suffer from Dry Eye Syndrome (DES). A large body of scientific evidence indicates that supplementing with fish oil benefits eye health in people of all ages. Why? Because fish oil is the best source for EPA and DHA, the omega-3 essential fatty acids that are among the most important nutrients for the development, function, and maintenance of eye tissue throughout life. DHA is particularly important for the eyes, and attains its highest concentration within eye tissue, where it accounts for approximately 30% of the total fatty acids.

Why Nordic Naturals?

- Surpasses International Pharmaceutical Standards are the vitamins you recommend tested?
- Exceptional Freshness Levels oxygen-free manufacturing process maintains product integrity
- Great Taste Utilizes proprietary enzymatic & oxygen-free processing technology to ensure no fishy smell, taste, or aftertaste.
- True Triglyceride Form produced in true triglyceride form to ensure optimal absorption and results
- Research-Based Innovation - as a research-driven company, where processing technology and product formulations are driven by the latest scientific developments

ProOmega® - Concentrated fish oil formula gives you more omega-3s in fewer soft gels. Supports tear production, protects eyes from oxidative stress and damage, and provides concentrated EPA and DHA for healthy ocular tissue. Contains 2200mg EPA+DHA per 4 soft gels. Enfused with a great lemon taste!

ProDHA Eye® - Concentrated DHA fish oil with FloraGLO® Lutein for additional eye support. Supports healthy tissue moisture and tear production and provides natural pigment to ocular cells. 788mg DHA, 190mg EPA per 2 soft gels. 20mg Lutein and 4mg Zeaxanthin in each serving.

Visit www.bernell.com to learn more information about these professional grade supplements. Help your patients, while also helping your bottom line!



BLEPHADEX[™] - SCIENCE BASED EYELID CARE Unlike many of the store bought lid wipes that only offer a cleanser, Blephadex™ offers a combination of Tea Tree Oil, Coconut Oil and a gentle lid cleanser. Blephadex™ can effectively offer relief while removing excessive oils and debris that may cause many symptoms associated with Blepharitis and Demodex. Available options for In-Office and At-Home treatments include: Blephadex[™] Eyelid Foam Cleanser Blephadex[™] Eyelid Applicators hade Renh (For In-Office Use): (For At-Home Use): 0 Offers a simple and easy in-office solution for The convenience of a bottle of foaming BigBiephadeipha debridement of Blepharitis and Demodex. Comes Blephadex cleanser. Foam comes as a one-month supply. in a box of 20 individually packaged applicators for ease of application and treatment of patients. **BPXFM1** (Foam Cleanser) **BPXAPP** (Applicators) Blephadex[™] Eyelid Wipes and Warming Wipes Blephadex (For At-Home Use): **FLID WIPES** Wipes come in a box of 30 individually packaged wipes. Blephadex Warming Eyelid Wipes bring you the same effective treatment and relief that you trust along with warming comfort to refresh your eyes. BPXWI30W (Warming Eyelid Wipes) BPXWI30 (Eyelid Wipes) FLEX PUNCTUM PLUG TWIN PACK Comes in six sizes ranging from 0.4mm to 0.9mm This is the best punctal plug I have found in the market. I have placed it in patients who have lost several of the other brands and have yet to lose one. The size must be correct, as always. Works well on punctas that are volcano shaped instead of flat and on lacrimal systems where the canal quickly turns close to the puncta. - Craig Andrews OD FDA now requires these to be sold in packages of two. E30724 (0.4mm Duraplug) E30727 (0.7mm Duraplug) E30725 (0.5mm Duraplug) E30728 (0.8mm Duraplug) **E30726** (0.6mm Duraplug) E30729 (0.9mm Duraplug) Studies show the most common complication with LASIK is dry eye. Solve this and you'll keep a good referral source! DURAPLUG[®] Extended Temporary Plugs These plugs are great for 60-180 days of treatment. Good for post-op dry eye prevention; treatment of seasonal dry eye/ allergies. Great for enhancement of ocular medication. Supplied as 2 sterile inserts in a package and 10 packages per box. These are very popular with LASIK surgeons to reduce the post-op complaints of dry eye following surgery. Be proactive and have a larger referral base for more LASIK procedures! Priced for 20 plugs (\$13 per Plug) E0077 (0.3mm Duraplug) **E0078** (0.4mm Duraplug) E0076 (0.2mm Duraplug) 5

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

SURGICAL & SANITARY


PART #

Т

ALGERBRUSH III - Corneal rust ring remover. Fast turning motor has clutch that is easily stopped by healthy corneal tissue.

PRICE

NAME	DESCRIPTION	PART #		
(A)	Comes with high-quality precision	ACAB5M (0.5mm Burr)		
AlgerBrush II	0.5mm or 1mm burr.	ACAB (1.0mm Burr)		
	Steinless steel hum Comes as a set of 2	ACBR5 (0.5mm Burr)		
(B) Replacement	Stainless steel burr. Comes as a set of 2. (Unsterilized)	ACBR1 (1.0mm Burr)		
	(0)	ACBRSET (1 of Each)		
Burrs	Stainless steel burr. Sold in package of 5.	ACBR5S (0.5mm Burr)		
	(Sterilized)	ACBR1S (1.0mm Burr)		
(C) Replacement	Stainless steel chuck with burr.	ABACHUCK5 (0.5mm)		
Chuck	Sold Individually.	ABACHUCK (1.0mm)		
(D) AlgerBrush	Chuck with 2.5mm Round Fine Grit Diamond	AC25RFD (2.5mm Round)		
Il Pterygium	Chuck with 3.5mm Round Medium Grit Diamond	AC35RMD (3.5mm Round)		
Burrs with	Chuck with 4.0mm Wheel Fine Grit Diamond	AC40WFD (4.0mm Wheel)		
Chuck	Chuck with 5.0mm Wheel Fine Grit Diamond	AC50WMD (5.0mm Wheel)		



STEEL TOOLS STAINLESS NAME DESCRIPTION (A) Miltex Ellis Curved Stainless Steel Foreign Т

Foreign Body Spud	Body Spud, 1.2mm Blade	MIL18380
(B) Miltex Lacrimal Dilator	Medium on one end, Lg. size on the other end	MIL18688
(C) Miltex Foreign Body Needle	Curved stainless steel.	MIL18370
(D) Miltex Davis Chisel- Shaped Foreign Body Spud	Straight shaft stainless steel construction	MIL18372
(E) Miltex Castroviejo Double-Ended Spatula	All stainless steel surgical instrument.	MIL18564
(F) Miltex Bailey Chrome Plated Eye Magnet	One end magnet, other nylon loop.	MIL18398
(G) Magnetic Foreign Body Remover (Dix and	Dix style = straight tip	MFBRDIX
Ellis Style)	Ellis style = curved tip	MFBRELS

STAINLESS STEEL FORCEPS

D.

Ε.

Α.

Β.

C.

Nаме	DESCRIPTION	PART #	PRICE
(A) Bent Tip Serrated Jewelers Forceps	Bent tip, serrated forceps. 4-1/2" long	V707B	
(B) Bent tip Non-Serrated Jewelers Forceps	Bent tip, non-serrated forceps. 4-1/4" long.	V707	
(C) Miltex Jeweler Style Forceps Style #3 (Narrow Fine Tips)	Swiss Jeweler Style Forceps, Non-Magnetic Stainless Steel, Styl 3, Narrow, Fine, 4-3/4"	MIL17303	
(D) Miltex Jeweler Style Forceps Style #4 (Very Fine Tip)	Swiss Jeweler Style Forceps, Non-Magnetic Stainless Steel, Styl 4, Fine, 4-3/8"	MIL17304	
(E) Miltex Swiss Cilia Forceps	3-3/4" long, with a raised angled platform & excellent grasping capability.	MIL181107	
(F) Miltex Beer Cilia Forceps	3-1/2" long, with rounded tips and smooth jaws.	MIL181102	
(G) Miltex Barraquer Cilia & Suture Forceps	4-1/2" long, with 5mm smooth platform tips.	MIL181111	

F.


G.



PRICE

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

SURGICAL & SANITARY



Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE

SURGICAL & SANITARY

BLUMINATOR® OPHTHALMIC ILLUMINATOR

The Bluminator's patented LED-based design revolutionizes examinations of the cornea. The Bluminator produces a bright, spectrally pure beam of blue light, which causes foreign bodies, abrasions, and other epithelial defects to fluoresce vividly during examinations with fluorescein. The blue beam produced is precisely the wavelength required to maximally excite the dye. Since the magnitude of the fluorescence is proportional to both the intensity and spectral purity of the light, the Bluminator yields a clearer, more detailed view of the eye than ever before. The view is enhanced by a 7X magnified image to facilitate detection and diagnosis of corneal injuries and irregularities, streamline contact lens fittings, and speed up treatment.



to find up-to-date pricing and stock status.

5ml weight, surgical-type latex gloves are touch sensitive, & provide

optimal dexterity & resiliency. Unlined, non-sterile, disposable gloves

QUANTITY

100 GLOVES

100 GLOVES

Examination and Procedure Gloves

D

DESCRIPTION

В

LATEX GLOVES

fit either hand. Non-Powdered.

GLOVE SIZES

SMALL, MEDIUM LARGE & XLARGE

SMALL, MEDIUM & LARGE

stine 50

Powdered-Free

Ε

OCONTEO

Due to COVID-19, our examination and procedure glove prices and stock status fluctuate. Call or go online at www.bernell.com

PRICE

ONLINE

ONLINE

G




ASSESSMENTS

EXAM ROOM NEEDS



LOW VISION

REFERENCE





BLU562 - (Blue LED)

BLU562CW - (White Light LED)

Α

B.

A. CAVICIDE

A ready-to-use disinfectant. Effective against TB and MRSA. Disinfect, clean, and decontaminate at one time. Safe on most instruments or on non-porous surfaces even around infants.

В

Cide

CAVWIPE (160 Towelettes)	CAV24 (24 Fl. oz. Spray Bottle)
CAVGAL (1 Gallon Bottle)	4+

B. STAINLESS STEEL INSTRUMENT TRAY

A stainless steel tray with removable stainless steel lid for soaking instruments. GF3254 (Tray Size: 8" x 2.75" x 1")

C. Sporicidin Germicidal Solution

SPC32 (Original Scent (32 oz.)

D. Sporicidin Disinfectant Solution

Gallon of solution for cleaning and disinfecting equipment or surfaces



8

F

Used for pre-operative prepping and irrigation of the ocular and periocular surfaces






NEV

EXAM ROOM NEEDS

ASSESSMENTS

RETINOSCOPY/DIAGNOSTIC LENSES



with magnetic attachments to fit the Bernell Spot Retinoscope head, allowing quick and easy card changes. The magnetic attachments will also fit most newer model streak and spot retinoscopes with magnetic plates. The set features four cards including age appropriate color pictures, black/ white symbols, alphabet card with intentionally missing letters, and a Dolch Word target, all designed to aid in dynamic retinoscopy. This retinoscopy fixation card set is will make checking accommodation easy for practitioner and fun for the patient!

Schematic Eye

This retinoscopy training eye has a sliding axis scale in diopters

from 3D hyperopia to 5D myopia. Cells for three trial lenses.

Adjustable retina position for precise calibration. The Schematic Eye can be tilted to enable observer to precisely sight on optic axis.

Made of plastic

BC2174SB (5cm Base)

BC2174 (8cm Base)

Normal fundus photo included.



These cards cleverly attach magnetically to the magnetic logo of your Welch Allyn retinoscope head, allowing quick, card changes. easy, They provide age appro-

W18250


priate symbols and words based on the Pierce MEM card, for dynamic retinoscopy. These retinoscope fixation cards make checking accommodation easy. The card set includes four washable fixation cards which fit most newer model streak / spot retinoscopes.

Welch Allyn 3.5V Standard **Ophthalmoscope** Head

Halogen light for true tissue color and consistent. long-lasting illumination. 6 apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter. 28 focusing lenses with a range of -25 to +40 diopters.

W11710

Additional Welch Allyn items are available on the web at www.bernell.com.



EXAM ROOM NEEDS

Available In:

5cm Base

8cm Base

RETINOSCOPY/DIAGNOSTIC LENSES


RETINOSCOPY RACKS

Two very nice retinoscopy bars. One is red and one is black. Each has 16 powers of lenses from 0.5D through 4D in half steps and 5-10 in one diopter steps and 12 & 15D. Comes with folding carrying case. Fits easily in a pocket.

DELUXE SKIASCOPIC LENS RACK SET

Contains two iron racks that come in a deluxe bamboo case to protect from dust. One rack contains 9 stationary 22mm diameter plus lenses (ranging from +1.00D to +9.00D) with two sliding 22mm diameter minus lenses with one at +0.50D and the other at +10.00D. The second rack contains 9 stationary 22mm diameter minus lenses (ranging from -1.00D to -9.00D) with two sliding 22mm diameter minus lenses with one at -0.50D and the other at -10.00D.

L122084 (Size 10" x 1.75")



DSLRS8



Gray Standardized Oral Reading Paragraphs

This spiral booklet of 12 pages of oral reading paragraphs in two different font sizes are used for Book Retinoscopy.



70% of the proceeds from this sale go to the Vision Leads Foundation.

ION DIRECTVIEW 4 MIRROR NF (No FLANGE) WITH 2 POSITION HANDLE

Perfect for fast and effective anterior segment examinations! This single element ALL GLASS lens provides four 64° mirrors offer complete 360° view with minor rotation. Has a small 8mm



ALL GLASS, "no fluid" contact surface diameter. Magnification is .9x. Laser spot factor is 1.1x. The 2 position, ergonomic handle is perfect for straight and angled usage. Simply insert for your preference!

> Viewing Range: 4 x 64° Image Magnification: 0.9x Laser Spot Magnification: 1.1x

GRAY STANDARDIZED

ORAL READING

annunnunn

VLORP

PARAGRAPHS

1 Mirror SLT & SLT Flange

This single element, ALL GLASS lens has been developed specifically for angle examinations and Selective Laser Trabeculoplasty (SLT) procedures. The DirectView 1 Mirror SLT Lens has an offset top convex surface that places the image and laser application precisely over the internal reflection viewing mirror with a magnification and laser spot factor of 1.0x.

1 MIRROR SLT + ION070 1 MIRROR SLT FLANGE ION071 +

ION DIRECTVIEW 2 MIRROR



This single element ALL GLASS lens has a top convex surface that places the image and laser application precisely over the 2 viewing mirrors and application area.

ION023

"Flanged" 14.5mm contact surface is ideal for laser treatment.

ION DIRECTVIEW 3 MIRROR NF (No Flange) ION026

ion VISion 3 mirror lenses are of a traditional design, incorporating 3 mirrors angled at 60°, 66° and 78°. An 11.5mm contacting surface facilitates dynamic compression technique. Magnification and laser spot magnification factor is 1.0x.



Perfect for fast and effective anterior segment examinations! Four 64° mirrors offer complete 360° view with minor rotation. No flange has 8mm all-glass "no fluid" surface diameter. "Flanged" 14.5mm contact surface is ideal for laser treatment. Magnification is .9x. Laser spot magnification factor is 1.1x.



AL ASSESSMENTS

RETINOSCOPY/DIAGNOSTIC LENSES

View 20D

ion VISion, Inc.

All-Glass Optical Lenses

1 1 5D



Since 2003, Ion has been committed to providing the highest quality, innovative medical devices to their customers and deliver superior product performance, reliability and customer service.

Their products are developed based upon physician driven needs for innovative solutions.

Ion VISion knows how committed you are to helping patients manage their health. We share your commitment with a passion for providing lifelong solutions that can assist you in managing patients' medical

conditions. They are committed to providing the ophthalmic physician with the most advanced options for treating a broad range of ophthalmic clinical indications. ion VISion's mission is to stay at the forefront of the medical device industry by offering innovative and versatile solutions to today's medical community.

BIO Lenses "Great lenses! The ION BIO lens has excellent optical qualities for viewing the retina, with a wide range of pupil sizes." - Gregory G. Kautz, OD, Optometric Physician Portland, OR

All ion Vision lenses are composed entirely of precision instrumentgrade glass for superior diagnostic and treatment performance.

Their optimized design and exclusive enhanced broadband coating ensures maximum image integrity.

Just one look and you'll recognize the quality advantage that their ezView BIO lenses deliver.

or lens maintenance.

eZView 20D

ZVIan

ezview 220

There's a 20D in every physician's pocket. Our improved design delivers superior image clarity throughout the viewing field. Computer optimized lens surfaces eliminate optical aberrations. A true "20 diopter" design. Available in colors: Black, Blue, Purple, Red or Green.

eZView BIO	15D	20D	20DSS*	22D	25D	28D	28DSS*	30D	40D
Part #	ION006	ION002	ION002SS	ION014	ION007	ION008	ION008SS	ION009	ION010
Viewing Range	30°- 47°	46°-60°	46°-60°	56°-73°	52°-68°	53°-69°	53°-69°	58°-75°	69°-90°
Image Mag.	4.11x	3.08x	3.08x	2.68x	2.54x	2.27x	2.27x	2.15x	1.67x
Laser Spot Mag.	.24x	.32x	.32x	.37x	.39x	.44x	.44x	.47x	.6x
Duites									

eZView | Non-Contact Slit Lamp Lenses



eZView 90D

The standard slit lamp biomicrosopy lens. Their 90D design delivers superior clarity, stereopsis and outstanding small pupil capabilities. Compare it to your 90D and see the difference.

Available in colors: Black, Blue, Purple, Red or Green.

eZView 78D

Ideal general diagnostic and treatment lens. Available in colors: Black, Blue, Purple, Red or Green.



*Steam Sterilizable | SurgiView 20D & 28D ion VISion SurgiView lenses are composed entirely of instrument quality glass — there is no plastic to fog, haze or crack during everyday working conditions. This amazing lens coating is a TRUE steam sterilizable coating which doesn't require complex care instructions

eZView 124SP

"The IonVison eZView 124SP has excellent optical quality for viewing the retina and vitreous, especially in a small pupil. Usually this kind of lens is difficult to use because it can collect moisture easily, but eZView 124SP is very difficult to collect that moisture. This lens is one of the best non-contact fundus lenses through a non-dilated pupil!"
- Akira Negi, M.D., Professor & Chairman Department of Ophthalmology Kobe University Graduate School of Medicine

Enhanced Broadband AR Coating

Our exclusive anti-reflective lens coating is specifically designed to minimize ambient reflection, delivering a pristine retinal image! This unique AR formula is optimized for use with all laser systems.

ADDITIONAL ION LENSES AVAILABLE AT WWW.BERNELL.COM

				-					
	eZView 60D The high magnification offered by the 60D is ideal for detailed observation of the posterior	eZView Non-Contact Slit Lamp Lenses	60D	78D	90D	95D	115D	125D	124SP
w 600	pole. Superior stereoscopic properties.	Part #	ION005	ION004	ION003	ION012	ION011	ION001	ION013
		Viewing Range	81°	97°	89°	98°	116°	124°	124°
		Image Mag.	1.15x	.93x	.76x	1.0x	.76x	.57x	.45x
	NY10- 115	Laser Spot Mag.	.87x	1.08x	1.32x	1.0x	1.3x	1.75x	2.22x
New PS		Price							

eZView

25

REFERENCE

13

LOW VISION



i3D HiMaq

Using superior, medical optical high index-glass and a superior aspheric optical design, the i3D HiMag lens is a next generation, high magnification viewing lens for the Slit Lamp microscope. With great detail, imaging of the Macula, optic disc and pathologies within the central zone will be critically examined with excellent stereopsis under highly magnified conditions.

i3D WideField

The i3D WideField is a superior upgrade to the 90D lens by utilizing high index, high resolution medical grade glass along with a next generation, 3D optical aspheric design. The i3D WideField offers far reaching field of view while maintaining typical image magnification compared to 90D style lens performance.






i3D Ultra WideField

Utilize the i3D Ultra WideField to extend the field of view beyond 90D type lenses. The primary exam lens for routine retina exams on a daily basis. With the Slit Lamp microscope scan beyond typical boundaries compared to other exam lenses, plus this lens is easier to use with patient's having smaller pupils.

i3D UltraWide Non-Myd

This is the lens that answers the call for non-mydriatic pupils or for patient pupils difficult to dilate. May be used as the ideal exam lens for many eye-care professionals due to the superior, specific optical quality material and next generation optical design.





What is an Anti-Reflective (AR) coating? The anti-reflective surface coating is one of the most important components of lens quality. Anti-reflective coating is an extremely thin, invisible multi-layer (5 or more) coating of metal oxide that is deposited on both surfaces of the lens. This multi-layer coating reduces the amount of light reflected from the lens surface allowing more light to reach the retina and thus improving image transmission.

Why do you need AR coatings? Between 4-16% of light is reflected back or lost by a non-coated lens. By adding an AR coating, the overall transmission of the lens can be increased by up to 15%!

i3D Next Gen	i3D HiMag	i3D Widefield	i3D Ultra WideField	i3D UltraWide Non-Myd
Part #	ION091	ION092	ION093	ION094
Field of View	97°	117°	125°	125°
Magnification	1.00x	0.76x	0.57x	0.45x
Laser Spot Factor	1.00	1.32	1.75	2.22
Price				

All i3D Series rings come in your choice of ring color: Black, Blue, Purple, Red or Green.

i3D | Contact Slit Lamp Lenses



OmniView 85

Ultra-high magnification viewing of the posterior pole. Perfect for focal laser treatment of central retinal conditions.



i3D Contact Slit

Lamp Lenses

Part #

Field of View

Magnification

Laser Spot Factor

Price

OmniView 135

OmniView

85

ION019

84°

1.06x

0.94

Designed to deliver the ideal laser spot size for photodynamic therapy procedures. Provides excellent visualization of the choroidal neovascular membrane. Also well suited for focal/grid laser therapy.

OmniView

135

ION018

137°

0 66x

1.50x

OmniView

137

ION015

137

0.44x

2.27



OmniView

145

ION017

144°

0.51x

1.97

OmniView 137

Ideal for use on difficult to dilate patients at the slit lamp. This lens delivers an outstanding retinal view through pupils as small as 3mm.

OmniView 145

Ideal for wide field observation and treatment of the retina. Great balance of magnification and field of view.

OmniView 165

Ultra-wide field imaging to the ora serrata. The OmniView provides an incredibly clear view of the retina from the posterior pole to the periphery without the need of shifting or tilting the lens. Perfect for PRP.

"The OmniView 165 delivers a terrific wide field view for PRP and other procedures. The all-glass contacting element doesn't haze over time like other designs, which means I get a sharp, high contrast view every time I use the lens."

- Peter K. Kaiser, MD Director, Retinal Clinical Research Center Cole Eye Institute Cleveland, OH

Glass is better...

All imaging components in our contacting lens systems are made entirely from glass - even the contacting element - to ensure maximum image integrity. The glass contact element also allows for non-fluid applications, alcohol and peroxide cleaning, and long term durability. Say goodbye to the cracking, crazing and hazing that are typically associated with plastic optics.

These are the ONLY lenses on the market with this exclusive feature. No other manufacturer offers such a wide array of all-glass lenses. It's your money, why settle for less?



OmniView

165

ION016

165

0.5x


2.0x



3 Lens Cases:

4 Lens Cases:

IONCS01 (\$49.95): Lens 1: eZView 15D or 20D or 22D, Lens 2: eZView 95 or 60D or 78D, Lens 3: eZView 90D IONCS02 (\$49.95): Lens 1: eZView 15D or 20D or 22D, Lens 2: eZView 95 or 60D or 78D, Lens 3: eZView 115



OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

EXAM ROOM NEEDS	Occluders						
Ш Z	Long Handle Plexiglass	Occli	uders	LENGTH: 9-3/8 CUP SIZE: 60MI	🖁 Short (Handle	1
MO			10		PLEXIGLA	Iss Occ	Luders
RO				-	Length: 6-3/4	" Cup Size	: 60мм
AM	А. В. С. Д.	Е.	E	G.	(
EX,	2						
လ					H. L	J.	· .
L							
SM			5. 1				
SES							
ASSESSMENTS						•	
	DISPOSABLE OCCLUDERS	LETTER	Түре		ESCRIPTION	Part #	PRICE
AL	White cardboard occluders for use with patients	A	DOUBLE END		AND RED MADDOX	BC207	
OPTICAL	when contagion is a concern.	B C	DOUBLE END		ACK AND PLAIN RED	BC208 BC209	
OP		D	SINGLE END		ED MADDOX	BC209 BC205	
		E	DOUBLE END		20/400) & CLEAR MADDOX	BC207CF	
≻		F	SINGLE END		BLACK	BC204	
RAP	Available	G	SINGLE END		WHITE	BC210	
μ̈́	CALL FOR DETAILS	н	SINGLE END		OX SHORTER HANDLE	BC1200R	
Ē	HP12 (Set of 12)	J	SINGLE END		DOX SHORTER HANDLE	BC1200G BC1200C	
VISION THERAPY			SINGLE END		DOX SHORTER HANDLE		
VIS	MANY OF OUR ITEMS ARE AVAILABLE AT SPECIAL QUANT	TITY DISCOU	NTS OR WITH PF	IVATE LABELING.	CALL 574-259-2070 FO	R DISCOUNTED	PRICING!!
	Long Handle Plas	STIC C	CCULD	ers	Length: 9.5"	Cup Size: 60	
Z					Double End Occ	PART #	
Ĭ			· • • • • •	0	A Red Maddox &	Black	
Ę					B Red Maddo Occluder	X BC205	
NEURO REHABILITATIC	A. B. C. D.	E		FRONT	C Clear Maddo	BC207C	P
HA				F.	Single End Bl	ack	
RE					Cupped Occlu	ider	
IRO					E Black Cupped I Pinhole Occlu		Р
					Length: 7-5/16	"Cup Size :	50mm
~					F Pinhole Occu		Р
z		•	•		Pinhole can be covered w		
LOW VISION			Sho	rt Hand	LE VINYL OC	CULDER	3
Ň			Length: 6-		ATERIAL: VINYL	CUP SIZE: (
Š			Letter	Nаме	DESCRIPTION	Part #	PRICE
Ľ		(igle Pinhole Occluder	Black cupped occluder wit single 1mm pinhole in the ce		
			L Black	Multi-Pinhole	Black cupped occluder wit		
CE		-		Occluder	multiple pinholes.	BC1201B	
XEN				ck and White yl Occluders	Hand-held cupped occlude Available in: Black or Whi	rs. (Black)	
REFERENCE	😈 G. 💽 H. 💽 I. 🖉	•	VIN	y occuuers	Available III. DIACK OF WIT	(White)	
REF	15 Fax Orders to Save Time & Money at	t (574) 25	9-2102 or 2′	103 or Go To	www.bernell.com		

OCCLUDERS

Molded Mask OCCLUDER

Similar design to our mask occluder, but in a molded design. Has a multiple pinhole opening that can flip in or out of place and a longer handle.

RJ025

	Nаме	DESCRIPTION	PART #	PRICE	
А	Multi Pinhole Occluder	Used to determine refractive errors. Has metal spring clips to safe-guard patient's lenses.	BC3010	\$49.95	
в	Clip-On Occluder development Heavy, durable, black plastic 45mm occluder has metal spring clip to safe-guard patient's lenses.		BC3009	\$39.95	
с	Clip - On Translucent Occluder	Similar design to the BC3009, but with a translucent material. Patient can see outlines, but less than 20/400 distance.	BC3009T	\$44.95	


Occlud-A-Measure™

Size: 60mm long

BC7019


Inexpensive device that fulfills seven functions. Use as an occluder, measure inter-pupillary distance, palpebral fissure height, corneal diameter, pupillary diameter, near point of accommodation, and near point of convergence with dot or small clown target.





FRBNR

BERNELL STREFF BI-NASAL WEDGE[™] Updated: Ergonomic design on thicker vinyl material.



MASK OCCLUDER





This is the gold standard for a reversible occluder. Can be used as a mask or an occluder. It features a multiple pinhole opening which can be changed with a flip of a lever to create a full monocular opening.



Translucent Measuring Occluders

This is a set of three occluders made of flexible acrylic. Each is labelled with the level of acuity expected when looking through the occluder. Useful for determining the level of acuity required to hold binocularity in high exophores. Also good in detection of latent nystagmus.

Important Note: These occluders are designed for use when blurring at distance and are not for use in near acuity testing.

Set of 3: 20/100. 20/200. & 20/400 BC214SET





Translucent Occluder

Allows viewing of the covered eye by the doctor during cover test while keeping the distant acuity of the patient below 20/400.



NEURO REHABILITATION

LOW VISION

REFERENCE



OPTICAL

REFERENCE

LOW VISION

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com



This eye patch that lets the light in, but keeps the patient fully occluded even at near. This patch is for use when doing fields such as the Humphrey Field Analyzer or the FDT. This patch reduces the light adaptation difference.

BFEP (Bag of 6)



These patches allow the patient's eye to be seen during patching. The patient sees 20/400 or less at distance through the patch to allow peripheral fusion during patching.

BTEP (Bag of 6) BTEP12 (Bag of 12) Improved translucency & flexibility vs our original.



Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 18

REFERENCE

EXAM ROOM NEEDS


ASSESSMENTS

OPTICAL

VISION THERAPY


NEURO REHABILITATION

LOW VISION



19 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

PRISM BARS



	10, 20, 23, 30, 33, 40 & 43
B.) Extended Bar - AHB10	1, 2, 4, 6, 8, 10, 12, 15, 20 & 25
C.) Better Bar - AHB6	1, 3, 5, 10, 15 & 20
D.) Good Bar - AHB5	3, 5, 10, 15 & 20

Vertical Prism Bars

В

С

D

ΝΑΜΕ	PRISM DIOPTERS	PRICE	E
E.) Complete Bar - AVB15	1, 2, 3, 4, 5, 6, 8, 10, 12, 14,16,18, 20, 25 & 30		
F.) Better Bar - AVB10	0.5, 1, 1.5, 2, 3, 4, 5, 6, 8 & 10		A 1997 Print
G.) Good Bar - AVB6	0.5, 1, 1.5, 2, 2.5, 3		G

EZ VIEW VERTICAL AND HORIZONTAL PRISM BARS

PRISM BAR SET	
EZ View Horizontal and Vertical Prism Bar Set (with Case)	
EZHVSETC	

*Combination set shown

PRISM DIOPTERS	PRICE
Horizontal Sizes: 1, 2, 4, 6, 8, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, and 45 Diopters.	
Vertical Sizes: 1, 2, 3, 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 25, and 30 Diopters.	

The EZ View Prism Bar eliminates a problem that exists with other prism bars: the numerical magnitude of the prism has typically been located on the thick edge of the prism bar. This created a problem for the examiner. When using the prism bar to measure either the angle of deviation or to measure

convergence and divergence amplitudes, the examiner had to either turn the prism bar sideways to view the numerical magnitude or make a pronounced movement to the side to view the prism magnitude. When moving the position of prism magnitude with EZ View the examiner can now

see the numbers without the need to either turn the prism bar or change viewing position. The issue of turning the prism bar sideways during the assessment of vergence amplitudes could interfere with the accuracy of the measurement. These issues can be particularly problematic for less experienced examiners; EZ View's design allows clinicians to concentrate on the important aspects of the examination procedures and easily view the results.

Prism Bar Features: Improves vectographic results. Reduces washout. Eliminates reversal. Superior diopter accuracy.



POWERS PRINTED ON FRONT OF PRISMS INSTEAD OF SIDE!

INDIVIDUAL PRISM BARS	PRISM DIOPTERS	PRICE
EZ View Horizontal Prism Bar GOEZHP	1, 2, 4, 6, 8,12,14,16,18, 20, 25, 30, 35, 40, & 45	
EZ View Vertical Prism Bar GOEZVP	1, 2, 3, 4, 5, 6, 8, 10, 12,14,16,18, 20, 25, & 30	

REFERENCE



PEDIATRIC

ADULT

Pediatric Fixed PD's: 52, 54, 56, 58, 60

(Sold in Assorted Colors)

CNTFCH (Pediatric)

REFERENCE

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

(Color: Gray)

CNADCH (Adult)

Adult Fixed PD's: 60, 62, 64, 66, 68

Includes Case

<u>Available in:</u> +/- 0.25,+/- 0.37, +/-0.50, +/-0.75, +/-1.00 & +/- 2.00

Same Jackson Crossed Cylinders pictured above, but in larger case.

Includes: +/- 0.25D, +/- 0.50D, +/- 0.75D, +/- 1.00D

in foam lined aluminum case.

Jackson Crossed

CYLINDER SET TOM3030SET

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



22

REFERENCE

PUPILLARY REFLEX DEVICES



Bernell Pupillometry Exam Kit

This combination kit was developed to serve as an all-in-one assessment kit for pupil measurement. This handheld circular pupillometer device is used for measuring different size pupils. It was developed and researched by Dr. Jack Richman, O.D., as part of the Massachusetts Drug Evaluation and Classification Program (MDEP). This pupil measurement device is used by law enforcement officers across the country to assist in determining whether people are possibly impaired and/or under the influence of drugs.

The reverse side has expected pupil size in different lighting conditions for non-impaired people. Each kit includes (2) pupillometers plus a modified LED penlight that produces the same intensity light as an incandescent bulb. Scan QR code to read the comparative study on the modified LED penlight vs. an incandescent penlight, which was performed by Dr. Jack Richman, OD, FAAO, FCOVD.

Each Kit includes (2) of the MDEP Pupillometers. plus the modified Bernell MDEP LED Penlight.

Penlight Specifications:

Silver LED Penlight Reusable LED penlight clips to pocket. Uses (1)

AAA battery (not included). LED light provides

CLEAR PUPIL GAUGE WITH

This pupil gauge was designed by Sharon Alamalhodaei.

COMT. The clear sturdy material makes it easy to

measure pupils accuarately and quickly, and includes

a handy reference chart of abnormal pupil findings.

Material: Thick, clear plastic

LED100

consistent lighting and LED lasts for years.

On/Off button on back.

ŝ

ETTPC1

- \rightarrow Material: Black, metal housing \rightarrow Switch: Click On/Off button on end of device
- Size: 4-3/8" long by 1/2" diameter
- → Battery: 1 AAA Alkaline (not included)

Uses (2) AAA batteries (Included). SUPER BRIGHT LED PENLIGHT WITH CLIP This super bright LED penlight is much brighter than a standard penlight. This is an excellent light source for foreign body inspection when the Cobalt Blue Penlight (ALPENF) is added. It provides a bright target light when used with a Maddox Rod or similar test at distance. Lightweight reusable LED Uses 2 "AAA" Cell batteries. LED lasts for a lifetime. Push button switch located on back.

UV PENLIGHT This penlight is ideal for pupil examinations of a patient **NEURO REHABILITATION** with a dark iris. The UV light causes the lens of the eye Length: 6.21" Body Dia.: 0.375" (15.17 x 0.95 cm) Pupil Characteristics Chart LOW VISION Density Bar Has 6 neutral Density levels from 0.3 to 1.8 log units. Used in conjunction

prevent dust from gathering between the lavers. It comes with a handle and leather-like pouch for storage.

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

G15119

MEASUREMENT DEVICES

Pocket CFF Tester™ €

The **critical flicker frequency** testing device is a compact, hand-held device for determining the CFF values of a patient's vision. Critical fusion frequency is the point at which a light is flickering so fast that the visual system perceives it as a solid or non-flickering light.

In certain neurological conditions, the CFF of a patient will change. The measurements may be taken on the right eye, left eye or both eyes.

The CFF unit will give you three values:

- 1) The point where the solid light first starts to flicker is referred to as non-seeing to seeing or NSTS.
- 2) The point where you are seeing the flicker and the frequency increases till you no longer
- see the flicker, which is referred to as seeing to non-seeing or STNS.
- 3) The average of the two values above.

RAF (Royal Air Force Rule) Convergence Rule


RAF Rule (Royal Air Force Rule) provides a binocular gauge to measure Objective and Subjective Convergence as well as Accommodation in 1 mm increments. The RAF Rule consists of a 50 cm long rule with a slider holding a rotating four-sided cube.

RAFCR1

Each side has a different target. The first has a vertical line with a central dot for convergence fixation. The others provide a limited number of lines of near reading examples.

Placido Disk, Clip, & Handle

This economical keratoscope is used to determine the curvature characteristics of the anterior surface of the cornea.



Exophthalmometer

Made of clear plastic, exophthalmometer has millimeter scale on both sides to check for parallax. Luedde design. **3-1/4**" long. Fits in coat/shirt pocket.





Duane Age) with sliding cardholder.

Accommodation & Near Point Rule

The Accommodation and Near Point Rule measures

50 cm long, 10 mm square, it is engraved with Centimeters,

Inches, Diopter Scales and Duane Age (the mean of the

near point of accommodation at different ages according to

Dark Grey

The near point rule is the point closest to the bridge of the

nose where convergence and binocular single vision can

no longer be maintained as an object approaches.

Accommodation Convergence Rule
Two Cards printed on Styrene

GR50B

Kit Includes:

Slide Card Holder

Made of gray opaque

plastic with white

, markings.

Used to test & train monocular, binocular, accommodation, and rock. Trains convergence using physiological diplopia as the control on suppression. Features four scales; dioptric values, duane age, centimeters and inches. Made of durable plastic with engraved black numbers.



VISION SCREENING

• BC11931 - Snellen 20 Ft Acuity Chart

• BC1284150 - Hyperopia Flipper +1.50

• CVT1 - Color Vision Testing Made Easy

• BCDI - HOTV Distance Acuity Chart

· BCNIE - HOTV Near Card with Elastic

• BCVTPK2G - Pre-Kindergarten Test Book

• BC816AMAE - AMA Near Test Card w/ Elastic

· VA1010S - Butterfly Acuity Test

DEMNV - DEM Test and manual

• BC200 - Bernell-O-Scope

· BSROD - 9" Brock String Rod

• BCF402 - Recording Forms

• BCFIX - Fixation Stick

Manual

• BC204P - Occluder

BERNELL PRELIMINARY TEST KIT FOR YOUR OFFICE Save Over \$70 with this kit! BPTKIT



Often primary care optometrists are being forced into seeing more patients per hour or per day in order to survive economically with third party plans. Unfortunately, this sometimes means either foregoing certain tests or delegating them to technicians. Sadly, particularly with children, the assessment of the binocular system is one of the first things to be passed over in a routine exam. In order to free up the doctor's time and yet get the necessary information of the binocular skills for the doctor, Bernell has developed a preliminary room test kit which may be used by a technician. It is designed around simple free space tests instead of an instrument such as a stereoscope to allow a more realistic test of real world demands.

• BCBS548 - Worth 4 Dot flashlight with Letters which is used to measure fusion & convergence abilities. BCFIX - Fixation stick to test pursuits & accommodative skills & NPC.

 BCSALCD - Saladin test card quickly & easily measures vertical & horizontal phorias & fixation disparities. • BC1275VF - Vergence flipper with 6 Base-Out & 3 Base-In lenses equally before the eyes. BC1284150 - Hyperopia Flipper +1.50 for testing farsightedness

· DEMNV - DEM test of eye movements to assess visual processing & eye movement skills • Manual - Simple manual written for the technician

878

0

BERNELL DELUXE SCHOOL SCREENING KIT Save Over \$200 with this kit!

т

14444444444444444444444444

0

E D - 4

ECFD

z

. 3

5

BSSK

After years of working with different optometrists and ophthalmologists, Bernell has developed a kit with for screening of the major visual skills neces-

sary for the primary grade student. Tests cover visual acuity at distance and near, convergence and accommodation, stereopsis, color testing, saccadic eye movements, phorias and fusion, suppression,

and hyperopia. Kit comes with instruction manual that explains how and why of each test. Designed for use by school nurses or persons with some vision screening training.

BROKEN WHEEL TEST[™] FOR NEAR

This convenient test appears as a near rotary card for the phoropter, but is instead handheld. One can rotate the card to show two windows in a column. The patient simply selects which car has the broken wheel. The correct answer appears on the back of the card so that the administrator of the test can easily score the response. A chart that is included gives an approximate visual acuity based on the farthest distance at which the patient consistently and correctly responds.



SHERIDAN GARDINER TEST The Sheridan-Gardiner Test contains both

ADDITIONAL KITS AVAILABLE ONLINE



BWT

nt_111

::: ::

::: :::

A.G. A.A.

near vision and distance vision tests according to the Faculty of Ophthalmologists and Snellen Standards. This Test also includes a reduced Snellen Test. Testing depends on matching shapes rather than identifying or naming letters. All are reversible. The Sheridan Gardiner Test is primarily designed for children but is also usable for deaf mute or other patients unable to interact with traditional letter based visual acuity tests. Results are more reliable than illiterate E or similar tests. A series of orthoptic booklets in washable plastic finish are included with the basic kit to augment the range of letters in each size.

o b 00 <u>____</u> -

BROKEN WHEEL TEST^{***} 66 6 6 6 10 2

FOR CHILDREN WHO HAVE TROUBLE WITH THE TUMBLING "E" TEST Easily administered, valid and reliable test of visual acuity in children 2 to 5 years of age, the visually or mentally handicapped, the multiple handicapped, deaf, or brain damaged, of any age. Uses a clinical approach of testing for visual acuity by incorporating the Landolt C. The Broken Wheel Acuity Test utilizes cards that have a familiar, non-threatening symbol (car), presented in a forced choice response. The simple recognition of

the gap in a Landolt Ring is the critical feature. The acuity level in this test has been established to be equivalent to Snellen Letter optotype. No interpretation of the directional orientation of the optotypes is necessary, since both cards are the same except for the critical feature. Responses can be by finger pointing, hand pointing, or by eye pointing. Includes a pretest for response reliability. 18 cards. VA from 20/100 to 20/10.

Broken Wheel Test[™] - Advanced Guide

Detailed guide for utilizing the broken wheel test cards. Includes training exercises with accompanied illustrations. Letters, screening worksheets and glossary are also included.



OPTICAL

VISION THERAPY

25

LOW VISION

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

VISION SCREENING

OPTEC PLUS SMART VISION SCREENER FROM STEREO OPTICAL

Your practice requires a new generation of vision screening devices.

The Stereo Optical OPTEC PLUS all-in-one digital vision screener enables quick, thorough and repeatable results in the pre-test room, saving valuable time in the exam lane.

New features include a user-friendly interface, an expanded library of tests with limitless customization including auto-mode, randomization and protocols.


- Testing distances: Near, Far and intermediate
- . Monocular and Binocular
- Lighting conditions: day / night
- Constant and homogeneous luminance levels with con lighting in accordance with the international standard

- Multiple language options available •
- No training or certification needed
- Lightweight and portable
- Tablet, laptop, or personal computer (PC)

EASY TO USE AND INTUITIVE INTERFACE

 Control your test on a single screen 	22
 Choose your test conditions in 1 click: Night & Day 	1 E M
 Monocular LAP, Briocular Par, Internetiale, Maie datactor salectori Pandori 	23600
· Change your selection at anythre	
OPTEC PLUS PG enables running tests under	111 A
the plane condition	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1

THE OPTEC PLUS TEST LIBRARY IS DESIGNED TO FIT ALL YOUR NEEDS



restation, Visual accely/ latters / runitizers / Landstin / Tanalan / E-Tunising / ETDPE ing teacher Long of another Refractive Teats: American RestGener, Astgrister, Labort (s Non-refractive Teste: Coor Inds.) Contrast newsbydy test expla. Funcar, latenti and vertical prepis, D

Additional tests: HOTV tests, Road sign recognition, Alen test

tal portpheral vision tests Donal at 60°, 70°, 80°, 80°, 90°, 100° and 35° temporal, 43° c re receivery heat: Seconds after which subjects can see the dots post exposure to give IOPTEC PLUS PG.

CUSTOMIZATION: CREATE YOUR OWN PROTOCOL IS SIMPLE



AUTOMATIC MODE": INCREASE THE EFFICIENCY OF YOUR PATIENT FLOW AND SAVE TIME IN THE EXAM LANE



- Self administered by the polyint or
 Delegated to your staff.
- Text to speech technology guides patient through the text, hearing you from repetite
- Present auto-mode protocols
- The OPTEC PLUS records the patient

CONNECTIVITY: GET ALL THE ADVANTAGES OF A SMART DEVICE

- Dealer records in pdf like that can be printed. cled and exported
- · EMR compatibility to some your data transfer
- Operated and montorestor: PG Laptop Tablet
- · Wretes contectivity: Bluetoon and Will



ntrolled	and the second s

*Only available in North America

OPTEC+

Additional Vision Screeners are available online. VISIT WWW.BERNELL.COM FOR MORE INFORMATION.



CRECIFICATIONS

Dimensions:	19" x 11.15" x 17.95"	
Weight: Input: Dutput: Compliant to:	14.7 Lbs 100-240V , 50 - 60 Hz 24 DC, 2.1 A (50W) CE, FDA, FCC and ISO 16284 standards	
Classification:	Class 1 Medical device without measurement function	
Option available:	Headphone, ancillary lenses, peripheral mask	

equirements:

- ocessor Dual core 1,8GHz minimum
- crosoft windows XP Sp3 minimum, 7 or vista pro or nilial, 8 & 8.1 (x86 or x64) 10, 8 RT not compliant
- m 2Gb DDR2 minimum
 - reen 1280x768
- o hard disk for the software
- oort USB 2.0 High speed or Native Bluetooth

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.



NEURO REHABILITATION


ASSESSMENTS

OPTICAL

VISION THERAPY

26

REFERENCE



LOW VISION

REFERENCE

27

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

P a + MANNEL B - NO.

ORDERS OF 500+ CAN BE IMPRINTED ON EITHER PAD WITH YOUR NAME & PHONE. PLEASE CALL FOR PRICING.

5" X 6.5

50+





Sized to fit over our Adult Optokinetic Drum (DAL300) this vinyl sheet is used to test for the optokinetic reflex. The banner has child interest targets similar to our DAL301. Can be rolled up for easy storage. Clear end pieces hold banner from curling. Washable plastic material similar to window shade. Velcro on backside of the banner allows for easy attachment around the Adult Optokinetic Drum (DAL300).



Pelli-Robson Letter Near Contrast SENSITIVITY CHARTS (SET OF 3)

Optotype Options: PATTI PICS, PV NUMBERS, or SLOAN

Features Traditional Sloan Optotypes found in the Original Pelli Robson Chart! For use at 16 in (40cm) Includes three charts, each of which is scramble. Choose from traditional Sloan Optotypes, Patti Pics Pediatric Optotypes or PV Numbers Optotypes.

Patti Pics - Set of 3 PV5011PP Acuity Level 10m PV Numbers - Set of 3 (20/630) PV5012N at 40 cm (16 in). **SLOAN Letters - Set of 3** PV5013S

TEST

This easy and convenient test can be used to monitor a decrease on the contrast sensitivity function over time. This test also demonstrates the difficulty of seeing in poor lighting or at night. Conditions affecting Contrast Sensitivity include: Age, Glaucoma, Cataract, IOL, Keratomas, Diabetes, Macular edema, Optic neuritis, MS, Parkinsons, Papilloedema. Made of washable material.

RULES OF THUMB WITH THIS TEST:

. Level 5-8: Suffer significant loss and need contrast enhancement Level 9-12: Have noticeable loss and need lighting evaluation • Level 13-16: Have a near normal contrast sensitivity



L


HVCST

н E

HAMILTON-VEALE CONTRAST SENSITIVITY

• Level 4 or Less: Suffer from severe loss or are legally blind

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com



28

Dimensions: 16" × 16" × 1"

accuracy and reliability.

Chart Comes with: Test

Instructions, Recording Forms, Soft-tipped Pointer, 2 meter

cord Conversion Charts

0

1074, 3.074, 1.2074

zкрус

PV2180


The test chart is also bound and includes instructions and recording forms.

Each level uses traditional Sloan optotypes in an ETDRS format for an


equally difficult test from beginning to end. Every level uses a different

series of optotypes to help reduce any possibility of memorization. Also, each contrast sensitivity level is highly calibrated and will produce accurate


results. This flip book is an excellent choice for research because of its



ACUITY - DISTANCE



ACUITY - DISTANCE



Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com




EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY



	~ ~-		-			-	
PROJECTOR SLIDES (STEREO OPTICAL)							
Ι.	J.	Κ.	L.	М.	Ν.	0.	P.

Slides fit most projectors. Vectographic slides do require a polarizing screen. To learn more about the slides listed, you may visit our website at www.bernell.com. They are listed under Primary Care-Projector Slides & Screens.



EXAM ROOM NEEDS

OPTICAL

REFERENCE

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

RPS

RSLPNV1

BERNELL

EXAM ROOM NEEDS


ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



Set includes four paddles, three with gratings printed on both sides and one solid gray face for presentation.

FACE DOT TEST Developed by Jack Richman, O.D.

This tests for visual acuity in difficult-to-test pre-verbal or nonverbal infants and children with low cognitive or developmental functioning levels. The face is an easily recognized target and provides a strong attention stimulus for low functioning patients. The test employs the use of contrast sensitivity in a forced choice testing procedure. The results may be converted into an approximate visual acuity based on the testing distance.

The set includes three paddles that be used at different test distances to approximate visual acuity. Distance chart is on back of test paddle. Instructions for the different testing methods are also included.



ACUITY - NEAR



COMPUTER & MULTIFOCALS CARD - 2ND EDITION ALLOWS YOU TO CHECK IMPORTANT COMPUTER & MOBILE DISTANCES IN THE EXAM CHAIR

Specifically designed for patients requiring corrective multifocal evewear. Card may be placed at different heights and distances to simulate the complete computer screen. Card is printed on two sides - one side shows a computer screen with Microsoft Word document with different font sizes. Simulates a 15" computer monitor with shiny screen. It includes highlighted and colored words. The second side shows a mobile tablet and phone with text to stimulate real world applications.

BCCRT2 (8.5" x 11" Card)





"BEE" CONTINUOUS Text Near Card

Continuous Text Card for testing patients 10 and up. Simple sentences in proportionally spaced (logMAR) lines with Snellen equivalents. Acuity range 8.0M to .4M (20/400 to 20/20). Pediatric card tests children with elementary reading levels (grades 1-3). Two-sided, measures 18x23cm (7x9in).

LIBEEC



both adults and children. Range of print sizes with visual acuity notations in Snellen, Jaeger, Decimal, etc. Also shows indicated magnifications & reading ages. Includes a wide range of contemporary near tasks: Music, computer, phone book, newsprint, share listing,

"Bird" Continuous Text Near Card

Continuous Text Card for testing patients 7 and up. Simple sentences in proportionally spaced (logMAR) lines with Snellen equivalents. Acuity range 8.0M to .4M (20/400 to 20/20). Pediatric card tests children with elementary reading levels (grades 1-3). Two-sided, measures 18x23cm (7x9in).

.0.

0000

00000

00000

0 0

Chart Size: 12" x 9" x 1"

0000

LIBIRDC

PV2508

w

5

LEA SYMBOLS & LEA NUMBERS NEAR VISION CARD

This card combines a LEA Symbols® line test on one side and a LEA Numbers® line test on the opposite side. Proportionally-spaced lines from sizes 20/400 to 20/10 (6/120 to 6/3) equivalent. Attached 16-inch (40cm) cord ensures proper testing distance.




PATTI PICS NEAR VISION CHART

Memorization can alter testing results and create an inaccurate record for the patient. This near vision chart features multiple groups for repetitive testing help minimize memorization. This helps to create a more accurate and reliable vision chart that yields repeatable results. Each side is different to minimize memorization.

Near at 40cm Testing Distance.

Intermediate at a Choice of 66cm or 80cm Testing Distances. Both have cord for easily achieving proper testing distance. Near Vision Testing at 40cm: 8M (20/400) to .20M (20/10) Intermediate Vision Testing at 66cm: 8M (20/250) to .20M (20/6.3) Intermediate Vision Testing at 80cm: 16M (20/400) to .40M (20/10)



SLOAN MULTIPLE GROUP NEAR VISION TESTING CARD SINGLE SIDED

Two smaller groups 2M to .20 for repetitive testing (20/100, 6/30 to 20/10, 6/3 equiv.).

Memorization can alter testing results and create an inaccurate record for the patient. This near vision chart features multiple groups for repetitive testing help minimize memorization. This helps to create a more accurate and reliable vision chart that yields repeatable results.

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

ACUITY - NEAR

EXAM ROOM NEEDS


ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION




ALL CHARTS ARE PRINTED ON RIGID WHITE VINYL. Order Now by Phone at (574) 259-2070 or Go To www.bernell.com REFERENCE


Acuity - Near


NEAR POINT TEST CARD


Designed for use in performing the MEM (Monocular Estimate Method) test for determining lag of accommodation. The front of the card also has a cross-cylinder grid and a reduced Snellen chart. The back of the card provides additional reading material of various type sizes. The card is made of a plastic material which is easily cleaned.



BC11981 (4.5" x 5.5")







음

ensees → 3.5 € = V □ Value, cave fix, and protect your eyes. Make each and every one of those

mages as sharp and beautiful as they can be

37

LOW VISION

120052305055

÷.

SMD TEST CARD

A durable home monitoring test card for those patients needing to detect metamorphopsia. Especially useful with ARMD and PMF. The

back of the card contains near point reading material of various type sizes. This allows the patient to monitor possible changes in visual

acuity. The card is made of a washable plastic material.

REFERENCE

A DEPENDENCE OF A DEPENDENCE OF ii. 14 140 i. GENERAL INSTRUCTIONS old card at approximate out right and find, by co Ny 127 N vering left eye. the electron ν н ta are note BCNI (5.5" x 8") BCNIE (Elastic Cord) BC11980 (5.5" x 8") JAEGER TEST CARD A.M.A. Test Card Front side has 2 through 12, back side 14 through 20 and music line. 0.37 to 2.00 size on side one; 2.5 to 8 size on side two. Card has letters from 14/14 to 14/224 size in blocks of letters to test crowded conditions without word recognition. AMA RATING READING CARD the focal point can be advanced with a --N 10 10 10 10 convex lens 10 122 101 102 101 100 100 100 2111 concave -..... m 1111 1111 1111 ν 100 1111 1111 1111 HII den maleran a real faithile 1111 1111 1111 1111 and Filebert & F F JAEGER'S TEST TYPES BC816AMA (5.5" x 8") BC815 (5.5" x 8") SNELLEN X-CLY JAEGER CARD A.M.A. ENGLISH READING CARD Near point test card printed on both sizes. Front side tests for Jaeger 2 through 8. Back side tests Snellen and X Cylinder card is made of A.M.A. notation and visual efficiency 14/14 to 14/84 size on front. Distance meter notation 0.50 to 2.00 size on back. durable plastic and is washable. 500 1995 CARLANSIUS BAT 41.04 11 -11 1.010 14 100 . 1.00 47 th States of States of States BC1196871 (5.5" x 8") BC814SCJ (5.5" x 8") Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 38

HOTV NEAR CARD

This card is ideal for testing illiterate patients. It has HOTV letters on one side & Allen figures on the other. **BCNIE has elastic cord 16**"

т

long to verify testing distance.

BINOCULAR SKILLS

EXAM ROOM NEEDS


ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



REFERENCE 39

BINOCULAR SKILLS

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION





REFERENCE

LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



LOW VISION





REFERENCE

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

Introducing the upgraded Stereo Tranaglyph[™] Test Booklet, a comprehensive solution for screening binocular vision deficits in patients of all ages. With improved features like tabbed title pages, optimized test plate flow, and durable pages, this booklet makes testing simple & effective.

BC1025

RANDOM DOT

NEAR

Booklet Dimensions: 5.75" x 8"

Materials: Opaque Rigid Vinyl with Clear Spiral Binding



COLOR VISION TESTING



PIP represents pseudoisochromatic plates, the type of test images used to identify color vision defects and the fact that the book has 24 plates in total. The PIP24 screens for protan (red weak), deutan (green weak) and tritan (blue/yellow weak) color vision deficient individuals. This test is the perfect choice when testing both adults and children because it has a section to test each demographic. This product is accepted and used by the U.S. Navy, Army, Marines, Coast Guard and FAA.

Waggoner CVT1 **Color Vision Testing Made Easy (CVTME)** OLD NEW Woggener



Color Vision Testing Made Easy (CVTME) is the pediatric gold standard for identifying genetic (red/green) color vision deficiencies in children as young as 3. It consists of one demonstration plate and 9 test plates displaying circle, star and/or square throughout the plates and can be completed in under 60 seconds.



Waggoner Computerize d Color Vision Test (CCVT)

The Waggoner CCVT is the military-grade, color vision diagnostic application designed by Dr. Terrace Waggoner and validated by the Naval Aerospace Medical Institute (NAMI). Additionally, it is accepted by the U.S. Navy, Army, Coast Guard, and Marines for testing their pilots. It includes the pediatric Color Vision Testing Made Easy and D-15. The Waggoner CCVT is an all-in-one color vision testing suite to satisfy anyone interested in testing for color vision deficiencies ranging from school nurses to the U.S. Military.

Features:

- Extended Color Vision Test: Diagnostic
- Self-administered and self-scoring
- Adult and pediatric tests
- CPT 92283
- Approved codes: glaucoma, ARMD, cataracts, lupus, optic nerve neuropathies, diabetes, hypertension Cross platform versions: Windows, Android, and iOS



Waggoner CCVT Tablet with App Pre-Installed WCCVT1 Waggoner CCVT App Only						WCCVT				
	Genetic Deficiencies	Acquired Deficiencies	Screening Test	Diagnostic	Pediatric Test Included	CPT 92283	Cross Platform	Accepted by Militaries	Test Time (min)	Price
Waggoner CCVT Computerized Color Vision Test	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	3-6	
CCT ¹ Cone Contrast Test	\checkmark	\checkmark	\checkmark	\checkmark	×	\checkmark	×	\checkmark	5-6	
CAD Colour Assessment	\checkmark	\checkmark	\checkmark	\checkmark	×	\checkmark	×	\checkmark	10-12	
ColorDx CCT-HD Cone Contrast Test HD	\checkmark	\checkmark	\checkmark	\checkmark	×	\checkmark	×	×	5-6	
Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 46										

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

ASSESSMENTS

OPTICAL

VISION THERAPY

THE DYSLEXIA SCREENER

A practical & easy-to-use screening method for detecting accepted types of dyslexia (specific reading disability). Ideally suited for the regular classroom teacher, resource specialist, speech /language specialist, physicians & optometrists. This quick test provides a method for early detection of dyslexia so that further testing & therapy can be instituted.



DX1 (CHILDREN GRADES 2-9)

DX2 (Adult)




Dyslexia Determination Test (DDT) (DDTKIT (KIT)

Explains the language encoding and decoding paradox, determines learning disorders while identifying specific dyslexic patterns and offers approaches to therapy.

Replacement materials are available.

THERAPY IN DYSLEXIA AND READING PROBLEMS



By: Dr. John Griffin It features newly presented material for identifying deficiencies. Describes therapy to relieve reading failure. 105 pages.



Wold Sentence Copying Test

Standardized visual motor test designed to be used as a screening device in schools. Easily administered in a group setting. Used to determine if a child has the ability to rapidly and accurately copy a sentence from top to bottom of page.

• Pad of 50 forms

- Includes laminated instruction sheet
 - WOLDSC

QUICK NEUROLOGICAL SCREENING TEST 3R (QNST)




The QNST-3R includes updated norms (based on a nationally representative sample), along with revisions and clarifications to the Record Form, administration procedures, and scoring guidelines. The QNST-3R is an individually administered, empirically based assessment of the development of motor coordination and sensory integration (both of which have been shown to relate to learning as well as to general daily functioning), seen as neurological soft signs (NSS). NSS are minor irregularities that include poor motor coordination, sensory perceptual changes, and difficulty sequencing complex motor tasks. The presence of NSS can indicate neural

ideo dan indicate neural						
• 1 Manual	Kit Includes:	trauma and ar often harbinger				
• 25 Record Form • 25 Remedial Gui Forms (developme	deline	of learning dif- ficulties (without history of trauma).				

Product	Part #	Price
COMPLETE KIT	QNST3R-K	
Manual	QNST3R-M	
25 RECORD FORMS	QNST3R-F	
25 Remedial Guideline Forms	QNST3R-RGF	

TEST OF AUDITORY PROCESSING SKILLS, FOURTH EDITION (TAPS-4)

The TAPS-4 provides information about language processing and comprehension skills across three intersecting areas: phonological processing, auditory memory and listening comprehension. These areas underpin the development of effective listening and communication skills, and are critical to the development of higher order language skills, including literacy skills.



The TAPS-4 features new subtests along with revisions to subtests from the TAPS-3, fully updated norms, and an expanded age range. The TAPS-4 subtests were also reorganized into Index and Supplemental subtests, reducing testing burden and increasing flexibility in administration.

The TAPS-4 also features audio administration for the subtests in which proper pronunciation of speech sounds is critical, providing a greater degree of standardization and accuracy during the testing process.

The TAPS-4 has 11 subtests, organized into three indices: Phonological Processing Index, Auditory Memory Index and Listening Comprehension Index.

• 1 Manual • 25 Test Booklets • Administration CD	Kit Inc	ludes:
Product	Part #	Price
COMPLETE KIT	ATP2216	
Manual	ATP2217	
25 TEST BOOKLETS	ATP2218	
ADMINISTRATION CD	ATP2219	

Test of Visual Motor Skills - 3rd Edition

By: Nancy A. Martin, PhD



This TVMS assesses how well a person can coordinate visually guided fine-motor movements to copy a design while it is in sight. The TVMS-3 is used to determine whether there are any systematic distortions or gross inaccuracies in the copied design that could be the result of deficits in visual perception, motor planning, and/or execution. The test utilizes a simplified scoring system and no longer uses a protractor. A single test replaces the two levels of the 2nd version (TVMS-R and TVMS-UL), making the test easier to administrate and more economical. The TVMS-3 was normed on 2610 individuals.

Kit Includes: Manual, 15 Test Booklets & 15 Record Forms					
Product Part # Price					
COMPLETE KIT	EXKC519				
Manual	EXKC520				
(15) TEST BOOKLETS	EXKC521				
(15) SCORING FORMS EXKC522					

47

NOISIN WO

DEVELOPMENTAL EYE MOVEMENT TEST" (DEM")

Bernell's exclusive DEM[™] test provides an objective method of assessing fixational and saccade activity during reading and non-reading tasks. Incorporates a sub-test of number calling in a vertical array and provides the means to evaluate oculomotor function with numbers in a horizontal array. Developed by Dr. Jack Richman and Dr. Ralph Garzia. Comes with the Version 2 Examiner's Manual which includes instructions as well as 60+ pages of tables of standards based on age and grade level. Some case histories and clinical treatment considerations are presented. Also included is scoring software which displays the comparisons to standardized scores as shown in the 60+ pages of tables, plus provides a printable record. The software saves time and prevents errors in reading the many tables. Pad of evaluation sheets (25) included.



We sell so many more DEM[™] tests than the older saccadic tests that we have decided to no longer sell the others. See our website for studies which show why we think the DEM[™] is the best.


The Mobile Scoring Software available through Google Play or iTunes App Store.

If you are interested in having a DEM scoring software available on your phone or tablet, Bernell has partnered with MakeMind, a developer in Australia, to release a new scoring app for the DEM known as the DEM Calculator. This app is available for purchase through iTunes App Store or the Goggle Play Store. Visit the Bernell website for info on the new scoring mobile app available for iTunes and Android.



Product	Fart#	Price		
Developmental Eye Movement Test™ Kit Includes Scoring Software, Test Plates, Examiner's Manual & (25) Score Sheets	DEMNV			
DEM Examiner's Manual	DEMMANUAL			
DEM Scoring Software on USB (Windows Only)	DEMSTD			
DEM Replacement Test Plates	DEMTP			
DEM Score Sheets (Pad of 25)	DEM1			
Original plates may become grimy over time due to extensive use and can				

Original plates may become grimy over time due to extensive use and can be replaced. Replacement plates are for sale only to accounts who have purchased the full test. Replacements do not come with any instructions or test standards.



Test Kit: Manual, 25 Record Forms, 25 Design Completion Worksheets, 1 set of 5 Stimulus Cards, in portfolio

The VSA provides an easy way to screen for common visual skill difficulties that can impact academic performance and participation, including reading and writing tasks by way of five items that assess binocular, ocular-motility, and visual-motor skills. The VSA requires minimal equipment and can be administered as part of a school's annual vision screening program or with students who present with visual skill concerns to determine need for further optometric evaluation.

Product	Part #	Price
VSA-2 - Complete Kit	VSA2276	
VSA-2 Manual	VSA2277	
VSA-2 Record Forms	VSA2278	
VSA-2 Design	VSA2279	
Completion Worksheets	VGAZZIS	
1 Set of 5 Stimulus Cards	VSA2280	

STRUCTURED OBSERVATIONS OF SENSORY INTEGRATION-MOTOR (SOSSI-M) Assess proprioceptive and vestibular processing, and motor planning. Ages: 5 to 14 Testing Time: 20-40 Minutes

The SOSI-M is a new nationally standardized assessment of sensory-based motor skills, based in Ayres Clinical Observations, and developed by well-regarded experts in sensory integration.

The SOSI-M includes 14 sets of administered items that assess proprioceptive and vestibular processing, and motor skills.

<u>SOSI-M Kit Includes:</u> Administration Manual, Spanish Administration Manual, Quick Start Administration Guide, 25 Record Forms, 25 Comprehensive Observations of Proprioception (COP-R) Forms, one 15 cm and 20 cm square, tennis ball, pencil and pencil topper, and stopwatch. Tote bag and Airex foam mat also included.

L	Product	Part #	Price
L	SOSI-M - Complete Kit	SOSI2281	
L	SOSI-M Manual	SOSI2282	
L	SOSI-M Record Forms	SOSI2283	
L	SOSI-M COP Forms	SOSI2284	
	SOSI-M Quick Start Administration Guide	SOSI2285	
	SOSI-M Spanish Administration Instructions	SOSI2286	

IREST READING SPEED TEST

The IReST differs from prior text charts that used single sentences by employing linguistically standardized paragraphs. They close a gap in the diagnostics of reading performance. The paragraphs



Chart Size (Inches): 8.5×5.5 in

Chart Distance: 40cm / 16in

are designed to resemble everyday life situations of reading continuous text, such as in books or newspapers. The level of difficulty corresponds to sixth grade reading (10–12 years), which is comprehensible for teenagers and adults. Due to lower variance, measuring reading time of a complete paragraph rather than a single sentence or random words is more reliable. Furthermore, it can provide some information about fluency, fatigue, and mistakes.

- Easy to use
- Easy to handle
- Easy to measure (by stopwatch)
 Highly calibrated for reproducible results
- Easy to calculate (normal values and formula included)



Characteristics:

- · Set of 10 paragraphs of text for repeated measurements
- Comparability between languages for international studies
- Natural conditions: Newspaper print, meeting the demands of everyday life

Available in other languages. Go online to www.bernell.com for more info.

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

Wachs Analysis of	Cognitive Structures	(WACS) [®] Test Kit
Developed By: H. Wachs, O	D & L.J. Vaughn, MA	

WACSKIT COMPLETE KIT WITH MANUAL



WACSKITRB

(WACS)© Test Kit - Additional 15 Workbooks

Standardized inventory based upon Piaget's theories of measuring cognitive development in terms of body and sensory thinking in children ages 3 years to 5 years 11 months. Individually administered and easily scored in 30-45 minutes by assessment professionals or classroom teachers. Focuses on period during which child begins to acquire symbolic thinking and has been used successfully with both verbal and non-verbal children, as well as adults, whose development age ranges in the 3-6 age period.

A product of 30+ years of clinical research and 5 years of data collection, this inventory is a valuable tool for assessing appropriate type and level of activities to produce effective developmental and learning growth. Valuable instrument for all working in child development and especially useful with learning disabled, autistic, mentally challenged, culturally disadvantaged and as a predictor of school readiness. The test comes as a complete kit with manual. Blocks are made of Corian material (except for green) which is very rugged and safe for teething. Easily cleaned between testing unlike the older wooden blocks which could harbor bacterial growth.

The kit consists of: 15 clusters of task divided into subtext made up of a variety of manipulative, visual and sensory materials, as well as body movement tasks grouped into 4 areas of assessment: Identification of Objects, Object Design, Graphic Design, and General Movement.

(WACS)© Test Kit - Profile Sheets - Pad of 25 WACSKITPAD

Manual WACSKIT-M

Test of Information Processing Skills (TIPS) Executive Function, retention, recall, and memory capacity are easily detailed. Author: Raymond E. Webster, PhD

The TIPS provides clinicians with quick and reliable measures of how well a person processes information (letter strings) presented visually and auditorily. Although TIPS had its genesis in the Learning Efficiency Test (R. Webster, 1981, 1992, and 1998), it features new items, new subtests (Delayed Recall and Word Fluency), new scoring procedures, and new national norms (N=3,314).

Product	Part #	Price
TIPS - Complete Kit	TIPS8468	
TIPS Manual	TIPS8469	\$79.95
TIPS 25 Protocols	TIPS8471	
TIPS Stimulus Card Booklet	TIPS8751	\$60.95

- Visual & Auditory Memory/Sequencing
- Ages 5 through Adult (90+)
- Individual Administration
- Nationally Norm-Referenced on 3,314 individuals
- Administration Time: 30 minutes · Qualification: Level B

Test Kit Includes:

- Manual
- Stimulus Card Booklet
- 25 Protocols
- **Test of Handwriting Skills-R (THS-R)**

ce

Visual-motor integration difficulties can foreshadow learning problems. Author: Michael Milone, PhD The THS-R is an untimed, clinical assessment of neurosensory integration skills evident in handwriting (both manuscript and cursive) that are often disrupted in students with learning difficulties. The THS-R provides standardized assessment of children's handwriting ability for both manuscript and cursive styles that can be used with any of the popular handwriting programs (such as Palmer, Zaner-Bloser,

•		
Product	Part #	Pri
THS-R - Complete Kit	THS8439	
THS-R Manual	THS8440	
THS-R 15 Manuscript Test Booklets and 15 Record Forms	THS8444	
15 Cursive Test Booklets and 15 Record Forms	THS8445	
THS-R Training Video	THS8441	

- · Visual-Motor (sensory) Integration
- Ages 6-0 through 18-11
- · Individual or Group Administration
- Norm-Referenced
- Testing/Scoring Time: 25 minutes
- · Qualification: Level B
- Test Kit Includes:
 - Manual

D'Nealian). The THS-R is not a classroom assessment of penmanship.

- 15 manuscript test booklets
- 15 cursive test booklets
- 30 Record Forms Training Video

HANDBOOK OF COPY FORM CARDS (BC11985 (Set) Diagnostic Tests These cards provide a special test for eye-hand and visual motor coordination, visual organization and form reproduction. Written By: R.W. Lowry, Jr., O.D. This book is also referred to as the Monroe III Test. RL9513 Perceptual skills program tests for outstanding char-REPLACEMENT SCORE PADS acteristics in children ages 2-1/2 to 10 years. Shows Includes: 7 cards and pad of 50 score sheets manual skill, symmetry, proportions and orientations.

Extras and/or

Replacements:

LOW VISION REFERENCE

TEST OF VISUAL PERCEPTUAL SKILLS - 4TH EDITION

The TVPS-4 is the latest update of the standard comprehensive assessment of visual analysis and processing skills. The TVPS-4 remains an easy-to-use assessment for determining visual-perceptual strengths and weaknesses. Norms are based on a nationally representative sample. Additional lower-level items were added to address the needs of younger or more impaired individuals, and norms now extend through 21 years of age, making the TVPS-4 useful for a wider range of examinees.

The TVPS-4 now includes 18 items in each of seven perceptual areas. Ceilings minimize fatigue and ensure that more difficult items do not unduly tax younger students, while allowing for reliable and accurate assessment across the age range. The TVPS-4 is consistent with current methods of visual information processing used in occupational therapy and optometry.

MOTOR FREE VISUAL PERCEPTUAL TEST - 4TH EDITION

The MVPT-4 is the most recent revision of the only non-motor visual perceptual assessment that can be used throughout

the lifespan. The MVPT-4 provides a quick, reliable, and valid measure of overall visual perceptual ability in children and adults. The MVPT-4 includes 45 items from the MVPT-3 which have been reorganized and grouped for easier administration. Stimuli are comprised of black-and-white line drawings and designs, with answer choices presented in an easy to record multiple-choice format. No motor involvement is needed to make a response, making the test particularly

Price

An updated and streamlined version of the classic visual-perception test.

Part #

ATP04XA

ATP042M

ATP042XD

ATP042XB

The	TVPS-4	also	measures	the	following
Catt	ell-Horm	-Carro	oll theory n	arrov	v abilities:

By Ronald P. Colarusso, EdD/ Donald D. Hammill, EdD

- Visualization
- Flexibility of Closure
- Visual Memory

Product

COMPLETE KIT

TEST PLATES

REPLACEMENT MANUAL

25 RECORDING FORMS

Memory Span

Product Price Part # COMPLETE KIT **TVPS4-K** SCORING FORMS TVPS4-F TEST PLATES TVPS4-TP MANUAL TVPS4-M

useful with those who may have motor disabilities.

The MVPT-4 is one of the most widely used visual perceptual assessment for recertifying adult drivers

The MVPT-4 is designed to be used for screening and

after head injury or stroke.



(pack of 25), Test Plates, and Manual

OPTICAL

ASSESSMENTS

EXAM ROOM NEEDS

NEURO REHABILITATION

REFERENCE

50

	_
	4
	C
	-
	S.
	-
	-
	5
	~

JORDAN 3RD EDITION FEATURES:

10

- ► Visual Reversals
- Individual or Group Administration
 - ► Ages 5 thru 18 ► Norm-Referenced
 - ► Testing/Scoring Time: 20-30 minutes
- - ► Qualification Level B

THE REVERSALS **FREQUENCY TEST**

Written by Richard Gardner MD. Useful in diagnosis of children with learning disabilities. It offers an objective assessment of letter and number reversal frequency. The test is helpful to differentiating neurologically based and psychogenic learning disabilities. It is standardized on both normal and learning-disabled children.

Three tests: reversals execution: reversal recognition: reversals matching.

GRFT

JORDAN LEFT-RIGHT REVERSAL TEST (3rd Edition)

research purposes by psychologists, occupational therapists, educational specialists, optometrists, and others

who may need to determine a person's overall ability to discern and understand visual stimuli.


The Jordan Left-Right Reversal Test-3rd Edition helps clinicians identify students who have difficulty with reversals, and what types of reversals are problematic: objects, letters, numbers, words, or letter sequences. Reversals are a major stumbling block for students learning to read. The Jordan-3 may be beneficial to add to a battery of test when assessing AD/HD.

The Jordan-3 has two parts (children ages 5-8 take on Part 1) consisting of five subtests, in developmental sequence. The first and last subtests are new. Included in the Test Kit are a Laterality Checklist (an informal survey to determine preferential use of one side of the body) and a Remedial Checklist that provides suggested remediation activities.

Kit Includes: Test Kit (Manual, 25 Record Forms, 25 Remedial Checklists, 25 Laterality Checklists)

ATP2538A (KIT) ATP2538BU (RECORD FORMS) TORDAN LAPT-REART REVERSAS TEST or Street Manhood Manhood of Taxa Text ORDAN LEFT-RIGHT REVERSAL TEST Emmon





Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 52

SUPPLIES & EQUIPMENT



Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

53

OW VISION

OPTICAL


NEURO REHABILITATION

METAL DISTOMETER WITH CASE BY GULDEN OPHTHALMICS

This product is used to accurately measures the vertex distance between a patient's cornea and the lenses in eyeglasses or trial lenses. It is sold in a case with a conversion wheel for calculating lens power relative to the distance from the eye.

A change in vertex distance of eyeglasses compared to that of the trial lenses results in a change in the effective power of the eyeglasses. This becomes especially important when the eyeglass prescription is greater than a plus/minus four or five diopter correction.

DST127



PROGRESSIVE LENS IDENTIFIER

Very important tool that all optical labs need for fitting and testing progressive lenses. Greatly assists in the guick and easy location of the various manufacturers' marks placed on progressive lenses by using the specially designed background contrast and lighting provided. Also useful for checking the surface quality for scratching or bubbling when coating is applied to lenses.



PGLID3



DIGITAL ULTRASONIC CLEANER

(MULTIPLE FRAME CAPACITY)

A five cleaning cycle model with a larger tank capacity to accommodate two or more frames at the same time. Removable lid keeps cleaning solution free of dust when not in use. Stainless steel tank. One year warranty.

Power: 110/120 V AC; 60 Hz Ext. Dimensions: 9.0" x 2.7" x 6.3 Tank Capacity: 3 pints (1,400 ml) Cycle Settings (Seconds): 90, 180, 280, 380, 480 Weight: 3.7 lb (1.7 kg)



ULTRA BLUE (ULTRASONIC CLEANER CONCENTRATE)

Ultra Blue is a new generation water soluble liquid ultrasonic cleaner concentrate.

cleaning action without the use of any solvents or other volatile organic compounds (VOC's). The Ultra Blue formulation is rapidly biodegradable and safe for the environment. An 8 ounce bottle of concentrate makes 4 gallons of cleaning solution.



AM1011



AM1007

Provides accurate millimeter measurement for all types and sizes of lenses quickly. Tape is durable, easy-to-read White plastic with 0-300 calibrations. ABS plastic housing.



USOLS

USOR

REPLACEMENT TAPE

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 54 ASSESSMENTS



BRAIDED EYEWEAR RETAINERS

ER20GR+

**Specify Color When Ordering




Revolutionary new eyewear retainers will help you get a "Gripz" on your expensive sunglasses and prescription eyewear. Gripz eyewear retainers are engineered from high-grade braided silicone that will stand up to the elements and toughest use and abuse. In fact, they've been lab tested to hold the weight of a brick, earning the seal of "Brick Tested Tough." Whether you're biking, fishing, rafting or participating in any other vigorous outdoor activity, the patent-pending Gripz retainers uniquely combine friction and a contracting force to provide an extremely secure hold that will keep your valuable glasses safe.

- · Lab Tested to Hold the Weight of One Brick earning the seal of "Brick Tested Tough"
- Engineered from High Grade Braided Silicone Which Allows the Gripz to be very Heat Resistant
- Fits MOST Frames
- Designed to have a secure hold and protect your eyewear under the harshest conditions.

AVAILABLE COLORS: BLACK ASH, THINK PINK, LAVA ORANGE, SAFARI GREEN, SURF BLUE, & DUSK GREY

EXAM ROOM NEEDS

ASSESSMENTS





The prism lenses in glasses angle vision allowing patients to lie on their backs to read or watch TV without eye fatigue or strain. Frame design permits use with corrective glasses. Approved by optometrists, ophthalmologists, and orthopedic physicians.

Wt: 5.4 oz. Comes with case and cleaning cloth.



Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

CONTACT LENS PRODUCTS

JUMBO CONTACT LENS CASES (PKG. OF 50)

Ensure that your office is stocked with extra contact lens cases for your patients! These jumbo cases for soft contact lenses have wells 10mm deep x 22mm diameter. They are available in White, Turquoise or Yellow. We also sell a smaller Light Blue case that is 7mm deep.

Part #	Dертн	Color
NLSCLCY	10mm	Yellow
NLSCLC	10mm	WHITE
NLSCLCBB	10mm	BRIGHT BLUE
NLSCLCB	7mm	LIGHT BLUE

QUANTITY	PRICE
5+	
2 - 4	
1	

Please specify color when ordering



These tweezers are 3.15" long

non-returnable item.

and work great for handling soft contacts. Comes sterile

in a plastic sealed bag. For hygiene reasons, these are a

These tweezers are 2.25" long and work great for handling soft contacts. Comes sterile in a plastic sealed bag. For hygiene reasons, these are a non-returnable item.

LOW VISION

57 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

White polypropylene tweezer 5" long x 3/4" wide.

For hygiene reasons, these are a non-returnable item.

OPTICAL ASSESSMENTS

VISION THERAPY

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



contact lenses.

VWPLT24

polishing. Unit last for years.

M1800

VPSOAP

CLPOL

Readers & Sunglasses

EXAM ROOM NEEDS

ASSESSMENTS


OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE




FLIP-UP MAGNIFYING LENSES

CLIP AND FLIP MULTI POWERED CLIP-ON

The Clip & Flip magnifying lenses turns any pair of regular glasses into a "flip-down" hands free magnifier. Flip the magnifier up for normal viewing or down for a super magnified image. The Clip & Flip magnifying lenses are perfect for reading, fly-tying, hobbies, crafts, model building and so much more! The crystal clear acrylic lenses will give you a sharp, distortion free view that allows you to see the finest details with ease. Each Clip & Flip magnifying lens comes included with a protective soft pouch to protect the lenses from dust and scratching.



Magnifications	Part #	Price
1.5x power (+2.25 diopters)	OD10CML	
1.75x power (+3.00 diopers)	OD12CML	
2x power (+4.00 diopters)	OD14CML	



HALF EYE READERS

"See Wees" Infant Sunglasses

Designed to fit infants and toddlers up to approximately 4 years of age. Soft, flexible no screw frames provide comfort and safety. Polycarbonate lenses offer 100% UVA and UVB protection. Certified lead free.

Frame Size (mm): 42/17/105

AM1821 (Pink) AM1824 (Navy)

Wee2Cool Weezers™ Toddler Suanglasses

Designed to fit children approximately 2-6 years old. They have soft rubber frames for safety and comfort. Polycarbonate lenses offer 100% UVA and UVB protection. Certified lead free.

Frame Size (mm): 48/13/105

AM1035



Computer Eyes

Designed for clear viewing of video display terminals. Converts ordinary bifocals into glasses that relieve eye strain and relieve neck fatigue associated with VDT operations. Eliminates the need for additional glasses. Made of quality optical plastics. Plastic coated clip prevents scratching. Flip up to see distance without removing.





These slip-ins were developed by an optometrist for use of his patients in VT programs. They work well for computer work or as spare readers for those on the verge of bifocals. The slip-ins are similar to post-myd slip-ins, but have a prescription of +1.12D power. Made of polycarbonate with UV protection and a scratch resistant coating. Come with handy display that can also be attached to monitor to hold peekers when not in use.

Suggested resale price to patients is: \$29.95 (PCR125)

Make-up Magnifying Glasses

Make your patients happy by selling this device. Now women with presbyopia can apply their eye makeup/mascara accurately. Lens comes +1.50 or +2.00. Lens swivels at nose from side to side. Comes in Gold.



"It amazes me how many we sell of these each year in my clinic! You just have to ask and they will sell!" - Craig Andrews, OD

3M Press-On[™] Bifocals

Similar to Press-On Prisms. Come as a set of two bifocals. No adhesive necessary. **Can be re-used**.

GREAT FOR TRIALS OF STRONG BIFOCALS ON

Low Vision Patients



PART #	PRISM DIOPTERS	PRICE	
B+ (Power)	+1.00, +1.50, +2.00, +2.50,		
D+ (FOwer)	+3.00, +4.00, +5.00 & +6.00		
Description of the second s			

*Specify diopter size when ordering such as B+1.50 for +1.50D plus lens.



MANUFACTURING THE WORLD'S SAFEST AND MOST FLEXIBLE FRAME!

We take great pride in boasting about all the safety features in our frames. FlexFrames offers 17 sizes, 13 models, and 21 colors of our Flexible and Safe frames; ensuring the largest selection from which your patients may choose. Below are important facts that you should know about our products:

- **Proposition 65 Compliant**
- Can hold a +/- 28 prescription
- **BPA Free**
- **Phthalate Free**
- **Rubber Free**
- Latex Free
- **Metal Free**
- **Holds Prisms**
- We cover EVERY age group

For FlexFrames, there's no such thing as prodigal spending when it comes to children's safety. Our Flexible and Safe line has many SAFETY attributes. The novelty about these frames is that there are no metal parts and the frames are so soft that children and adults can LITERALLY sleep while wearing them without causing injury to themselves or ruining the frame. And with these frames, no-one can pinch themselves even while sleeping. Many children and adults love that they can see as soon as they wake up. They don't have to worry about searching for their frame to be able to see. The anatomically designed bridge eliminates the need for nose pads and the inherent risk associated with them. The elastic band behind the users' head between the ears assures the frame stays on the users face. We designed the tensors so that they take up any elastic slack, and there is no leftover elastic to cut or unravel.



Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

iccolo

ROUND SHAPE

Available sizes 33 (Item #: FFPI33 + Color Code)





PICCOLO 33 (33 - 11 - 112)

Ages: 8 Months - 2 Years





BC Clear Pink



Coming Soon

FC Clear Beige



IN Orange



JC Clear Grey



L Lavender

Coming Soon

S Neon Purple







HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate



REFERENCE

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



RASCAL 39 (39 - 14 - 121) RASCAL 42 (42 - 14 - 125) Ages: 2 - 5 Years

Available sizes 39 & 42 (Item #: FFRA39 or FFRA42 + Color Code)





60

EC Clear Blue



l Red





KP Burgundy Pearl





BC Clear Pink



D Dark Blue



FC Clear Beige



IN Orange



JC Clear Grey



L Lavender

Coming Soon

S Neon Purple





DS Navy Blue



HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate





LOW VISION

BAMBAM 31 (31 - 11 - 112) BAMBAM 34 (34 - 15 - 112) Ages: 3 Months - 2 Years

OVAL SHAPE Bambam

Available sizes 31 & 34 (Item #: FFBB31 or FFBB34 + Color Code)







EC Clear Blue



l Red



J Dark Grey



KP Burgundy Pearl





BC Clear Pink



D Dark Blue

Coming Soon

FC Clear Beige



IN Orange



JC Clear Grey



L Lavender

Coming Soon





DS Navy Blue



HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate



REFERENCE

Nt0

OVAL SHAPE

PLUTO 38 (38 - 17 - 117) PLUTO 40 (40 - 18 - 120) Ages: 2 - 6 Years

Available sizes 38 & 40 (Item #: FFPL38 or FFPL40 + Color Code)



B Pink

CP Royal Blue



EC Clear Blue



I Red



J Dark Grey



KP Burgundy Pearl





BC Clear Pink



D Dark Blue



FC Clear Beige



IN Orange



JC Clear Grey



L Lavender

Coming Soon



BS Magenta



DS Navy Blue



HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate



LOW VISION

P Plum **S Neon Purple** Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

EXAM ROOM NEEDS

TAZ 37 (37 - 14 - 115) TAZ 44 (44 - 16 - 125) Ages: 1 - 7 Years

OVAL SHAPE

Available sizes 37 & 44 (Item #: FFTZ37 or FFTZ44 + Color Code)







EC Clear Blue



l Red



J Dark Grey



KP Burgundy Pearl





BC Clear Pink



D Dark Blue



FC Clear Beige



IN Orange



JC Clear Grey



L Lavender

Coming Soon





DS Navy Blue



HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate




EXAM ROOM NEEDS

ASSESSMENTS

66

S Neon PurpleVC Clear GreenOrder Now by Phone at (574) 259-2070 or Go To www.bernell.com



EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE



Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

68

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE

FLEX FRAMES

EXAM ROOM NEEDS

ASSESSMENTS


OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION





Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

FLEX FRAMES



SPORT SHAPE

GROOT 42 (42 - 16 - 118) Ages: 3 - 5 Years

Available sizes 42 (Item #: FFGR42 + Color Code)



B Pink

60

CP Royal Blue



EC Clear Blue



l Red



J Dark Grey



KP Burgundy Pearl





BC Clear Pink



D Dark Blue



FC Clear Beige



IN Orange



JC Clear Grey



L Lavender





BS Magenta



DS Navy Blue



HP Yellow Pearl



JS Black



N Frost White



M Milk Chocolate



P PlumS Neon PurpleVC Clear GreenFax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com


71

NEURO REHABILITATION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL



72 Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

FLEX FRAMES

EXAM ROOM NEEDS **MEASUREMENTS** FlexFrames **PATTERN FLEXIBLE & SAFE** ASSESSMENTS A: HORIZONTAL ED: EFFECTIVE DIAMETER DBL: DISTANCE BETWEEN LENSES C: CIRCUMFERENCE B: VERTICAL **Frame Style** Eye Size А В ED Temple Bridge DBL Circumference Length Size Aspen 48 48 49.96 31.02 51.46 138 16 14.47 132.68 Aspen 50 50 51.98 34.36 54.78 142 19 17.26 142.45 Aspen 53 53 56.82 34.76 58.50 142 19 17.45 150.48 **OPTICAL** Bambam 31 31 35.44 28.37 35.90 112 15 12.99 106.32 Bambam 34 34 37.87 29.89 38.26 115 15 12.99 106.32 23.53 17 Brainy 39 39 42.13 43.40 124 15.15 114.49 44.9 Brainy 42 42 26.25 46.50 130 16 14.75 123.35 45 48.04 27.51 49.72 17 15.03 Brainy 45 126 131.03 47 49.41 28.56 17 15.95 Brainy 47 51.24 128 134.41 **VISION THERAPY** Cameron 51 51 52.31 40.14 57.64 140 19 17.34 152.40 Groot 42 42 45.00 29.59 46.64 118 16 14.89 121.96 Piccolo 33 33 35.72 32.56 36.52 112 11 13.17 109.39 Pluto 38 38 41.56 34.04 41.60 117 17 14.49 119.70 Pluto 40 40 43.78 34.73 43.86 120 18 16.09 124.63 41.29 41.29 Rascal 39 39 30.33 41.88 117 14 116.84 Rascal 42 42 43.8 31.6 14 12.67 123.23 44.40 117 Sherlock 39 39 42.06 27.05 43.92 124 14 12.57 115.58 15 Sherlock 42 45.43 28.87 48 125 14.23 124.20 42 **NEURO REHABILITATION** Sherlock 47 47 50.43 36.27 53.24 128 17 15.58 142.51 37 40.29 14 Taz 37 29.74 41.04 115 13.17 110.88 Taz 44 44 46.03 34.55 47.04 130 16 14.57 127.71 17 Yogi 40 40 42.52 33.72 45.46 112 15.20 125.66 17 Yogi 42 42 45.03 35.71 48.42 130 14.46 132.89

17

15.54

137.76

STRAPS

37.01

49.84

130

46.59

44

Yogi 44



BUILT UP BRIDGE FRAMES

Currently Available in the following styles and colors

Bambam (33 eye size) Brainy (39, 42 & 45 eye size) Vogi (38 eye size) Taz (37 eye size)



REFERENCE

LOW VISION

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com



Erín's World Frames


SPECS4US was created and founded in 2004, by Maria Dellapina, a single mother of four with more than 30 years of experience in the optical field. The name SPECS4US stands for Superior Precision Eyewear for Children who are Special.

Maria's youngest daughter Erin was diagnosed with Down syndrome at birth. After several frustrating years of searching for an eyeglass frame that would fit Erin properly, and finding out that nearly 87% of children with Down syndrome will need glasses by the time they enter preschool, Maria decided to combine her passion as a mother, with her knowledge of eyeglasses, and SPECS4US was formed.

Erin's World, is a frame line designed to accommodate and fit the features of children and adults with Down syndrome, as well as, other individuals with unique facial features. This innovative eyeglass design includes a bridge which is adjusted to properly fit a low nasal bridge and temples (or "arms" of the glasses) which are modified to keep the glasses from constantly slipping down.

The SPECS4US mission is to improve the sight and quality of life for children and adults with Down syndrome and other unique facial features by providing custom frames that fit their unique needs. Erin's World frames are now available in a variety of sizes and styles.





How is the Erin's World frame line different then a typical frame?

Erin's World frames are custom-designed with shortened temples, extra wide frame fronts, and a lowered bridge. Temples are either 180° spring hinges or TR90. Specs4Us recognized the need for a frame line that gives a solution to individuals with low nose bridges.

Our frames have been proven to increase the quality of life as they are the first thing our customers put on in the morning which is why our motto is "See Something New"! Having a frame that fits properly and is comfortable for all day wear is why we continue to be one of the best frames in specialized eye wear.



EXAM ROOM NEEDS

ASSESSMENTS

NEURO REHABILITATION

VT STARTER KITS

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE



• BC1011S - Red/Green Barrel Card (1) • BCVT - Guide To Vision Therapy™ Manual

VT STARTER KITS

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

Bernell Expanded VT Starter System (Better Starter Kit)

Designed with a larger assortment of products. Utilizes the most advanced techniques in vision therapy in one complete system.

Save Over \$350 with this kit!

INCLUDES

- · BC50K, BC500K, & BC600K Tranaglyph Kits
- BC801 & BC925 Tranaglyph Trainers
- BC1050BK Aperture-Rule™ Trainer
- BC1270 +/- 1.00, 1.50, & 2.00 Accommodative Flippers
- BC11484 & BC11488 4D and 8D Adult Prism Training Goggles
- BCVTA Vision Test Book
- BCBS538 Worth 4-Dot Flashlight
- BC405IO Mirror Stereoscope
- BC1003 Dual Polachrome with Lamp
- BC200VPT Bernell-O-Scope™ w/ Variable Prismatic Trainer
- BC906 6.5" x 10.5" Red/Green Stereo Trainer
- BC63RG Large Red/Green TV Kit
- BC200LLI & BC200LLO Base-In and Base-Out Landscape Line Training Cards • BCHC - Hart Card Set
- BC10901 Single Phy-Dip Cord Ball Trainer
- BC275TRAS & BC550TRBS Transparent Slides (1 Each)
- BC1011S Red/Green Barrel Card (1)
- BCVT Guide To Vision Therapy™ Manual



Pro Starter System (Best Starter Kit) Bernell VT

Utilizes the most advanced techniques in vision therapy in one complete system. Contains all products in "Better Start Kit" as well as many additional popular therapy tools.

SAVE OVER \$850 WITH THIS KIT!



- · BC50K, BC500K, & BC600K Tranaglyph Kits
- BC801 & BC925 Tranaglyph Trainers
- BCB6-10 Black Background Tranaglyphs with R/G Flipper
- BC930K Variable Tranalgyph (930 Series)
- BC11 Necklace Variable Tranaglyph
- BC275TRAS & BC550TRBS Transparent Slides (1 Each)
- BC1003 Dual Polachrome with Lamp
- BC906 6.5" x 10.5" Red/Green Stereo Trainer
- BC63RG Large Red/Green TV Kit
- BCBS538 Worth 4-Dot Flashlight
- VS1-5 Cheiroscopic Tracing Pads (VS1, VS2, VS3, VS4 & VS5) • VTP-1060FGD - VTP Edition Floor Rotator with Geo Pegboard • BC200VPT - Bernell-O-Scope™ w/ Variable Prismatic Trainer

INCLUDES

• BC200VI & BC200VO - Base-In and Base-Out Character Training Cards

• BC1270 - +/- 1.00, 1.50, & 2.00 Accommodative Flippers

• BC113810 - 10D Child Prism Training Goggles

BC200LLI & BC200LLO - Base-In and Base-Out Landscape Line Training Cards

• BC11484, BC11488 & BC114812 - 4D, 8D & 12D Adult Prism Training Goggles

- BCVTA Vision Test Book
- MARSBL2RG Bernell Soft White Red/Green Marsden Ball
 - Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 76
- BC1050BK Aperture-Rule™ Trainer • BC10901 - Single Phy-Dip Cord Ball Trainer
- BC1011S Red/Green Barrel Card (1)
- BCHC Hart Card Set
- BC675C Macula Integrity Tester Version 1

• BC405IO - Mirror Stereoscope

LOW VISION

VT PATIENT KITS

ERNELL VT PATIENT KIT **BVTPKIT1** 12+

Help bring vision therapy to your patient's home by providing them with this comprehensive kit!



We are excited to introduce the Bernell VT Patient Kits! These kits were designed to offer a variety of essential vision therapy tools in an economical package and to act as a starter kit for offices to distribute to new patients starting a vision therapy program. Elevate your patient experience with the Bernell VT Patient Kit!

This vision therapy kit includes an expertly-curated selection of VT equipment that your patients can use at home to complement their inoffice vision therapy. These kits come packed in a convenient carry tote bag with additional room so that practices may add additional items to the kit bags as their patients progress through therapy. Also included on each bag is a clear vinyl bag tag so that office may feature their business card or logo to brand the bags for their patients.

KIT INCLUDES:

- · 10 ft. Physiological-Diplopia Cord (Brock String)
- · Marsden Ball with Hook
- · Free Space Fusion Cards with Holder
- Translucent Occulder with Toggle Adjustment
- · Peripheral Awareness Chart with Letters/Figures
- +/-1.50 Accommodative Flipper
- 65mm Plastic Lens (-3.00)
- "LifeSaver" Cards (Translucent and Opaque)
- · Standard and Wide Straw with Red Pointer

ACCOMMODATIVE PROCEDURES

· Soft Carrying Case with Bag Tag


- · Red Pointer
- · Flashlight
- Saccadic Pencils (2)
- · Bernell Bull's Eye Target
- · Custom Scribble Page
- · Custom Hart Chart
- · Parquetry Shapes

KIT

- · Parquetry Workbook
- · 4 Chart Drill

77





WACHS KIT - HOME VT KIT BY DR. HARRY WACHS

Industry standard materials

BERNELL

WACHSKIT The Dr. Harry Wachs Vision Training Kit is a comprehensive vision

Affordable and accessible

X

therapy tool developed by Dr. Harry Wachs and updated and refined by Bernell. It is based on a hierarchical sequence of training that has been proven to be effective over many years. The kit includes patient and doctor manuals of instructions and is suitable for primary care doctors who do not practice visual therapy and do not have easy access to a doctor who does. Home-based vision therapy can help many patients with problems that cannot be corrected with glasses alone and it does not require a lot of time.

KEY FEATURES

Convenient and comprehensive Integrated and instructional

OPTICAL

EXAM ROOM NEEDS

VT PATIENT KITS

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



Petrosyan Vision Training Kit

The Petrosyan Vision Training Kit is a collection of essential vision therapy equipment recommended by Dr. Tamara Petrosyan, OD, an associate clinical professor at SUNY College of Optometry. This kit is suitable for patients starting vision therapy and as they progress through therapy. It is only available for purchase by vision care professionals but can be shipped directly to patients. The kit includes a variety of items and for more details, customers are encouraged to contact the company's customer service team.

Petrosyan Vision Training Kit Includes:

- **Progressive Freespace Fusion Cards**
- Petrosyan Oculomotor Workbook
- Petrosyan Visual Search and Scan Workbook
- Red/Green Reading Unit with Narrow Bars
- Red/Green Goggles with Elastic Band
- Free Space Fusion Card Translucent & Opaque
- Life Preserver Card Translucent & Opaque
- +/- 1.00 Accommodative Flipper
- +/- 2.50/-3.50 Accommodative Flipper
- Four Corner Hart Chart
- Hart Chart Peripheral Awareness Chart
- Barrel Card
- 20ft Brock String
- 4PD Loose Square Prism
- Bernell Bull's Eye Chain
- Translucent Eye Patch
- Penlight
- Marsden Ball





Introducing THE NEW

S2

COMPUTERIZED BINOCULAR HOME EYE EXERCISE SYSTEM

Developed by Dr. Jeffrey Cooper MS, OD, FAAO & Rodney K. Bortel

Use HTS2 WEB-based Eye Exercise Program to Adjunct an In-Office Eye Exercise Protocol or Utilize HTS2 as a Stand Alone Home Eye Exercise Program



The HTS2 Program is Fully Automated and Easy to Use

Exercises include:

- Pursuits
- Saccades
- Base-In/Base-Out Vergence
- Jump Ductions
- Jump Ductions Random
- Base-Up/Base-Down Vergence
- Accommodative Rock

The New HTS2 runs on any iPad, Android, Chromebook, Windows 7-10, Mac OSx, and any other device/tablet that is 7.5 inches in width or greater.

HTS

HTS2 maintains a toll-free tech support line to assist your patients with any computer related issues they may encounter. The ease of operation, proven efficacy, automatic monitoring, and affordable price all allow you to immediately offer an extremely profitable new service to your patients.

BVA Binocular Vision Assessment*

Helps you identify those patients that will benefit from the HTS2 Program. Measures Horizontal and Vertical Phorias, Fusional Ranges, Accommodative Facility, Worth 4 dot, Pursuits and Saccades, plus an Asthenopia Survey.

* Windows only

HTSBSP

Contact Bernell for pricing or go online at www.bernell.com to order.

Your first HTS2 order will include a complimentary doctor's kit containing everything you need to immediately start offering HTS2 to your patients. ONLY ORDERS FROM VERIFIED DOCTOR'S ACCOUNTS WILL BE HONORED.

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com


*ONLY ORDERS FROM VERIFIED DOCTOR'S ACCOUNTS WILL BE HONORED. ORDERS WILL DROPSHIP FROM THE MANUFACTURER.

79

LOW VISION


The New HTS2 WEB-based Program Allows You to Manage Your Patients' Program via The Internet





Proven Procedures

HTS2 is a web-based eye exercise program using proven methods of behavior modification in an operant conditioning paradigm. The program uses random dot stereograms and other objective targets. HTS2 "knows" if your patient has performed the prescribed exercise session correctly.



Patients who cannot appreciate an RDS may be assigned a host of Flat Fusion targets. The patient may select the target they wish to use for each eye exercise session.

		HT	\$2		
		- ABAR	ette 1		
	culler :	1000	alle.	Sugartanet Innes	-
Dense	Petiant	8 Bectwertal Ne	100		
Lain	Bully	5 New Section New			2
Courtney	Hawkins	8 Dear Drawlant Ver			
Duana	regene	2 Bear Street and the		0	2
Notes	Separ	ii Rey Disalat Ma	-		
Amatia	Handhall	8 Days Briterium Ver			
Insika	Moone	A Desc Dros Carl Uni		0	
Candace	Cierk	A Sectional As		0	1
Michael	Lynch	8 Peer Transiant the	- 100		
		.0		-00	Ð.

Instant Review

Simply log onto HTS2 to see a list of all of your patients. You may view their results, check compliance, modify their program, or send notes right from your computer.



Dispense the HTS2 program in the Auto Mode and the program will automatically advance patients through successive levels of exercises as goals are achieved. The Manual Mode allows you to assign your own custom exercise protocol.



Provides Motivation

HTS2 presents your patients with bar graphs showing exercise goals. Performances from their last and current sessions are displayed. HTS2 virtually "pats the patient on the back" encouraging them to try harder.



Easy Monitoring

The HTS2 program provides graphs depicting each session's performance. The program displays date, time, duration and performance of each exercise session. HTS2 denotes when your patient has performed the exercise session and also if it was performed correctly.

NOTE: HTS2 DOES NOT TREAT OR DIAGNOSE ANY SPECIFIC MEDICAL CONDITION OR DISEASE.

Contact Bernell for pricing or go online at www.bernell.com to order.

*Only orders from verified doctor's accounts will be honored. Orders will dropship from the manufacturer.

ASSESSMENTS

OPTICAL

Windows 7 through 10 Compatible

PTS II is a home-based doctor dispensed computerized perceptual processing program, that has been designed to address a variety of visual perceptual/information processing domains.

These include simultaneous processing, sequential processing, speed-of-information processing, visual temporal processing and rapid automatized naming.

THE PTS II PROGRAM... DESIGNED SPECIFICALLY TO ADDRESS THE FOLLOWING PROBLEMS:

- Ordinary reading disability
- Word decoding difficulty
- Problems in reading comprehension
- Spelling problems
- Memory disorders
- Laterality-Directionality difficulties
- Slow Speed of information processing

PTSII Screening Program

Measures and quantifies patient's abilities in the following areas:

Simultaneous, Sequential and Speed of Processing. Results displayed at the completion of testing include Raw Score,

Non-verbal learning difficulties

Percentile, Scale Score and Descriptors

Mathematics difficulty

COMPUTERIZED PERCEPTUAL PROCESSING SYSTEM

Developed by Sidney Groffman, O.D., F.C.O.V.D.

COMPUTERIZED PERCEPTUAL PROCESSING

PTS2



Procedures

Eleven perceptual processing procedures that are highly motivating, interesting and successful techniques for all ages.

	THE REPORT OF	Contract of the local division of the local	Long I	
Villagenergy	5	1.1	00	
	1	1.	39	
		100	30	
faire they bearing		F TE	28	
		1	30	
Million Malarsh			30	
Statistics biographic	100	100	251	
Participal Stimulation	100	100	30	
	1.0	100	39	
	100		100	
			199	
Read			Annual	

Complete Control

You may select Auto or Manual mode operation, Starting and Ending Procedure Levels plus Session Duration



Easy Monitoring The Program Records Each Session's

Data for Your Review

PTSSP

NOTE: PTSII COMPUTERIZED PERCEPTUAL PROCESSING SYSTEM DOES NOT TREAT OR DIAGNOSE ANY SPECIFIC MEDICAL CONDITION OR DISEASE.

Contact Bernell for pricing or go online at www.bernell.com to order.

*Your first PTS II order will include a complimentary doctor's kit containing everything you need to immediately start offering PTS II to your patients.

*Only orders from verified doctor's accounts will be honored. Orders will dropship from the manufacturer.

REFERENCE

from HTS SANET VISION INTEGRATOR Developed by Robert B. Sanet, OD, FCOVD & Rodney K. Bortel

Using a 50" HD TV and interactive touch screen, The Sanet Vision Integrator (SVI) is a multi-sensory system incorporating eye hand coordination with cognitive processing and balance to help improve these abilities for a wide range of visually related learning problems.


The key features that make this "all in one" instrument so useful are the variety of programs and many ways that visual abilities can be enhanced. The programmable instrument incorporates features of a saccadic trainer, virtual rotator, tachistoscope and programmable metronome. The SVI instrument actually "speaks", instructing the patient to respond to verbal commands, improving auditory-visual integration and memory.

SVI can be used to enhance the following visual abilities: Pursuits, Saccades, Fixation Stability, Eye-hand coordination, Visual reaction time, Speed and span of recognition, Automaticity, and Contrast sensitivity, plus Visual and Auditory Sequencing and Memory.





Tachistoscope Stimuli includes numbers, letters and words.




Rotator Two rotator sizes, variable speed and direction, with a large range of stimuli parameters.



Five Categories Each category features several specific procedures.



Selectable Parameters Allows you to quickly modify the visual demand of any procedure.



Numbers/Letters/Word Saccades Helps build skills in saccadic accuracy and helps develop automaticity.



Programmable Metronome Provides immediate feedback using both auditory and visual stimuli.

NOTE: SANET VISION INTEGRATOR DOES NOT TREAT OR DIAGNOSE ANY SPECIFIC MEDICAL CONDITION OR DISEASE.



View our Video Demonstrations at www.svivision.com

The VTS 4 - HoloLens Computer Orthoptics Holographic Vergence Exercise System

Exercises Include:

- Smooth Vergence
- Rotations
- Jump Ductions
- Multiple Choice Vergence
- Pursuits
- Saccades
- Arcade
- Brock String

Measurements of the following:

- Phoria
- Fusional Ranges
- Pursuits
- Saccades
- Worth 4 Dot
- Motor Field
- Fixation Disparity
- Cyclo
- Aniseikonia

VTS4HOLO

The HoloLens Augmented Reality Technology Provides Detailed Holographic Vergence & Stereo Targets



Simultaneous Perception

Offers various vergence exercises to help

coordinate the eyes so that they may better work together as a team.

VTS 4-HoloLens Computer Orthoptics[™]

NOTE: THE VTS 4 - HOLOLENS COMPUTER ORTHOPTICS SYSTEM DOES NOT TREAT OR DIAGNOSE ANY SPECIFIC MEDICAL CONDITION OR DISEASE.



Flat Fusion



Stereopsis





REFERENCE



tem consists of 12 activities designed to make visual and sensory training fun and immersive, while providing effective treatment and metrics. Each game allows customization for contrast and brightness for suppression, convergence and divergence therapy, and monocular fixation in binocular field. Our VR activities are designed to work on eye-hand coordination, reaction time, visual memory, visual discrimination, peripheral awareness, sports vision, pursuits, depth perception, vergences, and measure eye alignment.

VR

remium

Oculus Rift S (Virtual Reality Headset and Tracked Controllers)

Software Custom Settings Allow:

ptics Trainer Starter Package (V2)

- Selective Eye Targeting
- Vergence and Torsion · Amblyopia settings to control Contrast, Brightness, and Blur for each eye

12 activities designed to make visual and sensory training fun and immersive!

OPTICS TRAINER SOFTWARE

Optics Trainer	VR Basic	P
Vision Therapy Settings (MFBF, Vergence, and more)	\checkmark	
Patient Prescriptions and Monitoring	\checkmark	
Stereo 3D Displays	\checkmark	L
6-Degrees of Freedom for Movement and Interactions	×	
7 Basic VR Activities	\checkmark	L
12 Premium VR Activities	×	

Dotics Trainer VR Basic

Optics Trainer Virtual Reality Basic is designed for home use and complements our Optics Trainer VR Premium Suite for the office. This software is specifically for headsets on the market today including:

- Oculus Rift, Rift S, Go, Quest, Quest 2
- · HTC Vive, Vive Pro, Cosmos
- Valve Index
- · HP Reverb and Windows Mixed **Reality Headsets**

This is a stand-alone headset that does not require any additional hardware. With the same tools as our clinical systems for amblyopia (lazy eye), strabismus (crossed eyes), double vision, and convergence insufficiency (CI).

Get the most sophisticated VR treatment system today!

Get the Home Therapy Software for your clinic and enjoy 7 of our fun activities!

OTVRH1 (1 Month Subscription)

*Software modules require a doctor's prescription for patients to purchase directly from Bernell. Hardware not included with a subscription purchase and is not available from Bernell, but retails for \$199

SPECIAL OFFER FROM BERNELL: TRY VR BASIC FOR 60 DAYS FOR FREE

Contact our Customer Service Team today to learn more details about this offer!

84

LOW VISION

Revolutionary Eye-Tracking Technology



FDA-cleared, objective eye-tracking biomarker that improves lives & performance.



Light-weight testing kiosk with research-grade eye tracking



Easy-to-understand patient assessment and test results

	Portel Countrees Countrees Indeed
37) (1) (1) (1) (1) (1) (1) (1) (1) (1) (1	X view rokens were locking were roken. Describe an eventscher here
	w () Kannan Ankan Maran Mara Maran Maran Mar
	A Marana and A Marana Angela and A Marana an



LOW VISION

EXAM ROOM NEEDS

LOW VISION


o visit

86

Objective • Repeatable • Quantifiable • Automated

Revolutionary Eye-Tracking Technology that Tracks, Measures & Analyses Dynamic Vision

- FDA-cleared medical device and break-through device designation
- Measure, record, analyze, and report minute eye movement patterns at a level of detail not feasible through observation-only exams.
- Uncover hard-to-identify eye movement impairments in real-time
- Transform the traditional 'follow-my-finger' test into digitally tracked ocular movements
- More than 20 digital, automated eye tracking tests
- Create patient-friendly reports with eye movement illustrations and video playbacks
- 5 different software modules: Functional Vision, Sports Vision, Reading Skills, Sensorimotor Exam, and complete customizable tests.
- Seamless practice integration!





Non-invasive, fast, and easy-to-use

Over 4 Million Patient Tests Taken!

PACKAGES: 1 YEAR, 3 YEAR, 5 YEAR OR LIFETIME SUBSCRIPTION

Unlimited Assessment Package

REYUAP

Dynamic Vision Package

REYDVP

Reading Skills Package

REYRSP

Sport Vision Package
REYSPVP



Scan QR code to visit

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

SENAPTEC SPORTS VISION EQUIPMENT

SENAPTEC SWIFT TOUCH

A modern light board designed to work on eye hand coordination, center/peripheral vision coordination, go/no-go, split attention, visual memory and reaction time.







SWT-721

INCLUDED WITH PURCHASE

- 4x4 ft. durable carbon fiber board with Kevlar finish (35 lbs.)
- · Android mobile controller device with customizable four training modules
- Standard wall plug
- · Static wall mount
- · One year manufacture's defect warranty

ADDITIONAL MOUNT OPTIONS AND ACCESSORIES AVAILABLE

Call or visit www.bernell.com for more information.





PRODUCT DETAILS

The Swift Touch has 72 LED buttons arranged in an offset concentric pattern on a 4x4 ft. carbon fiber board with a Kevlar finish. Large enough to engage the full range of upper extremity motion and the peripheral field of view. The button color, geographic frequency, response time, and feedback are all adjustable settings. The buttons can be pressed with the fingers or the palms, but do not respond simply by swiping across them.

The unit runs a cloud-based application that stores user data and can show trends/ comparison over time. The system is controlled by an on-board smartphone (included). There is no need to attach peripheral computers in this complete solution. Training modules includes; Eye Hand Coordination, Go/No Go, Split Attention, and Eye Hand Coordination w/ Decay.

Designed with mobility in mind, it weighs less than 40lb and the board folds up for easy transportation. Mounting options include a portable tripod, wall mount (static and motorized) and a motorized pedestal.

SENAPTEC SYNCHRONY

SAT-201 (Standard)

SAT-502-PRO (Pro)

The Senaptec Synchrony simulates motion to represent various paths and/ or objects in motion with LED lights. The mobile app can activate the series of LED lights consecutively to simulate movement for testing and/or training of anticipation or coincidence timing. The app can also activate the LED lights in a non-consecutive manner to provide stimuli for testing and/or training for dynamic vision, accommodation, and vergence. Users can respond with the included button to get results in milliseconds. Additional accessory options available online.



FEATURES

Curved liquid crystal lenses for a full 180° field of view. Most opaque flickering lens available - great for outdoor use Easy to read OLED display to show level, mode, and battery life Easy to operate with 2 buttons to adjust modes and difficulty level Designed to keep moisture out Elastic strap to keep eyewear tight on head during movement Rechargeable integrated battery via USB Free mobile app to control and customize the Senaptec Strobes Example drills included in the

mobile app for various sports

SENAPTEC STROBE Sensory Performance Training Eyewear Enhance an athlete's visual skills in the training room or on their field of play!

The Senaptec Strobe is designed to train the connections between an individual's eyes, brain, and body. Using liquid crystal technology, the lenses flicker between clear and opaque, removing visual information and forcing the individual to process more efficiently. The Senaptec Strobe can be integrated into existing sports training drills and exercises, or to vision therapy protocols as an uploading technique.



This app provides additional value and functionality for your strobes at no additional cost!

See pg 86 for info on the using strobes for amblyopia therapy.



SENAPTEC QUAD STROBE

The Quad Strobe offers all the capabilities of the regular Strobe, though with more precision, as each lens is divided into four flickering quadrants. From the companion app or directly on the eyewear, the user can activate each quadrant independently and customize the occluding (vision blocking) patterns to best address their specific training or medical rehabilitation needs.



Includes: Case, charger, and instruction manual



EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

BTRACKS[™] Assess Balance System A Computerized Balance Assessment and Training System





Balance and Posture are important indicators of health and well-being!

The BTrackS Assess Balance System gives health professionals an all-in-one solution to assess and train balance by administering computerized protocols with gold-standard force plate technology. BTrackS is used for balance assessment, fall risk analysis, balance improvement programs, concussion management, Medicare wellness, chiropractic care, physical therapy, vision therapy and balance research. Protocols are easy to administer and results include comparisons to normative data where possible. Since being introduced in 2014, BTrackS has set a new standard for accurate, reliable and affordable balance testing and training. The Assess Balance System is lightweight and portable and does not require AC Power or Internet Connectivity to operate.

BTrackS[™] has set a new standard for accurate, reliable and affordable balance assessment and training.



✓ Accurate ✓ Computeriæ d 🛛 ✓ Affordable ✓ Portable 📹 Objective

BTRACKS[™] BALANCE PLATE + BTRACKS[™] Assess BALANCE Advanced Software with Foam Pad

The BTrackS™ Assess Balance Advanced System includes twelve protocols (i.e. seven testing and five training) that provide the most complete solution for identifying balance issues and training balance. All testing protocols include The system comes with a foam pad for vision and vestibular testing. The BTrackS[™] Balance Plate includes a two-year warranty. The testing and training protocols include:


- Balance & Fall Risk (Testing)
- Weight Distribution (Testing)
- CTSIB (Testing)
- Limits of Stability (Testing)
- Cervical Challenge (Testing)
- Single Leg Stance (Testing)
- Custom Design (Testing)
- Random Targets (Training) • Target Tracking (Training)

Left/Right Targets (Training)

• Diagonal Targets (Training)

• Front/Back Targets (Training)

BTRACKS




Reporting and Tracking

Balance Training and Rehabilitation

BTrackS Assess Balance comes with multiple training protocols that allow focused training of postural control. This includes Left/Right, Front/Back, Diagonal and Random Target acquisition protocols, as well as a Target Tracking protocol. BTrackS Balance Training helps individuals learn to control postural sway using biofeedback. During training, the real-time location of a person's center of pressure is represented on an image of the BTrackS Balance Plate. The individual must move their center of pressure into targets of varying difficulty.

Each of the five training programs run for default of three minutes but can be longer or shorter as desired. Left/Right and Front/Back deliver uniaxial targets in various sizes and locations. Diagonal delivers biaxial targets in various sizes and locations. Random delivers multiaxial targets that appear randomly in various sizes while target tracking delivers a circular target that is moving around the screen and center of pressure should reside inside the circle as much as possible.

One of the most important features of Assess Balance is the tracking and reporting of results. Each time an individual is tested, the results are stored and available to compare against prior results. A concise and detailed report is always available to share with the individual, their family and their health providers. These same reports can be used for billing and if needed, for personal injury validation. All the reports are stored locally on the computer and can also be stored in local or cloud based servers.



BALANCE BOARD (VTE) WITH 4 SUPPORT BASES

Board Size: 49 x 49 cm.

- 30; 40; 60; 80 mm, diameter

Cylinders Sizes:

- 40 mm, thickness.


Balance board comes with 4 cylindrical supports, easily to interchange and disassemble. Board has an anti-slip surface; handle is cut in the board.

Balance board is made of aluminum material. supports are made of metal material.



Aluminum Balance Board

Aluminum balance board with hand holds makes this easy to use and easy to store. Comes with four round fulcrums of differing diameter to allow different levels of support. Weather stripping on one side makes a non-slip surface. The other side is a raised tread patterned, non-slip surface. Keep your two favorite fulcrums on and just turn it over to use either side. Painted gray powder coat. Fulcrums screw by hand to change.



LOMAN BEAM & WALKING RAIL FROM VTE

The Lowman Beam is a blue base made of two wings of 10 degree slope separated by a central divide. The central Walking Rail is a 4" wide unit which snaps onto the central divide. Each unit is 19.7" long and interlocks with other pieces to form as long of a beam or rail as you desire. Made of a strong plastic with a skid-resistant surface. The light weight and easy assembly allows it to be stored or moved.

E/LBWR1 (Single Unit - 19.7" Long)

E/LBWR2 (Pkg. of 6 - 9.84" Long Total

Dizzy Disc THERAPYWIDE 24"

Dizzy Disc Therapy Wide is geared for older children and or adult populations that provides the same benefits of sensory input and or balance/vestibular.

- · 24" wide platform
- · Built in handles
- · Smooth non-slip surface pad

DDTW24W



TURNUNIT 16" 3-D ROTATIONAL BALANCE TRAINER WITH VARIABLE TILT

The TurnUnit builds incredible balance and the ability to make controlled rotational movements up to a 3 dimensional plane of movement. It has a 16" Top Platform which offers a progressive tilting function with a built in lever to change the angle. The dynamic axis of rotation makes it far more

interesting than traditional flat rotational platforms. The Turn Unit is designed for people that are serious in not only their workout regime, but more importantly their results!

DDTU16B

TACTILE WALKING BEAM - SET OF 8

TWBS8

This set includes 8 interlocking walking beam sections designed to provide sensory stimulation, thanks to the hundreds of smooth, contoured bumps, molded into its top surface. Size per section is approximately 19" x 5" x 3" in size, which provides over 12ft of walking distance!

· Colors may of individual walking beam sections may vary.

Set can be easily cleaned with

- light soap and water.
- Maximum Weight: 176 lbs



Maze Balance Board

MAZEBB

The Maze Balance Board adds excitement to balance training! This balance board includes a special inner track that has a ball. Patients can improve balance, coordination, and eye-hand coordination by controlling the movement of the ball through the Maze Balance Board. The challenge is to move a ball around the top of the balance board. This Balance Board is an intermediate level board. It helps to develop balance and coordination at a higher level compared to a circular board.

Includes: 1 board and 2 balls. Colors may vary. Measures approximately 21" long x 16" wide x 1" high. Lightweight and easy to move, carry or assemble. Weight Limit: 176 lbs.





BALANCE/REACTION/SPORTS VISION

STANCE CHECK IZONE VISION TRAINING AID

These specially designed goggles were created to train athletes in golf, baseball and softball, but may be applied to other sports, as well as vision therapy exercises. The goggles utilize a specially designed nose bridge which acts as a septum, forcing the athlete to reposition their stance so that both eyes are trained on the target. The Stance Check iZONE Training Aid trains the athlete to realign



Includes: Goggles and 2 Eye Shields

their stance so that both eyes stay on target, which helps lead to improvements in performance by training proper positioning of the eyes in the stance so that binocular function is not restricted.


HECOSTIX Revolutionizing Hand Eye Coordination Training Red-White-Blue or Red-Green-Blue

WHAT ARE HECOSTIX? We've taken a simple game of playing catch and added a twist! Each leg is colored in a fashion that allows participants to call a color or hand combination and then throw the HECOstix to the receiver. The person receiving must process the information while the HECOstix are in the air and catch the corresponding color or hand combination called. Seems easy right?



PERIPHERAL AWARENESS CHARTS LETTERS AND FIGURES

Peripheral awareness chart develops the ability to maintain fixation over time. It also improves peripheral awareness and helps to avoid peripheral vision restrictions and to reduce visual stress. Made of cardstock. **Size:** 12" x 12"



VECTOR® BALL - COGNITIVE VISION TRAINING TOOL

Vector® Ball is the most dynamic and versatile Cognitive Vision Training (CVT) device available. CVT are drills that exercise the eyes and brain to improve an athlete's Speed of Reaction.

The Vector® Ball emits a short burst of light, randomly changing the color of the ball (Red, Green, or Blue) each time it is struck, bounced, caught, etc. This provides athletes with unpredictable information that enables, live action, CVT drills.



Rainbow Bean Bags

Bean bags are a great addition to every therapy room. This set of six durable,

canvas square bean bags comes in a rainbow of colors. Great for developing gross motor skills and practicing eye/ hand coordination. Includes suggested activities for play.

RBBS6



PITCH RECOGNITION BALLS

How better to train "Keep your eye on the ball" than with different color and shape marked "soft" baseballs? Comes with six printed balls and manual or make up your own games--bunt if it is a square; swing if it is green; don't if it is red.

(Printed - Not stickers. Wont rub off.)


CBALL



TRACE THE 8s[™]

TR8ST

Tracethe8s™ is a fun, engaging, therapeutic activity designed to enhance: Visual-motor skills, Visual-perception skills, Attention and focus, Working memory, Laterality, Oculomotor skills, Scanning skills, Figure ground skills, Eye-hand coordination, Sequencing skills, Speech, and Brain integration.



REFERENCE LOW VISION

SPACE FIXATION CLING SLIDES & CLING TRANAGLYPHS For Sports Paddles™ & Bernell Space Fixator™ Our most popular tranaglyphs TM are printed on cling material in full tone inks used when cancelling is not as important as quick convergence or divergence training such as sports training. Also, two of our tranaglyphs are printed in half tone inks for use when cancelling must be maximized (HFIX909 & HFIX520). Check out the New Peripheral Awareness Cling Letter Chart (PER3A)! HFIXEC1 (6" x 5.5" HFIX601HD (10" x 5.5" HFIX52HD (10" x 5. Peripheral Awareness Cling Letter Chart Ε ALL CLINGS ARE PRINTED ON B Μ TRANSPARENT н MATERIAL, WHICH ADHERES AND Ρ G 7 HFIX520 (6" x 5.5" **REMOVES EASILY** FROM GLASS OR ACRYLIC SURFACES С S BY ADHESION FORCES. Ο BSFIXD (23" x 19") PER3A (12" x 12" HFIX520HD (6" x 5.5" SPORTS VISION FIXATION PADDLES Uses targets with cling material, which adheres and removes easily from the front or back of the paddles by adhesion forces. The In-Office unit (HFIXIO) is a two-piece model that comes with comfortable plastic handle. The Home unit (HFIXH) is one piece, clear polycarbonate. Size: 12" x 12" Size: 12.75" x 12" + (Not Including Handle) (Not Including Handle) **IN-OFFICE UNIT + 7** HFIX909 (9.375" x 10") ↑ HOME Clings & Dither Pattern SHOWN WITH N-OFFICE SHOWN WITH UNIT HEIXEC1 HFIX909 UNIT Save Over 15% with this kit! HFIXIOSET (Kit) HFIXIO (In-Office) HFIXH (Home) ernell Space Fixator SPACE FIXATOR **IMPROVED SWIVEL!** The unit now has the ability to tighten the U Channel if the BY VTE ITALY acrylic target ever becomes loose. Also, new bushing allows for effortless 360° swivel. It is used to train effi-Brass — Bushing For Effortless 360° Swivel The unit quickly adjusts for any height of child or adult. The rectangular acrylic panel easily ↑ ciency in saccadic rotates showing the too hard of a touch by the user. Bernell sells cling tranaglyphs that can BSFIXD on movements and much be used with the fixator. Cling material can be cut with scissors and partial targets used one BSFIX2 more. It integrates eye or more lines at a time as well as mixed as shown to right. You can set the unit in the middle movements with spatial of a room and do training while looking through it at a distance chart or even a Marsden Ball. awareness as well as visually guided move-KIT INCLUDES · These charts can be cut with scissors to limit clutter. ments in a coordinated · Some lines are Red/Green sensitive. & efficient manner. (8) CLINGS + FIXATOR All cling slides are · Clings can be reused or rejoined on the fixator. Rounded plastic with pictured above. • BSFIXD (Dots) · Stick on as many lines as you want at a time. permanently attached HFIX52HD (BC52 FULL INK) Arrows can be turned 90° to change the pattern. fixation dots. Has flat HFIX520 (BC520 HALFTONE) metal base and rotates in the "Z" axis easily. HFIX520HD (BC520 FULL INK) BERNELL SPACE FIXATOR KIT Improved design Adjusts easily for height HFIX601HD (BC601 FULL INK) BSFIX2KIT (Size: 24" x 32") for greater stability. differences. HFIXEC1 (EC1 SLIDE) HFIX909 (TRAINER HALF TONE) BERNELL SPACE FIXATOR ONLY

PER3A (PERIPHERAL AWARENESS)

E/SFIX

LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 92

BSFIX2 (Size: 24" x 32")

EDS	ROTATION T	RAINERS			
EXAM ROOM NEEDS		Rotat	ion Trainer	 Handheld Digital Remote with Magnetic Motor Att Magnetic Quick Connect Disc Attachment One-Handed Height Adjustment Clutch Sleek Anodized Aluminum Motor Enclosure Speeds of 3 to 50 Rotations per Minute 	achment
ROO		Floor &	Table Model Kits		
EXAM		Manageria E Riskers	perceptual and space awareness, dyna	fessionally designed and engineered to test and enhance eye/hau amic visual acuity and dynamic fusional training. This new rotator ha ors. The motor and handheld controller have been redesigned from t	as a multitude of
	Bicrosoft Economics		connects via a standard ethernet cord. T increase and decrease the speed by ho	k rectangular anodized blue aluminum housing. The handheld This new rotator can get down to 3 rotations per minute. The cont olding down the speed + or - buttons. It also features a clockwise or or desired rotation direction button and the disc rotation will reverse. The	troller can rapidly counterclockwise
ASSESSMENTS			started, stopped and turned on and controller's digital screen. The con	nd off via the corresponding buttons on the controller. All information is ntroller attaches to the motor via magnetic connection when not in use	displayed on the e.
WSSE	Massarine Expensional Africanser	Hackerry Gury Conner Service for the conner	No more unscrewi attachment adapter	w attached to the hub via Bernell's new quick connect ma ing discs from the hub when swapping discs. The discs now hir r that snaps into the magnetic hub. Additional front magnetic attach for other discs for the VTP Edition Rotation Trainer (VTP-HUBFR) f	ave a magnetic ments with knob
ASSE	e à) =- 	over when using mu The floor stand fe	ultiple discs. Plastic discs are 201/2" in diameter. Pegboards inclu eatures a one handed clutch adjustment. This allow the pole he	de tees. ight to be easily
	VTP EDITION FLOOR MODEL	VTP EDITION TABLE MODEL	COMPLETE SY cord, digital cont	and downwards to the desired height. Approximately 3½ to 5½ fee STEM INCLUDES: Stand with heavy-duty base, moto troller with ethernet connection cord, deluxe multi-position	r with power
OPTICAL	DISC OPTIONS:			ment hub and disc adapter, and disc.	
.dO		Туре Рак VTP ЕDITION FLOOR К/т VTP-106	T # PRICE PART #	PRICE	
≻	Disc	VTP EDITION TABLE KIT VTP-100 STANDARD FLOOR KIT BS1060			DRTS DISC
RAP	Pegboard Discs	STANDARD TABLE KIT BS1060 DISC ONLY BS106			
VISION THERAPY	Acrylic Peg			STANDARD ROTATION TRA	INERS
SION			Professionally designed	ABLE MODELS d and engineered to test and enhance eye/hand coordina ness, dynamic visual acuity and dynamic fusional trainir	
5		Floor Model	system features electro ible, variable speed (5	5-32 RPM), heavy duty motor drive, deluxe multi-positio d heavy-duty base. Floor Stand Height: Approximately 3'	only), revers- mal tilt head,
NO	Disco Dial	One-Handed Height Adjustment Clutch		in diameter. Pegboards include tees.	/2 10 0 /2 1001.
ITAT	Stpanded Vinyl Dis		Floor 8		Table Model eight Adjustment
IABIL			Tilt	Rings	
NEURO REHABILITATIOI	Disc Thickness: 0.125"			The following items come standard in each kit: c, (1) Pole Tilt, (1) Control Box, (1) Rotation Motor, (1) Hub for Rotator Stand/Base, (1) Control Box Mounting Rings (Floor	
EURG	PLASTIC GOLF	- Tees 🧼 🗍 Tai	BLE TOP MANUAL F		
z	Red & Green (Pro		PEGBO		
NO	22	Size: • Length: 68mm	L95TRP		
LOW VISION		• <u>Diameter</u> ; 3.8mm	1		
LO	Wooden			A Cho	
щ		angle (al table top rotator. Rotates either direc (5 settings from 90 - 0 degrees) table	top rotator. Non- to last. Has large ball bearings	s to easily
REFERENCE	BRIGHT COLORS (PKG. OF 60)	various	zed so the rotation is under therapist co s colored shapes that attach to the rota on of desired placement for any 0.19 sh	ator for easy iden- Comes with mix of bright colored	
FEF	BCTEES		MANUAL ROTATORS ARE GREA	AT FOR WORKING WITH STROKE OR TBI PATIEN	TS!
ш	93 Fax Orders to	Save Time & Woney at the	74) 259-2102 or 2103 or (Go Io www.bernell.com	

ROTATION TRAINERS

Vinyl Attachment Training Discs and **Clear Mounting Disc Adapter**

Complete Set of 4 Vinyl Attachment Discs with Clear Mounting Adapter Disc BS1060DASET

Designed for use with the Bernell Rotation Trainers (both the VTP Edition and Standard Rotator). We now offer 4 front/back vinyl discs. Dynamic Acuity, Anti-Suppression and two different versions of Simulated Depth (aka Wallach Ring) Discs.



The discs are 20-1/2" in diameter and made CLEAR MOUNTING ADAPTER DISC of washable vinyl. They are made to fit on any of the rigid acrylic or expanded vinyl discs that Bernell sells (Part #'s BS160AD, BS1060GD, BC1060RD, BC1060ED and BC1060SD). They also fit our Clear Mounting Adapter Disc (Part # BS1060CMDA1) which allows for quick change out of the discs.

> The clear mounting disc adapter is made out of clear acrylic and has 3 adjustable knobs to secure the disc. The adapter works with discs sized 20.5" down to 18.5". The URI Training Discs (Part # JWRDS) will work on this adapter as well. Add on an Individual Front Magnetic Attach for the VTP Edition Rotation Trainer to the clear mounting disc adapter (Item #: VTP-HUBFR) for easy on/off attachment or attach the clear disc to the Standard Rotator Attachment Hub (Item #: BSHUB) for use with the Standard Rotation Trainer.

> > DYNAMIC ACUITY This disc features Snellen letters

and numbers with sports balls on

one side and animal characters with

numbers and letters on the other.

The sports side was designed for use

in dynamic acuity testing and training,

which makes it a valuable addition to

any sports vision practice for use with

athletes. On the opposite side are

colorful and engaging animal images



OPTICAL

marine characters on the other. These discs work with Bernell red/ green goggles (not included) so that the red images are seen with the red filter and the green images VISION THERAPY

and are an engaging tool for anti-suppression training while also working to improve one's ability to accurately direct their eyes from one point to another point with ease and efficiency.

ANTI-SUPPRESSION

This disc features letters, numbers,

shapes, and lines on one side and

are seen with the green filter. White images act as binocular locks. The discs have excellent cancellation

This disc has two of the more

popular simulated depth visuals.

One side features the classic

Plateau spiral, which is used for

the Elliot Forrest myopia reduction

procedures. On the opposing side



FRONT/BACK

BCVDAD1

is often a good way to change the depth that is seen. It is best if the patient discovers this or other ways to shift the depth pattern.

If the patient finds this too difficult of a task, then you may want to shift to the opposite side which features a single continuous spiral. The depth is often more readily observable with this version of the Wallach rings.

SIMULATED DEPTH 1

that are designed for training eye tracking with pediatric patients.

The use of shapes, letters and numbers along with the animals

allows for numerous activities in the therapy room.

On one side of the disc is a special design known as a Wallach ring. It is broken up into three zones with the depth in each zone being the opposite of the adjacent zone. This is often a good place to start to see how well a person can organize visual information. Shifting attention from looking hard at the small circle to looking soft at the outermost region

generally more difficult to see the depth in the rings like the one at the 12:30 position.

Individual vinyl discs include four binding clips for attachment to any Bernell Training Disc that is 20.5" in diameter.

VTP Edition and Standard Bernell Rotator Extras and Replacements

	Product	PART #	Price
Α	CLEAR MOUNTING DISC ADAPTER	BS1060CMDA1	
В	ROLLING BASE FOR FLOOR ROTATOR STAND	BS1060RB1	\$184.95
С	DELUXE POLE TILT	BSRPT2	
D	FLOOR MODEL ROTATOR STAND	BS1060FS1	\$149.95
Е	TABLE MODEL ROTATOR STAND	BS1060TS1	
F	VTP Edition Rotation Trainer Standard Kit (Motor, Remote Control, & Hub)	VTP-1060SK	\$969.95
G	INDIVIDUAL FRONT MAGNETIC ATTACH WITH KNOB FOR VTP EDITION ROTATION TRAINER	VTP-HUBFR	
н	REPLACEMENT COMPLETE HUB WITH Magnetic Attachment for use with VTP Edition Rotator Only	VTP-HUB	\$42.95
Ι	HUB FOR STANDARD ROTATION TRAINER	BSHUB	



Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 94





FRONT/BACK

BCVASD1

BCSDR2

easily in the more tightly banded rings, such as the one seen in the 2 o'clock or the 6:30 position. It is

of the disc is a figure used for an activity called "High-Level Wallach Rings". Here eight different Wallach Rings are presented around the

considerably. Generally, patients will see the depth more

outside part of the disc. The thicknesses of the bands vary

MARSDEN & SACCADIC TRAINING



EPORT[®] II Vision Training System Clinically Proven, FDA Cleared, Doctor Recommended.



This EYEPORT® II is now battery powered for portability. The unit is powered via 4 AA batteries (not included). No more cords or chargers. The training device also features an updated tripod, new Deluxe Reversible goggles with special red/blue lenses (Elastic Version) and custom box.

The EYEPORT® II is the only vision training tool to be cleared by the FDA for the treatment of poor accommodative and vergence facility, convergence insufficiency, and large accommodative lag (in non-presbyopic subjects).

What makes the EYEPORT® II unique is its patented use of alternating red and blue lights which help patients to achieve dramatic results faster and more efficiently when training. The EYEPORT® II is a great tool for both in-office and at home vision training and is a gateway product for practitioners interested in adding Vision Therapy in their practice!




EYEPORT TRAINING SYSTEM INCLUDES:

- Eyeport II Battery Power Unit Requires 4 AA Batteries (Not Included)
- Deluxe Reversible Goggles with Special Red/Blue Lens Elastic Version (Includes Glasses Case)
- Versatile Tripod with vice functionality to allow for clamping to tables, doors etc.

Instruction Manual





EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

 Pre-programmed workouts Multiple speed settings

Protective storage case

Special glasses to train each eye individually

ALPHABET PENCILS

A-Z on one side and Z-A on the other. Great for all sorts of saccadic or pursuit training techniques and then as a reward for good work. Patients will take this to school and tell their friends where they got it.

LBPENCIL (Pkg. of 12)

Bernell Saccadic Sticks™

LOW VISION

Designed to be a multi-function training tool. These can be used for saccadic or pursuit training, physiological diplopia (as a type of Brock string with letters as targets holding the rod flat or on its angle), accommodative rock tool, training for "A" or "V" syndromes as well as sports vision training. These are a great addition to any therapy room. Instruction manual includes pictures of all 4 sides of each stick for quick reference.

Length: 25" (Including hook) Material: Clear Acrylic Weight: 8.4 oz. Each



BSTICKS

STICK TRAINER BY VTE

Round wood rods with black characters and figures on yellow backgrounds. For home training. Easy to use for saccades, pursuits, and accommodative rock.

> Rounded design allows for different targets: Random line letters on 1st side; Random line numbers and arrows on 2nd side; Line of figures for preschool children on 3rd side.


Available in: Package of 6 (3 Pair) or Package of 14 (7 Pair) E/PT E/PT14 (Pkg



95

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

REFERENCE



BERNELL SPORTS LETTER BALLS - BASEBALL SIZE

Baseball sized (9" circumference; 2.86" diameter) but **soft** like a softball. Available with four font sizes (12pt, 18pt, 24pt, & 30pt) or two font sizes (12pt & 24pt). As student/athlete improves, he/she can do similar tasks with smaller fonts. The letters are printed on the ball--not just stickers-- so they won't come loose & fall off. For Sports Vision Training - why not train with a ball that looks like a baseball instead of a play ball? It is softer than a real hard ball. Same font sizes are available with softball sized (12" circumference; 3.8" diameter). The balls come attached to 12ft string. Adhesive ceiling hooks are included, with replacements available (**MBHK**....\$1.25).





Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 96

LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

BROCK STRINGS



COMFORT GRIP HANDLE FOR BROCK STRINGS - PACK OF 3

Plastic handle for your brock string. The knot in your brock string slips over the groove in the handle and is held securely in place. The handle gives the patient a nice grip when training. Installs in less than 5 seconds.



GTVT BROCK STRINGS

Some doctors prefer the size and feel of these string units over the Bernell ones above. These have Red, Green, Yellow, Blue, and Orange beads.

PRODUCT	PART #	PRICE
12ft - Package of 12	FR12	
12ft - Single String	FR12ONE	
6ft - Package of 12	FR6	
6ft - Single String	FR6ONE	

Bernell Physiological-Diplopia Cord™

Inexpensive home use instrument for near training as well as to extend phy-dip training from near into distance. Utilizes physical diplopia for training suppression, binocularity and spatial localization. Heavy cord with colored beads.

(Available in 10ft (3 Beads) or 20ft (6 Beads) Lengths, Single or Dozen)



)	PRODUCT	PART #	PRICE			
I	10ft - Package of 12	BC1090				
I	10FT - SINGLE STRING	BC10901				
I	20FT - PACKAGE OF 12	BC1092				
J	20FT - SINGLE STRING	BC10921				

SPORTS BROCK STRING

Home use instrument for near training as well as to extend phy-dip training from near into distance. Utilizes physiological diplopia for training suppression, binocularity and spatial localization. Sports beads add additional level of focus when training.

Includes:

- Heavy Black Cord (10ft)
- 5 Sport Beads
- Red Handle for Grip

BSPDTC10 (SINGLE)



Phy-Dip Trainer[™]

Ideal home training device for non-strabismics, esotropes with near centration point, and patients with other binocular vision problems. Three movable plastic beads on hand-held copper rod allow extensions of physiological diplopia techniques from 2" to 20".

TNARSET

TONDEL'S ARROWS

BC1091

Grazyna M. Tondel, MD, MS, PhD

Consist of 6 pairs of colorful arrows. The arrow point from each pair touch each other along the middle, vertical line of the card. The card folds along this line and is shaped into a wedge. When used like a Brock String, the card enables training both smooth and jump vergence movements. Cards come pre-cut to easily fold into shape.

TNARSET5 (5 Cards)

SEE PAGE 83 FOR INFORMATION ON THE EYEPORT II VISION TRAINING SYSTEM Bernell C/L Lighted Phy-Dip String™

ADD SOME COLOR AND FLASH TO YOUR TRAINING!

The Phy-Dip or Brock string has long been a mainstay in vision therapy. It is a simple, but effective tool for teaching binocularity and anti-suppression. The CL Lighted Phy-Dip String™ is an improvement to the old stand-by. Kenneth J. Ciuffreda, O.D., Ph.D. and Diane Ludlum inspired the concept for the C/L Lighted Phy-Dip String™.

The advantage of the CL Lighted Phy-Dip[™] is the self-illuminating strings, which can flash at two different speeds or remain constant. The red/blue dual unit when used with red/green goggles allows the suppressing eye to be stimulated. The ability to light the string itself, not only adds interest and enhances the ability to see the string, but also helps break down suppression, making it more difficult to shut-off the non-dominate or lazy eye. The CL Lighted Phy-Dip™ is not intended to replace the original, but rather be a modern approach or technique.

CLSTRINGS (SINGLE - (1) White Lighted String with 5 Beads) CLSTRING (DUAL - (1) Red & (1) Blue Lighted Strings with 5 Beads Requires AA batteries (not included) (2) for Single and (4) for Dual string.

VISION THERAPY



For testing retinal correspondence between the two eyes, hold the device in a vertical position and flash it onto the macula of one eye with the other closed, then flash a horizontal image onto the macula of the other eye with the opposite eye closed. The patient can then receive visual feedback for proper retinal correspondence.

Includes: Handheld After Image Generator, After Image Enhancer Fan & Product Guide

Order Now by Phone at 1-800-348-2225 or Go To www.bernell.com

98

LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

All Bernell goggles come with temples or elastic strap. The goggles are available in packages of 6 or individually. *Bernell Red/Green & Red/Blue Glasses can be ordered without lenses glued to frame. Go online to order empty frames.


EXAM ROOM NEEDS

VISION THERAPY

NEURO REHABILITATION

LOW VISION

	Мси		RNELL MODE	
New New	FILTER TYPE		PART #	
		Dkg of 6		
	RED/GREEN	Single	BC11891	
Model	RED/GREEN			
-	WITH ELASTIC		BC1189E1	
Add NG to part numbers	POLARIZED	Pkg. of 6	BCPGNF6 BCPGNF1	
to order with lenses "Not Glued" to frames.		Single Pkg. of 6		
Not Glued to frames.	POLARIZED WITH ELASTIC	Single	BCPGNFE1	
		Pkg. of 6		
STIT.	RED/BLUE	Single	BC1189RB1	
1000	RED/BLUE		BC1189RBE	
	WITH ELASTIC	Single	BC1189RBE1	
	-			
~ ~	-	1		240
		6		7
			X	¢
	20	-		
			-	
		_	_	
ffers a similar size fr				
unded bridge and ter	nples for a	more co	omfortable	feel.
ur sturdiest model!				



nes offers a lightweight, attractive design similar in size to Model, but with rounded edges for a more comfortable feel.

, the filters are kept in place through a secure sonic weld. ed in ordering frames without the weld, which will allow easy the frame.

ORIGINAL	Original Bernell Model	Modern	Modern	Model	
	FILTER TYPE QUANTITY PART # PRICE		FILTER TYPE QUANTITY	PART # PI	RICE
MODEL	RED/GREEN Pkg. of 6 BC1180	MODEL	RED/GREEN Pkg. of 6	BC1185	
	RED/GREEN Single BCRG1 RED/GREEN Pkg. of 6 BC1181		Single	BC11851 BC1185E	
Hard temples can be bent to fit on	WITH ELASTIC Single BC1181E1	These goggles	WITH ELASTIC Single	BC1185E1	
and the second se	Polarized Single PCPC1	are larger than the	Pkg of 6	BC1185P	
patient's face.	Single BCPG1	original frames.	POLARIZED Single	BC1185P1	
Add NG to part numbers	POLARIZED Pkg. of 6 BCPGE		POLARIZED Pkg. of 6	BC1185PE	
to order with lenses "Not Glued" to frames.	WITH ELASTIC Single BCPGE1	Add NG to part numbers to order	Pkg of 6	BC1185PE1 BC1185RB	
	RED/BLUE Single BCRB1	with lenses "Not Glued" to frames.	RED/BLUE Single	BC1185RB1	
ATT.	RED/BLUE Pkg. of 6 BCRB6E	cluce to names.	RED/BLUE Pkg. of 6	BC1185RBE	
	WITH ELASTIC Single BCRB6E1		WITH ELASTIC Single	BC1185RBE1	
				6	
	e neversidle wodel		d Reversibl		ΞL
with larg	v design is a deluxe reversible frame e lenses and curved temples for a more ble and no slip feel.		goggles are reversible. The er than the old model & fold		v
- Contraction of the second	DELUXE REVERSIBLE BERNELL MODEL	March 1	CHILD REVERSIBLE	BERNELL MODE	EL
	FILTER TYPE QUANTITY PART # PRIC		FILTER TYPE QUANTITY	PART # PI	RICE
and the second se	RED/GREEN Pkg. of 6 BC1190		RED/GREEN Pkg. of 6	BC1186	
	Single BC11901	the second second	Single	BC11861	
	RED/GREEN Pkg. of 6 BC1190E WITH ELASTIC Single BC1190E1		RED/GREEN Pkg. of 6 WITH ELASTIC Single	BC1186E BC1186E1	
	Dkg of 6 BC1100DB		Pkg. of 6	BC1186RB	
	RED/BLUE Single BC1190RB1		RED/BLUE Single	BC1186RB1	
	RED/BLUE Pkg. of 6 BC1190RBE			BC1186RB6E	
*	WITH ELASTIC Single BC1190RBE1		WITH ELASTIC Single	BC1186RB6E1	




EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



LOW VISION



101 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

Letters kit includes all 3 charts & New Age Red/Green Goggle (BC11891)

Large Chart

Kit - All 3 + Deluxe Reversible Goggles

RGFL

RGFKIT

Figures kit includes all 3 charts & Deluxe Reversible Red/Green Goggle (BC11901)


EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



E/TVT (2 PAIR)

с (

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 102

*Cling pictured on iPad scre

BC65RG

LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

DPTICAL

VISION THERAPY

NEURO REHABILITATION



This system is similar to the Red Red Rock technique sold in the past by Mast-Keystone without the need for instrumentation. For use in the office as well as home VT. Useful for training anti-suppression using R/G goggles. Cards are made of same material as our Visual-Tactile Cards. Comes with a pair of New Age R/G goggles & an instruction manual. Provides therapies for accommodation, bilaterality, cognition, reversals & suppression.

LOW VISION

Includes: 12 laminated charts (front & back), R/G Glasses (New Age Model),


2 dry erase markers, eraser, and instruction manual.

coordination, Sequencing skills,

Speech, and Brain integration.

HEADLINE Anti-Suppression Charts

These chart packs (Bordered with Red) provide superior cancellation with Red/Green glasses. The Single-Color charts are fantastic for the ever important MONOCULAR IN A BINOCULAR FIELD training, and the Multi-Color charts provide exceptional suppression monitoring during binocular vision training procedures. All of the charts come on heavy duty 8-1/2" x 11" card stock.



BERNELL DUAL LED SLIDE ILLUMINATOR[™]- VERSION 3 (MVP3)

The unit holds two sets of vectograms, vectographs, or tranaglyphs. This includes BOTH Stereo Optical and Vision Assessment Company vectos! The unit utilizes LED lighting so that it can be used worldwide for thousands of hours. LED light intensity can now be dimmed via a rotational knob on the power cord. The cord is also now threaded through the unit for ease when adjusting tilt. We have redesigned this model to feature adjustable top and bottom bars to ensure any vecto or tranaglyph fits perfectly! Constructed of high rigidity plastic making it light weight, yet very durable.

The unit acts like a Ganzfeld so that the user has a minimal reference frame when viewing Tranaglyphs™, Vectograms, or Vectographs. This allows the target to better define the task.

HIDDEN CORD

HOLDS STEREO OPTICAL OR VAC VECTOS!

DUAL POLACHROME[™] ILLUMINATED HOLDER LED AND STANDARD VERSION

BACK SUPPORT FOR BETTER STABILITY

DIM SWITCH

Features individual units that easily snap together to create a dual unit, or can be used individually as hand-held slide holders. This is the perfect instrument for training with strabismics, non-strabismics or with pursuit and jump ductions at true distances. For in-office presentation of Tranaglyphs[™] or Polaroid Vectograms. Holds two pair of slides.

ADJUSTABLE BARS

LED Version are back-lit by three LED bars that are controlled by a dimmer switch. The Standard Version uses two 40-watt bulbs (included). Adjustable tilt legs. White ABS plastic construction. (Slides not included).

HORIZONTAL SLIDE HOLDER[™]


Shipping Note: This unit can be

shipped overseas! Works on 110/240V with plug adapters (Included).

Add this high quality, attractive

unit to your training room today!


Tough transparent plastic with diffuser plate, and doublesided white/ghost free card. Designed for use with vectograms or with all Bernell series tranaglyphs. Can be hand-held or used with legs. Size: 5.5" x 14". Fusion slides sold separately. Tranaglyph not included.





BC1003

STANDARD VERSION



BC1001 (1 Holder)

Side View

GHOST-FREE

BACKGROUND

CARD ONLY

DUAL HORIZONTAL SLIDE HOLDER WITH HANDLE

The lightweight Dual Horizontal Slide Holder is an economical tool that allows the patient to train in free space, shift accommodation through the paddles and back to the target. By turning the paddles around you can easily shift from base-in to base-out demand, as well as change the 3D target from coming towards you to going away. Perform jump duction training with two different base-in demands or two base-out demands

or have top base-in and bottom base-out and "jump" back and forth between the two images.

BSDSH1

lote: Cling targets, shown on page 80, can used on slide holder. The horizontal slide holders can be ed from the acrylig via their Velcro attachments.

The clear fixation paddle is 12" x 14" made of clear acrylic. It utilizes (2) horizontal slide holders that are attached to the paddle via Velcro® connections. These slide holders fit Stereo Optical Vectograms[®] as well as Bernell Tranaglyphs™

BP40T10IF

FREE SPACE SLIDE HOLDER & ILLUMINATOR



DFLS

ASTFSH

This free space slide holder is adjustable to accommodate slides of any size. Can be used with fixed or variable Tranaglyphs, Vectograms, Vectographs (not included). This slider holder may be used with the Dynamic Fusion Light Source (DFLS) or in free space across the room looking at the light from the DFLS. LED lights will last 20,000 hours. DFLS has universal transformer so it can be used outside US with an adapter which we will provide.

> Unit without Slide Holder DFLSW

LED VERSION

BC1003LED

LOW VISION

*Does not include pictured Tranaglyph and Vectogram

105

EXAM ROOM NEEDS VECTOGRAMS, VECTOGRAPHS & TRANAGLYPHS FIXED & VARIABLE VECTOGRAMS BY STEREO OPTICAL (EACH INCLUDE VIEWING GOGGLES) 6 CLOWN - VARIABLE SOV6 SPIRANGLE - VARIABLE SOV5 С Α в QUOIT - VARIABLE SOV2 FIGURE 8 - FIXED SOV3 COMPASS POINTS - FIXED SOV4 C 0 DX0 ASSESSMENTS **S**TEREO **O**PTICAL REPLACEMEN SLIDE HOLDER DADT # **BASIC FUSION - FIXED** SOV9 SOGUIDE STEREO TEST - VARIABLE SOV11 CHICAGO SKYLINE MOTHER GOOSE SPIRANGLE WITH BASEBALL- VARIABLI SOV7 SOV10 SOV14 SOV12 VARIABIE VARIARI I ACUITY SUPPRESSION - VARIABLE VECTOGRAN INFO: Polarized vectograms ensure testing and training at actual distances. Vectograms stabilize fusion and stereopsis, eliminate suppression and develop simultaneous vision. Used in mono and prevent anomalies involving projection and hand-eve coordination. Ideal for use with Bernell's Dual Polachrome Trainer and the Bernell Dual Slide Illuminator. Replacem cular, binocular, pursuit and saccadic training, vectograms help correct ent Polarized Viewers Available on pg. 100. **OPTICAL** VARIABLE POLARIZED VAC REPLACEMENT SLIDE HOLDER ECTOGRAPHS FROM VAC VA1050 VAC Vectograph Sets Include: Vecto Guides, Doctor & Patient Instruction Manuals, and Standard Polarized Viewers. 111 ADDITIONAL 00 VAC VECTOGRAPHS AVAILABLE ONLINE. VISION THERAPY VA1038 VA1039 CHESSBOARD VA1042 Zoo VA103 PANDA FARM VA1040 VA1041 PLANETS RACETRACK FRUIT BOWL VA1043 VA1044 VA1045 POOL BALLS VA1048 SPORTS BICYCLES MR. EGG VA1046 VA1047 PUMPKIN PATCH **NEURO REHABILITATION** POOL TABLE VA1049 DART BOARD VA1051 BALLOONS VA1052 BOOKSHELF VA1053 VA1054 VORTEX VA1055 AQUARIUM (B) (C) m E (E) (E) VA1056 TEDDY BEAR VA1057 VA1060 GEM W/ TARGET VA1060PL CHAIN Frog VA1058 POPCORN VA1059 GEN Vectograph KIT -0000 00 BERNELL NI 0000 000 0000 00 **LOW VISION** The new Bernell Mini Vectograph kit is an excellent device for training fusional **Bernell Mini** õ convergence and divergence ranges in true space. The mini-vectographs Vectograph Kit 0 6 60 Vertical orientation great cancellation and ability to evoke the perception of SILO (Smaller In/ BC275HV Larger Out), makes them a valuable asset in visual training. The BC275HV VA1022 VOLUME Bernell Mini Vectograph is used for training of lateral (eso and exo) phorias and training vertical phorias. Pair (Replacement) VAC Fixed The new clear holder has PD markings for both horizontal and vertical Kit Includes: Polarized Vectographs training. The Clear Holder в Mini-Vectograph Pair new sizing New Age Polarized Goggle of the miniā 0 (B) 0000 Instruction Manual 000 00 vectographs 00 0 80 0 õ 8 work in both 0 training orien-由 8 tations.

BC275VA

VA1020

VA1021


OPAQUE VOLUME

CLEAR VOLUME

UPDATED PRODUCT!

REFERENCE

106



REFERENCE

107 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

VECTOGRAMS, VECTOGRAPHS & TRANAGLYPHS


Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

Necklace 3D Variable Tranaglyph™

Outer Space Variable Tranaglyph™

not included)

Available as tranaglyph only or as part of a kit. Kit includes tranaglyph, horizontal slide holder & card (BC550), New Age red/green glasses & manual.




108

Vectograms, Vectographs & Tranaglyphs

KITS (50 SERIES & 60 SERIES) DEMANDS

Ideal for testing or training in true space at the office or home. Trains suppression while also developing fusion and stereopsis. Non-variable prismatic Red/Green slides build fusional convergence & divergence. 2PD increments from 0 to 30PD BI or BO. 10" x 5.5" Red/Green targets on frosted vinyl. Can be used with Polachrome Trainer.




109 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

BCB6-10 (B6 through B10 Slides)



ASSESSMENTS

OPTICAL



Some international flippers are priced lower then ours, but watch out! They will charge you international shipping and taxes. Save time and money and buy from Bernell!

111 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

RISM FLIPPERS

Professional training with two yoked prism pairs. Beginning with the highest powers through which a patient can fuse and focus, prism flippers are exchanged for the next higher powers through the training sequence. Great for early presbyopes with high exo with new bifocals! Please specify diopter power when ordering. Prism's base and power are marked. Prisms can be rotated, but are delivered as BI and BO.

Higher power flat lens flipper combinations are available, please contact us for more info.

FLAT LENS FLIPPERS

PART #	DIOPTERS	PRICE
BC1275+	0.25, 0.50, 0.75,1-7D	
BC1276+	8-10D,	
BC1276+	12D & 15D	

CORRECTED CURVED LENS FLIPPERS

Part #	DIOPTERS	PRICE
BC1275+CC	0.25, 0.50, 0.75, 1-7D	
BC1276+CC	8-10D	

Specialty Order Prism Flippers

Special Order Prism Flippers available in any four powers above. Pricing based on highest power. Please allow two weeks for delivery. Go online to download the Special Order Flipper Diagram.

co onine to dominoud the opecial order i hpper Diagram.						
PART #	DIOPTERS	PRICE				
BC1275S	Flat: 0.25, 0.50, 0.75, 1-7D					
BC1276S	FLAT: 8-10D					
BC1275SCC	CORRECTED CURVE:					
BC12755CC	0.25, 0.50, 0.75, 1-7D					
BC1276SCC	Corrected Curve: 8-10D					



Spherical lenses are 48mm in size and allow easy awareness of periphery. Ideal for perceiving SILO effect as well as other activities. Easy use for right or left handed persons.

Big enough for trial use of adds or to balance a prescription!

Available Powers:

+ 0.25, 0.50, 0.75, 1.00, 1.25, 1.50, 1.75, 2.00, 2.25, 2.50, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00, 4.50, 5.00, 5.50, 6.00, 7.00, 8.00, 10.00D.

- 0.25, 0.50, 0.75, 1.00, 1.25, 1.50, 1.75, 2.00, 2.25, 2.50, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00, 4.50, 5.00, 5.50, 6.00, 7.00, 8.00, 10.00D.

MFLIP + (Power)



Note: All prism flippers now include black leather imprinted cases.

FLAT PRISM FLIPPER SET Save Over 10%

Professional quality six-piece set includes: 2, 3, 4, 5, 6, & 8 diopter flat prism flippers with colored handles covers. Additionally, you receive a deluxe black carrying case (BERFC - \$30.95 value! See page 98.)

BC1275SET

Example for ordering: 5 PD flipper = 10 PD range

(5PD lens + 5PD lens)

ECONOMY PRISM FLIPPER SET

Comes with any 4 flippers of 0.25, 0.50, 0.75 or 1-7D power flat prisms, colored flipper covers and Acrylic Flipper Holder (BCPFS9 - \$32.95 value! See page 99.)

Save Over 15%

OUR MOST POPULAR KIT!

Corrected Curve Prism Flipper Set

Six-piece flipper set comes with same prism powers at BC1275SET but in corrected curve lenses, plus deluxe case and colored flipper covers.

BC1275SETCC

BC1275ESET

SAVE OVER 10%



Please specify diopter power when ordering. Prism's base and power are marked.

Prisms can be rotated, but are delivered as BU and BD.

Part #	DIOPTERS	PRICE
BC1278+	0.25, 0.50, 0.75, 1.00, 1.50, 2.00, 2.50, 3.00, 4.00D	
BC1278ESET	Set - Pick 4 Flippers plus Clear Acrylic Flipper Stand	

OPTICAL


VISION THERAPY

NEURO REHABILITATION

LOW VISION

Accommodative Rock Cards & Kit

Designed by Dr. Frank Terranova, card series measures & trains accommodative facility. Prevents patient memorization. Five 3-part sets include single letter, tumbling E, 3, 4 & 5 letter words.





The Accommodative Rock Cards Kit includes the 15 cards plus a pair of New Age red/green goggles and the new red/green bar reader (BC1323 -- pg. 90) that works specifically with the accommodative cards.

REFERENCE

SEE TO PLAY: GAZE STABILIZATION EXERCISE[©] (EYE BALLER) STPEBE

This sports vision exercise trains athletes to keep eyes up for "heads up" for distinct gaze control while performing eye, hand, body and mind activities. It improves stabilization of gaze while the body is performing athletic maneuvers, core body stabilization and trains decision making with separate right brain and left brain activities while filtering out visual noise. This exercise is also used in vision concussion rehabilitation as a clinical and home exercise. Charts are printed on white washable vinyl.

Includes: 3 Charts (2 - 8.5" x 11") and (1 - 5.5" x 8.5") & Manual

Fo	UR	Co	RN	ER	┣	łA	R	T	С	H	A	RT	Set
D	Y	М	L		H		J		U		M		NEET,
С	E	X	т		S		D		0		۷		
z	в	Α	v		X		A		W		G		Size:
U	S	к	G		Y		S		E		1	4	" x 3.75"
F	R	N	0		т		С		R		L		
-	L	F	В	J	-	R		K		С		Е	
	1	Q	X	z		D		M		X		0	P
	D	Т	κ	U		1		Y		F		Ν	Two Sided
	0	w	Ρ	H		w		н		۷		J	Front/Back
	N	S	A	۷		в		A		т		z	
	-	-		e-leet			-	-				- faits	


Charts designed for training saccadic movements, but that can be utilized for training techniques in the therapy room, such as accommodation training. These cards are printed on thick, lightweight flexible plastic material to allow the user to easily position and move them to various distances or configurations in the therapy room. Each set includes 4 charts, each one with 36 points letters. To allow for addition therapy techniques, each chart has two sides, one with colored letters and one with black.

Thick flexible and washable plastic material.


BC4DC36

Includes a QR code with link to various activities.

SIGHT WORD ROCK **USES DOLCH SIGHT WORDS**



This product lets you take accommodation therapy to the next level. After successfully completing the Hart Chart Therapy, the Sight Word Rock charts may be used while working on sight words. The charts are based on the Dolch sight word lists which vary in difficulty from pre-primer to third grade level. Sight Word Rock also trains the students to move their eyes from word to word instead of from letter to letter to better develop reading skills. There are five levels of Dolch sight words with two sheets of each level.



Set Includes: (1) 8.5" x 11" chart on white card with black print, (2) 3" x "5 laminated white card with black print, & Red/Green Goggle

Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 114

R

2

M

E

J 5

G

ASSESSMENTS

OPTICAL

LOW VISION

REFERENCE

HAIDINGER BRUSH EQUIPMENT

Macula Integrity Tester™ Version II

WITH MUCH IMPROVED EASIER-TO-SEE HAIDINGER BRUSH

MIT2 (Set)




BENEFITS OF THE MIT2[™]OVER THE ORIGINAL MIT[™]:

<u>Much</u> easier to see the Haidinger Brush

- Ability to move blue filter between 3 locations to vary difficulty of viewing the brush.
- Control speed of rotations and direction
- Pursuit training with LED flashlight
- Exo training using a distant light source

UPDATED LANTERN to reduce light dissipation with covered back.

The MIT 2[™] is an improved version of the original Macula Integrity Tester[™] for training and testing eccentric fixation. The MIT2[™] consists of a separate rotating poloroid unit and electronic control box. This allows the filters and light source to be placed behind the rotating poloroid which makes



Now Includes: CE wall plug

with international adapters.



Schiborr MIT[™] Training **MACULA INTEGRITY TESTER**" VERSION II - UNIT ONLY SLIDE SET - 11 SLIDES · wa Includes: Control Box, Rotating Developed By: Katja Schiborr, M.Sc. in Vision Science Polaroid Unit and Blue Filter Slide This expansive training slide set was designed by Katja Schiborr, M.Sc. in Vision Science, as an accessory set for the Bernell Now Includes: Macula Integrity Tester (MIT™) for use in testing and training of CE wall fixation using Haidinger's brush visual feedback. Various slide plug with designs are utilized in the set to keep patient interest, as well as nternational increase their visual skill level as they advance through therapy. adapters. The slide set includes instructions for use with the 11 slides which include: Target Number Letter Line and Maze slides KSMITSET MIT2U Made of clear, durable acrylic.

Macula Integrity Tester[™] Version I

Compact, quiet, illuminated instrument features Haidinger Brush for training and testing those with eccentric fixation.



System Includes: MIT[™], multiple dot slide, multiple ring slide, Snellen slide, MIT[™] aperture cards, blue/black goggle, cellophane reversal plate, blue cross slide, all-blue filter slide and two red pointers. Non-breakable plastic. 40W bulb included (not available in 220V). All replacement parts sold separately.

EXTRAS/REPLACEMENTS	PART #	PRICE
MIT [™] Tester Only (NO SLIDES/GLASSES/APERTURES)	BC675INC	
REPLACEMENT BLUE/BLACK GOGGLES	BC1035	
REPLACEMENT COLORED BLUB 40W	RP-675B	
BLUE FILTER CROSS SLIDE	BC1041	
BLUE FILTER SLIDE	BC1042	

MIT[™] TRAINING MAZES 8 SLIDES BC1044

This set of 8 MITTM Training Mazes has black lines with increasing degrees of difficulty. They are designed to fit in either the original MITTM or the new MITTM 2.

The object is for the patient to be able to move the Haidenger Brush through the maze, controlling the location of the brush with the fovea. This advanced technique should be done as one of the final phases of MIT[™] therapy. **Made of clear, durable acrylic**.

SINGLE OBLIQUE MIRROR STEREOSCOPE

Precision instrument used with either MIT[™] to train for macular fixation under binocular conditions.

Includes: 8 pairs of fusion cards, removable 42.5mm +6.00 lenses and 2 red pointers.

(Shown with MIT[™] tester **sold separately**, see BC675C.) Also works with new MIT2.




RAINER

RULE

A standard instrument in most Vision Therapy programs. Features progressive fixed demand. Polarized or colored lenses are not required. Helps improve low fusional and accommodative reserves, poor stereopsis,

ASSESSMENTS

NEURO REHABILITATION



APERTURE

Replacement Sports/Action Card Set (Included with Kit)

BINOCULAR TRAINING

Bernell-O-Scope™ Models

A stereoscopic/cheiroscopic instrument can be used for assessing binocular function as well as training. Features cheiroscopic drawing to break down suppression while enhancing binocular stability. Additionally, Base-In and Base-Out training improves binocular reserves. Carrier is adaptable for flat, curved or variable slides used for 200mm focal length skills test, Base-In and Base-Out as well as cheiroscopic training. Also adaptable for full-size non-variable book-form stereograms. High impact ABS plastic, the instrument is durable, lightweight and can be hand-held or set on table. +5.00D spheroprism lenses included. When inserted, these lenses make testing/training without near demand (appears as distant target). Movable occluders and wells for 38mm trial lenses. Has a headrest. Bernell-O-Scope™ is height adjustable. Testing and training cards sold separately. Original Model instrument now comes in die-cut, corrugated case for dust-free safe storage.



REFERENCE

cheiroscopic tracing forms.

BC200ESH

35

BASE-IN AND BASE-OUT SLIDE SETS FOR BERNELL-O-SCOPE™

Add More Fun and Challenge to Training with our New Landscape Line Jump Duction Base-In and Base-Out Card Sets That Come with Multi-level Suppression Checks. Another Added Option is the New Random Dot Stereogram Base-In and Base-Out Card Sets that feature hidden shapes in the random dot patterns. These hidden shapes only appear if both eyes are working together for proper fusion.



Bring interest and variation to your VT. The sets of eight or nine coated cards are bound together for use with the Bernell-O-Scope[™]. Each slide is marked with convergence demand and starts with a card with zero demand. Buy in volume and save with our discounts! Mix and match sets, or buy multiple cards sets for your therapy room.

LETTER	ΝΑΜΕ	BASE-IN	DESCRIPTION	BASE-OUT	DESCRIPTION	Volume
A & B	Landscape Line Jump Duction	BC200LLI	PD Range (0-14) 8 Cards	BC200LLO	 PD Range (3-27) 8 Cards 	DISCOUNTS
C & D	Random Dot	BC200RDI	PD Range (0-14) 8 Cards	BC200RDO	PD Range (0-17) 9 Cards	
E & F	3D Characters Jump Duction	BCBIB	PD Range (0-14) 8 Cards	всвов	 PD Range (3-17) 8 Cards 	
G & H	3D Characters White Background	BC200VI	PD Range (0-14) 8 Cards	BC200VO	 PD Range (0-17) 9 Cards 	
I & J	3D Characters Black Background	BC200FI	PD Range (0-14) 8 Cards	BC200FO	 PD Range (0-17) 9 Cards 	
K&L	3D Pictures	BC200CSBI	PD Range (0-14) 8 Cards	BC200CSBO	 PD Range (0-14) 8 Cards 	
M & N	Original 3D	BIRVT	PD Range (0-16) 9 Cards	BORVT	 PD Range (0-24) 9 Cards 	
0 & P	Clip Art Sets	BCBIC	PD Range (0-14) 8 Cards	всвос	PD Range (0-21) 8 Cards	

Bernell-O-Phore™

The Bernell-O-Phore^(The) is a unique mirror stereoscope, similar in design to a synoptophore, that can easily and quickly assess phorias, depth perception, and fusional abilities with one card set. The Bernell-O-Phore^(The) comes with a new randot card set of 26 (13 pairs) full color cards (Pictured on the right). Each pair contains either a unique sport figure, shape pattern or numbered sequence. Images appear when

fused. The remaining cards in the set are used to do base-in and base-out training. The testing maybe easily performed by the staff prior to the examination as a screening, alerting the doctor that more through testing may be required during the examination. Without the lenses placed in the trial lens holder in the eyepiece, it is test at a distance of 33cm or

+3.00 of accommodative demand. With the addition of +3.00 lenses (included in the set), there is zero accommodative demand; making it a distance test. Additional card sets will be offered in the future.

BC300 (Instrument & Cards)



BINOCULAR TRAINING Bernell Functional Binocular Assessment Test Booklet (FBAT™)




The product includes 12 test plates in booklet form with instructions, two pairs of Red/Green Glasses (Adult and Child versions), and two red pointers.

BCFBAT1

The Bernell Functional Binocular Assessment Test Booklet, or Bernell FBAT™, was designed to act as both a convenient and portable way to assess the functional binocular vision skills in many different clinic populations. It is a wonderful test when working with non-verbal patients and the pediatric population. Toddlers, or nonverbal children, will often grab in the air in front of a plate, showing good use of their binocular system and localization skill. When used with strabismic patients, the FBAT allows for grading of the decree of binocular dysfunction. There are flat fusion and full 3-D fusion targets, which take up different amounts of the visual field and with different fusion locks or reference targets. The 12 test plates were designed by clinical optometrists, Dr. Gregory Kitchener and Dr. Paul Harris, based on their experience working with these populations.



Included in the booklet for each test plate are descriptions of the binocular elements tested, expected or common responses, and what information the test provides about the visual process. All of this is printed on the side of the test facing you, while the patient is looking at the plate being described. This fun, yet informative test booklet is a great addition to any practice that sees children, non-verbal patients, or people with strabismus and amblyopia.



MORGENSTERN COLOR FUSION CARDS (DURABLE NON-GLARE FINISH, 10 CARDS PER SET) - FOR BERNELL-O-SCOPE

Mixed Convergence DIVERGENCE



As therapy proceeds along for the patient - give challenges - go into flexing from a convergence task to a divergence task - build broad ranges. SMCDC

PRE FUSION




Start every patient with this set and keep going back to it. SMPF

BASIC FUSION



The "workhouse" of building fusion, alignment and ranges.

OEP SELECTED KIT



OEP recommended kit for binocular training. A selection of 10 of the most popular cards from Morgenstern Color Fusion Cards.

BABOSET

We offer additional Morgenstern products online. To learn more information, as well as view additional products, visit our website at www.bernell.com

DIVERGENCE B



RGENCE



Basic divergence ranges of therapy - diverge on target and flex ranges

SMDA



SMBF

CONVERG

Intermediate convergence therapy, ranges of convergence, variety for interest building.

Convergence A



Basic convergence ranges of therapy converge on target and flex ranges.

SMCA

119

OPTICAL

LOW VISION



BINOCULAR TRAINING

VTP Edition Complete Training Kit Includes:

VTP Edition Mirror Stereoscope VTP Edition Adjustable Easel 12 Pair of Fusion Cards 18 Drawing Targets	• Two Red Pointers • Ruler • Manual
--	---

VTP Edition Mirror Stereoscope features magnets that are flush on the backside of the device for easy, yet strong attachment to the VTP Edition Easel. The Mirror Stereoscope and Easel are made of rigid ABS plastic and has hidden rubber hinges for smooth movement. The VTP Edition Easel is easy to adjust height and is has a heavier base than the original, which provides better stability in training.

The stereoscope is used during intensive therapy & for retention therapy. It is as important to retain the new visual abilities as it was to train them initially.

This is true space fusion trainer with a 40PD Base-In to 50PD Base-Out range. By changing the configuration, the instrument serves as a single oblique mirror stereoscope or cheiroscope. Trains for simultaneous vision, elimination of suppression, ARC, fusion, stereopsis, cheiroscopic drawing skill & development of fusional convergence/ divergence reserves. Contains front surface mirrors for less distortion of the image or angles.

Heavy-Duty Adjustable Easel: Designed to allow the user to easily adjust the angle to fit the height and comfort setting of the patient. The easel is especially useful when using the mirror stereoscope for cheiroscopic training. This easel features a heavier base as well as two strong magnets for easy, yet secure attachment of the VTP Edition Model. Included on the easel is a guide for easy monitoring of the prismatic demand of the patient when training.

3D FUSION CARDS SET

FOR THE MIRROR STEREOSCOPE


3D modeled base-in and base-out Mirror Stereoscope Fusion Cards. Set of 26 (13 pairs) full color cards. Unique 3D animals and backgrounds. Fusion cards work with the Mirror Stereoscopes (VTP-405 & BC405) models.

NEEDS

ROOM

EXAM

ASSESSMENTS



(INCLUDES BC405IO SET PLUS EASEL)

- WITH TARGETS & CARDS

BC405ES1

Home Training Kit

Easel with Adjustable Angle

Made of Black Foamcore. Lightweight. Includes BC405H plus Measuring Easel.

BC405HSET

Α

HOME **MODEL** (BLACK FOAMCORE) BC405H - WITH TARGETS & CARDS Δ

For home vision training during intensive therapy and for retention therapy. Doctors now realize that it is as important to retain the new visual abilities as it was to train them initially. Many successful vision training patients have been lost to follow-up because they thought they were cured. It is for this reason that Bernell is develops vision training equipment for the patient to buy from their optometrist and take home during or after their training. The home model is great for patients to do compliance and retention therapy at home!

BC405IOSET

B. BC405IO

IN-OFFICE TRAINING KIT

Similar design to the Home Model stereoscope, but with

the durability necessary for constant in-office use.

Office Model Mirror Stereoscope

×

This is a Wheatstone-type stereoscope; a true space fusion trainer with a 40PD Base-In to 50PD Base-Out range. By changing the configuration, instrument serves as a single oblique mirror stereoscope or cheiroscope. Trains for simultaneous vision, suppression. ARC, fusion, stereopsis, cheiroscopic drawing skill & development of fusional convergence/divergence reserves. Utilizes front surface mirrors for less distortion of the image or angles. Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

Measuring Easel

DRAWING TARGETS

C. BC405IOEA

Fusion Cards

D. BC405FC24

E. BC405CD18

LOW VISION




WSCH (CHEIROSCOPE & 3 TESTING PADS)



This instrument was designed by Bruce Wolff O.D. as a practical instrument for performance testing and visual therapy procedures. While incorporating new technology, this updated design remains as practical & durable as the original model. Simple design & low maintenance assures years of affordable use.

Easy adjustment makes it practical to administer cheiroscopic testing to patients from age 4 through adult. Testing provides information about eye-hand coordination, bi-manual symmetry, and qualities of binocularity such as suppression and binocular stability. Hard copy results can be a valuable tool for understanding functional visual difficulties.

Included is an instruction manual prepared by the Behavioral Vision Project (BVP) that provides a structured guide to a recommended testing sequence. The manual also includes discussions of observations available in testing, analysis and use of testing information, and therapy procedures utilizing the cheiroscope. The Standup Cheiroscope testing and training is taught in seminars held by the Behavioral Vision Project. An initial supply of 3 pads of testing forms (25/pad) and recording forms are included.



K.) BACT GRAPH - 160MM

L.) BACT GRAPH - 170MM

M.) PADULA DIAMOND DRAWINGS

VSBACT4

VSBACT5

VSPDD1



EXAM ROOM NEEDS

OPTICAL

121

LOW VISION

control is enhanced by using brightly colored pencils for tracing and

by using a "rubbing" motion like that used when erasing. The patient

puts a slash line at the location where the suppression has occurred.

BINOCULAR TRAINING

BAC1

BSBACTS + (1-9)

Bernell Advanced Cheiroscopic Trainer™

The Bernell Advanced Cheiroscopic Trainer[™] (BACT) was designed as a modern replacement for the original cheiroscopes used for in-office diagnostic & therapeutic applications. The unit is used with the Bernell Star Projection Forms and other cheiroscopic Tracing Forms. Cheiroscopic training on the unit helps isolate perception problems. The unit is also an excellent instrument for diplopia, fusion and accommodation training as well as for anti-suppression training. Our instrument is backlit with a focusing window that uses a dimmable light to control illumination during training. This allows for use of transparent and opaque materials to be used during training. Two knobs on each side of the unit allow for control of the angle for improved posture when training. A calibrated scale underneath the light window supports the training cards and sheets while also allowing for direct prism computations at various accommodation levels. The BACT comes with a second head set at 16" for distance/ near training. A truly professional update of a commonly used product.


This product is CE approved and works 110-240V and 50/60Hz. Comes with several interchangeable foreign electrical plug formations.

INCLUDES

- Bernell Advanced Cheiroscope
- 9 Durable Clear Vinyl Dry Erasable Slides (S1-S9)
 2 Red Pointers
 100 page padded booklets of Slides S2-S5 & S9
 Instruction Man
- (5) 100 page padded booklets of Slides S2-S5 & S9
 (VSPDD1 Padula Diamond Drawings)
 Extension Viewer

Dry Erase Marker Set

CHEIROSCOPIC TRACING SLIDES BEACTCEET



Васк

SIDE

VIEW

VTP-350

Dimmable LED

Light Switch

DELUXE CHEIROSCOPE[™] TRAINER

VTP Edition Deluxe Cheiroscope[™] Trainer was designed to provide the ideal weight needed to eliminate movement while tracing. Allows for use by right or left handed patients, as well as right or left eye viewing. The device features tilt control, which makes it easier for a patient to stay seated and trace while looking through the lenses. This can be adjusted using the black rubber knobs on the ends of the base. The device also has a mm scale on the front and back side of the base to indicate demand, as well as a smooth surface for tracing or drawing.

The instrument is made of rigid plastic with non-skid rubber feet. It also features a milled, clear acrylic viewer with +5.00 accommodation compensating lenses. The clear acrylic viewer helps to reduce visual distraction when using the device. The lenses are held in place by spring loaded threading. They can be easily removed without the use of tools for cleaning or replacement.

VTP Edition Deluxe Cheiroscope™ Trainer includes: 12 pairs (24 cards total) of special fusion cards with suppression checks and 36 tracing shapes (BC405CD18).

SINGLE OBLIQUE Stereoscope™ & Cheiroscope™



Helps train anti-suppression; builds simultaneous vision and binocular fusion at tropic angle and fusional movements back to ortho and beyond. For left or right handed patients. Removable 42.5mm +6.00 lenses.

INCLUDES: Trainer, tilt legs, eight pair fusion targets, that can also be used for tracing, two red pointers, and manual. Works with BC405CD18 targets, also.






Side View

122

BC160C (REPLACEMENT 16 CARD SET) BC160 (SET) Order Now by Phone at (574) 259-2070 or Go To www.bernell.com

BINOCULAR TRAINING



ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION


Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY



on fine slipping or fixation disparity of torsion. The combination should help symptomatic patients willing to work to improve their ranges.

Includes: 3 Laminated Cardstock Torsion Training Cards on two sheets and instruction manual.





Available in either Opaque (White Vinyl) or Transparent (Clear Plastic).

BINOCULAR TRAINING



8


XD

0

Updated set with improved targets.

FREE SPACE RANDOM DOT CARD SET BCFSRC1

The random dot images will give a unique 3-D Image when properly fused. The card set consists of 12 images that range from 4 PD to 26 PD in 2 PD increments. Some images should stand out toward the viewer and some will go away when properly fused. Which images stand out and which go away, will depend on whether the subject is crossing their eyes in front of the targets (crossed or base-out fusion) or crossing the eyes behind the target (uncrossed or base-in fusion). The standard viewing distance is 40cm.



NCLUDES Material: Cardstock, non-laminated PD Increments: 4 PD to 26 PD in 2 PD increments (12 Cards) Shape Response Panel: Cardstock, laminated

The answers (Hidden Images) of each random dot card is located on the backside of the card.

RANDOM DOT CARD SET BCFSET2

GOTTLIEB NEAR TRAINING

CARDS: STRING OF PEARLS

********************* - Harrison - Harts Hartson

These cards are designed to help with developing fine

convergence/divergence control and 3D awareness. This

series of 10 charts (5 Opaque and 5 Transparent) starts

with the largest pattern of fusional targets and decreases

The cards include instructions for convergence and divergence

training, as well as advanced techniques utilizing the cards.

in size with each chart for enhanced difficulty.

Developed by: Dr. Ray Gottlieb, OD, PhD.

0

Features 12 sets of random dot figures for use with a stereoscope or in free space. Provides variety and added difficulty to orthoptic techniques.

12 sets of random dot figures Printed on: Cardstock


Card Size: 3.125" x 3.125" Includes: Instruction manual and answer key



Consists of two copies of a card with seven rectangles with a pair of eccentric loops in each. Each loop is formed out of arrow shapes. Inside each rectangle are parts of a simple arithmetic sentence and one or two words. These are designed to help with management of convergence dysfunctions.

Detailed instructions and goals for using cards included. TNCDSET


RANDOM DOT PAD (Pad of 36)




PROGRESSIVE FREESPACE FUSION CARDS

Two white and two clear cards with eccentric fusion images that get progressively farther apart to increase fusion demand for eso & exo training. Each line has separation listed. Demand varies from 8-26mm.

PFF





GNTCSP

2

100


Perceptual





Game Includes: (25) one inch cubes of various colors, (25) 3/4" white cubes, (24) of each set of cubes are printed with A-X and 1-24, symbols, and dots. One of each set is blank. Comes with players guide and two Symbol Cards for matching. Magnetic box - 8" x 5" x 2". Get Two Extra Game Cards when you buy from Bernell!

Establishing a working memory provides for a more powerful approach when engaging in the activities of acquiring academic instruction. This game works well with patients in many different developmental levels. The sequence of cubes develops the "when" and the spatial pattern of the cubes helps develop the "where". The cognitive challenges incorporated into this activity provide skill in learning the "what". The game develops



CARL HILLIER'S MULTI-MATRIX GAME

better skills and automaticity with these three crucial aspects of visual information processing.

Skill in sequential and spatial information processing makes life and learning easier.

CREATIVE COLOR CUBES SUPER SET

This set allows you to work to improve perceptual skills by creating and matching three dimensional patterns and designs with this set of wooden cubes in six brilliant colors. Set includes 100 wooden cubes and 18 color cube pattern cards. Increase visual discrimination, fine motor skills and understanding of basic math concepts with this set of 36 pattern cards and 100,

VTMMG

training course

Creative Color Cubes offer hours of creative, constructive play that builds essential math and pre-reading skills:

- Patterns and Sequences	- Visual Perception
- Spatial Reasoning and Logic	- Geometric Concepts

Contents: 100 Wooden Cubes in 6 bright colors (cubes are 3/4"), 33 Design Cards, and Activity Guide



LOW VISION

PERCEPTUAL

WOOD PARQUETRY BLOCKS


PERBLKS

This training method uses wooden blocks to guide the student from the basic to the most complex perceptual skills. Improves visual analysis, handeye coordination, motor control, visual motor and visualization ability.

LOW VISION

REFERENCE

EXAM ROOM NEEDS



Parquetry Workbooks on pg. 131

MULTI-COLORED VISION BLOCKS (SET OF 32)



Now includes a 46 page activity book.

The Vision Blocks are used to create patterns, designs, sequences and explore angles and geometry with shapes designed to build essential math and pre-reading skills. These 20mm wooden blocks are similar to "Parquerty" blocks. They help develop matching and discrimination

abilities using single and 3 dimensional figures in oneto-one correspondence. They provide for: Visual discrimination, matching, spatial rotations, tactile feedback. L67BLK32



Set includes: (32) 20mm thick wooden blocks: (8) triangles, (8) squares, (16) parallelograms, activity workbook and a wooden box for storage.



CLEAR TRANSPOSITION FACILITATORS

Designed for use with parquetry or vision blocks, these clear acrylic facilitators allow for transposition of patterns while training.



Available in round or square. These come in two sizes: 10" or 12". SIZE/SHADE DADT # DDIC

SIZE/SHAPE	PARI#	FRICE
10" ROUND	L76R10	
12" ROUND	L76R12	
10" SQUARE	L76S10	
12" SQUARE	L76S12	



127 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com


PERCEPTUAL

WACHS ANALYSIS OF COGNITIVE STRUCTURES (WACS)[®] Test Kit COMPLETE KIT WITH MANUAL

Developed By: H. Wachs, OD & L.J. Vaughn, MA

Standardized inventory based upon Piaget's theories of measuring cognitive development in terms of body and sensory thinking in children ages 3 years to 5 years 11 months. Individually administered and easily scored in 30-45 minutes by assessment professionals or classroom teachers. Focuses on period during which child begins to acquire symbolic thinking and has been used successfully with both verbal and non-verbal children, as well as adults, whose development age ranges in the 3-6 age period. See page 49 for complete details.

The kit consists of: 15 clusters of task divided into subtext made up of a variety of manipulative, visual and sensory materials, as well as body movement tasks grouped into 4 areas of assessment: Identification of Objects, Object Design, Graphic Design, and General Movement.






This interactive teaching tool allowing patients to create patterns using their hands while only seeing the pattern laid out by their doctor or therapist. A great tool for encouraging thinking skills, including perception, prediction, analysis, and synthesis. Piaget's theories of measuring cognitive development in terms of body and sensory thinking is the brainstorm behind this ingenious box.

Openings for hands on both sides with curtains allow the patient to work inside the box without seeing their hands. Magnetic hinged back doors allow for self-check inspection while the open viewing area in front allows for monitoring by the doctor or therapist.

The box size is approximately 10"D x 18"W x 9"H inches. Does not include blocks.



A game of tactical tracking and speed! Use problem solving and a keen eye to find the path connecting 2 blocks. Then recreate the path as fast as you can with your own pieces. Watch out, though, your mind could be playing tricks on you! Three levels of play keep you challenged.

 Includes 44 pieces, 25 double-sided pattern cards and Activity guide

· Cards measure 8' 'L x 8' 'H

LER9279





Color-matching fun for everyone! Place any card to start. Each card includes different combinations of up to four colors, so be sure to look at your cards closely. How many colors can you match with your next play? Match the most colors, and win! This addictive, make-you-think strategy game builds matching, critical thinking, and addition skills.

Includes 40 cards and Activity Guide. Cards measure 3" x 3". For 2-6 players.

KANOODLE® GAME (KANP12

Pick a puzzle from the book, set up the challenge, and then place the remaining pieces to solve the puzzle. Hundreds of possible combinations, but only one completes the puzzle!

200 Brain-building puzzles

WACSKIT

- Tricky 2-D rectangular and mind-boggling 3-D pyramid puzzles at multiple levels of difficulty
- All pieces store inside portable, pocket-size case for easy cleanup and play anywhere
- For 1 player
- Includes 12 puzzle pieces and 44-page illustrated puzzle book

Teaching Tac-Tiles™

LER9075

Tactile learning is a bag full of fun! Colorful, engaging pieces with rich textures (smooth, bumpy, ridged, and rough) provide a great way to engage a child's senses and strengthen fine motor skills. Reach into the bag to find two matching shapes or textures only by touch, or follow the Activity Cards. Learn shape identification while building vocabulary, and early math skills.



Includes: 20 pieces in 5 shapes and 4 textures, and 10 Activity Cards. Shapes measure 2"H. Aaes: 3 - 7



BAXSTROM SIMULTANEOUS PERCEPTION CARDS

GEOMETRIC SHAPES HDBCGFV2

These cards are designed to be used with red/green anaglyphic filters (not included) and provide good retinal rivalry for enhancing binocular vision training. as well as a wide variety of training activities, which work well in any vision therapy program. These cards allow for training of Ocular Motility Training (Tracking Skills), Accommodation Training (Focusing Skills), and Visual Perception/Visual Training (Visual Information Processing Skills).

The deck provides 3 representations of each shape in solid lines and two levels of broken lines for training visual closure skills. These can be used for simple one tone matching, visual memory and form constancy activities. Developed by Curtis Baxtrom, OD, FCOVD and Thomas Headline, COVT



ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

MISCELLANEOUS



Mental Blox[®] Go!

Help kids hone their critical thinking skills even when you're on the go. The Mental Blox -On The Go set combines the boredom-fighting fun of a traditional travel activity set with the trusted educational challenges. Contained within the sturdy, travel-friendly carrying case is a 16-piece puzzle ready to tease brains and challenge critical thinking skills during long flights or short road trips. Each set comes with a puzzle book filled with 10 different types of challenges, including games and activities designed for both one and multiple players-over 30 challenges in all targeting different skills and subjects.

LE9286






Mental Blox[®] Critical Thinking Game

Put your critical thinking to the test and fire up the strategic problem solving! Stack, pattern and match the multicolor shapes. Don't be fooled by appearances, these endlessly enjoyable blocks mean brain-boosting business. Mental Blox help children practice following verbal instructions, asking questions and build critical thinking skills. Activity cards provide children with examples of what to build. Children can recreate what is on the card or can give verbal instructions to a partner or group. Pieces feature four different shapes and 3 different attributes for a real challenge!







ELECTRONIC CONVERTER

This voltage converter converts foreign electricity from 220-240V to 110-120V. This 2000 watts max. overseas adapter converter features 2 modes. Adaptor mode for dual voltage appliance such as 100-240V devices, as well as single voltage devices. (Please remember to select the proper mode for your device "Adapter (dual voltage) or Converter (single voltage)" before using this product). Comes with LED power indicator, plug lock system, over-current (overload), high-temperature, short-circuit and surge protections built-in.

Picture shows different adapter options on converter. Adapters slide out of device using switches on side.

CVTR (Includes Case)



Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

NEURO REHABILITATION

LOW VISION

REFERENCE

129

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION



offers nine height adjustment positions to provide a variety of

viewing angles. Folds flat for compact storage and portability.

Holds standard letter-size ring binders. Supports up to 10 lbs.

This ergonomic slant board has been designed incorporating the principles of visual hygiene. The Slant Board, with it's sleek and functional design, is easy to use on a desk, table, or in bed. All in all: Safe, resistant, light and attractive. When used for reading and writing, it serves to support proper body position and working distance and to discourage adverse postural adaptations. Adult Version is pictured to the right. Child Version has figures across the top.

VISUAL HYGIENE BOARD FROM VTE

Color: Platinum

Sheet Capacity: One Book/Pad

Adjustability Features: Height



E/HRD1A (Adult) E/HRD1C (Child) E/HRD2A (Adult - Pkg. of 10) E/HRD3C (Child - Pkg. of 10)

VISUAL EDGE[™] SLANT BOARDS IN COLOR AND ACCESSORIES

The new and improved Slant Board is a slightly thicker model than the old design and includes strips of Velcro backing for a more versatile paper catch and for easy attachment to our new Lap Pillow. The thicker board will be sturdier and will last longer while the Velcro backing makes the Slant Board compatible with our other products shown below.


Buy 6 for \$37.85 each and save 5% Buy 12 for \$34.95 each and save 12%



WORKBOOKS - DISCOUNT ON ALL BERNELL WORKBOOKS ON THIS PAGE

These workbooks were designed to look a part of a professionally run therapy program. You can reduce your costs

Each workbook is spiral bound with



131 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

WORKBOOKS ernell Activity/Therapy Workbooks **AVAILABLE INDIVIDUALLY OR AS A** SET OF ALL FOUR: BLACK AND WHITE BERNELL AMBLYOPIA WORD WORD SEARCH THERAPY THERAPY BERNELL THERAPY SEARCH THERAPY CMB (Color Maze Book) for the Ambiyopic Patier Progressive Word Searches for the Antalyopic Patient BWMB (B/W Maze Book) CWSB (Color Word Search Book) BWSB (B/W Word Search Book)

BERNELL is excited to offer these workbooks designed to make vision therapy a little more fun and interesting. The word search and maze books are extremely useful in the treatment of amblyopia and suppression. They are fun and interesting for the child to do, yet designed to get progressively more challenging with each page. The word search and maze books have puzzles that get smaller and smaller to train the weaker eye to see better and better. Both workbooks are available in a black print version for use with active patching or a Red/Green version that allows the subject to do monocular training in a binocular field with Red/Green glasses. Also can be used for Figure-Ground work.



Ann Arbor Workbooks



A). WORD TRACKING: HIGH FREQUENCY WORDS WORKBOOK AT8706

High frequency tracking of commonly used words in the English language. Provides visual training in tracking, left-to-right directionality, spelling, recognition, matching & selecting. Writing exercises using these words reinforce vocabulary development & written language skills. 64 pages.

× BERNELL

B). SENTENCE TRACKING: WORKBOOK

8-1/2" x 11"

SOFT COVER

DETAILS:

ERNELL

AT7483

Designed to teach scanning a line of type to improve reading, word recognition and discrimination skills. Sentence Tracking aids in vision training by correcting problems with left to right eye movement, omissions and reversals. 3rd to 4th grade reading level. Type size ranges from 16 to 10 point, Reproducible,

C). SENTENCE TRACKING: HIGH FREQUENCY WORDS WORKBOOK AT9214

Uses the same model sentence approach as the time-honored Sentence Tracking program but includes the Francis-Kucera word list that is used in the best-selling Word Tracking: High Frequency Words. 2nd to 3rd grade level. Type size ranges from 18 to 10 point. Reproducible.

D). LETTER TRACKING: WORKBOOK

Letter Tracking is designed to improve visual discrimination in students. Problems that may be caused by the lack of good visual discrimination include reversals, skips, left to right directionality, omissions, substitutions, and excursions over the page with the eyes. For beginning readers, remedial cases at all levels, and adults. 24 to 19 point font size

E). SYMBOL DISCRIMINATION & SEQUENCING WORKBOOK AT85157

Designed to improve visual discrimination, visual memory for sequence, noting likenesses and differences, visual matching, and left-to-right progression. For pre-school/kindergarten level. Will help to develop important skills of visual memory and left-to-right progression.

F). SYMBOL/LETTER TRACKING WORKBOOK

Enables each learner to track symbols, numbers, lower case and upper case letters in groups of six through 20. The book is designed to improve visual discrimination, symbol memory, and left-to-right directionality while correcting problems with reversals, omissions, substitutions, and additions. Appropriate for non-readers and beginning readers. Type size ranges from 24 to 12 point. Reproducible

AT7351

AT7327

REFERENCE

ASSESSMENTS

OPTICAL

VISION THERAPY

SAVE OVER 15%

WSMSET (Set of All 4 Books)

NEURO REHABILITATION LOW VISION

SYNTONICS SHEARER COLOR SENSITIVITY TEST (SCST) SYSTEM SCST

SCSTG-B3

OPTICAL

SCSTG-GN3

The purpose of the SCST system is to both objectively and subjectively test for positive outcomes in three main categories of color sensitive patients. These groups consist of patients with concussion, reading disabilities and the severe headache or migraine sufferer. The percentage of patients suffering from visual discomfort after concussion is quite high; upwards of 70% in some studies. The percentage of these patients that can be helped with color filters or lenses is currently under investigation. However, some patients may get relief from their symptoms who do not respond to more traditional therapies.

EYELIGHTS

Includes Clear Lens

SCSTG-R3

1	SCSTG-Y3	SC	STG-GY3		
			Τιντ		DD
		Light	MEDIUM	Dark	IK
	GREEN	SCSTG-GN1	SCSTG-GN2	SCSTG-GN3	
	Red	SCSTG-R1	SCSTG-R2	SCSTG-R3	
	BLUE	SCSTG-B1	SCSTG-B2	SCSTG-B3	
	Yellow	SCSTG-Y1	SCSTG-Y2	SCSTG-Y3	

SCSTG-GY1 | SCSTG-GY2 | SCSTG-GY3

SCSTF5 (Flipper Set)

Some patients simply respond better to specific colors of tinted lens than simple gray lenses, which act as neutral density filters. Gray lenses simply reduce the overall amount of light reaching the eye. It is still unclear exactly how certain wavelengths of light work to achieve this, but it is thought that it changes the timing between the m-cell and p-cell system, or more commonly referred to as the central and peripheral vision.

133

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

SCST FLIPPERS

The SCST flippers contain the most common colors chosen by color sensitive patients. The use of color overlays is an alternative method for determining the

GRAY

optimum color, but studies have shown that the difference between colors chosen with overlays may vary significantly compared to the colors chosen when affecting the full field with lenses. Ideally the most exact way to choose a color is with a colorimeter, but the cost of the instrument and the expensive of custom tinting may be prohibitive for many practices.

> A much less expensive and yet effective way of testing and treating the color sensitive patient is the combination of diagnostic flippers and confirmation glasses of the SCST system.

Note: The red colored tint in the kit is known as FL-41 which seems to help many with migraines.


 Flipper Set Includes:
 5 Flippers & Protective Case

By stimulating the weaker side of the brain you allow it to become stronger, thus improving overall performance both mentally and physically. Eyelights can be used by anyone since everyone has one side of the brain that functions at a higher level than the other. Eyelights are being worn by children with learning disabilities in order to keep them more attentive and focused in the classroom. Athletes are wearing Eyelights to increase performance on the field or to push through a plateau. Anyone can wear Eyelights to deal more effectively with stress, depression, or neurological disorders.

EYELIGHTS HAVE BEEN KNOWN TO HELP WITH THE FOLLOWING:

- Cellular Immediate Early
 Gene Response
- ADHD
- DyslexiaStroke
- Multiple Sclerosis
- Autism

Eyelights target the weaker functioning side of the brain via the non-dominate eye. Light stimulation directly to the non-dominate brain causes and excitatory barrage to travel to the mesencephalon, the most metabolic area of the brain, where an increase in cellular function occurs.





CONFIRMATION GLASSES

RICE

Syntonics

DELUXE FCFTESTER SOFTWARE PACKAGE

The Computerized Functional Color Field Tester (FCFTester) software system provides for quick, efficient, and reliable testing of the functional color visual fields. The software allows for the measurement of Form Fields (functional campimetry), the Color Fields (Red, Blue, & Green), and Blind Spot.

A key benefit to this system is the vast number of settings and customizations which allow you to create uniform test conditions for your patients, as well as the ability to store test parameters to ensure testing accuracy and consistency with patients when retesting. A built-in calibration module allows you to run the software on monitors 15" up to 24". A variety of new enhancements now are included with the new Deluxe FCFTester software package, including the ability to monitor eye position during testing. For existing FCFTester customers, you can purchase this module as an upgrade. Contact us for details.

In addition to ongoing software enhancements, the viewing hardware for the FCFTester is now constructed of high-impact ABS, which provides a durable, professional appearance. Also, the viewer now features a trial lens well, as shown in the image to the right.

Save time and improve accuracy by upgrading from your manual color field testing equipment to this professional software system, which is quickly becoming the new standard in color visual field testing for both its ease of use and value in the exam room!

Interested in upgrading your old viewer hardware?

The new deluxe software packages will include the new ABS plastic viewer (shown in the picture). If you are a customer who purchased the original RFFT software with the original acrylic viewer and are interested in upgrading your viewer, then you can purchase the new deluxe viewer individually.

DFCFTVO (Viewer Hardware Only)

CERIUM PREMIUM OVERLAY TESTING SET CONTINUE

The Cerium Premium Overlay Testing Set consists of 2 sets of A5 overlays in each of the 12 colours equally positioned around colour space. The color is selected by process of elimination. Assessment should be carried out only in the manner indicated on the record sheet. Full instructions and assessment sheets are included. The pack is supplied in a high quality presentation folder.

The Wilkins Rate of Reading Test and DVD training by Professor Arnold Wilkins is also included. The Wilkins Rate of Reading Test was designed to objectively measure the benefits obtained by individuals who suffer with Visual Stress, viewing text through Cerium Coloured Overlays. It is not a measure of cognitive ability or intelligence but compares the effects of reading simple words with and without the selected overlay. The patient is simply required to read words aloud while reading times and errors are noted. The test uses simple words, familiar to young children. The words appear in each line of a block of text but in a random order as we are not assessing the child's comprehension skills but purely the reading speed.

Practitioners find the test easy and quick to administer and the results reliable and repeatable.







COLORED OVERLAYS

The overlays consist of a robust transparent sheet that has a non-reflective surface, but one that allows the print to be clearly seen. The overlays may be superimposed on each other, but no more than two overlays should be used any one time due to the increasing density of tint. The edges are rounded for safety.

CPO+ (SINGLE A4 FULL SIZE OVERLAYS)



DFCFTS

Customer Service Team to learn more information.







SOLAN BLUE FILTER SHEETS

Special blue vinyl sheets as used in the work on Scotopic Stress Syndrome by Dr. Harold Solan.

> Includes: Six Sheets Size: 8-1/2" x 10-1/2"



CPO+ (BUY 5 OR MORE)

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION LOW VISION

SPATIAL LOCALIZATION & PRISM TRAINING



LOW VISION

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



LOW VISION

Prism Training Bernell Prism Training Goggles

Bernell is Your Source for Prism Training Goggles

ADULT SIZE (48mm)

DELUXE GOGGLES

CHILD SIZE (38mm)

DELUXE GOGGLES

WITH NOSEPADS

WITH NOSEPADS

Adult Size (48mm) Standard Goggles

48mm Sizes have PD approx. 64

CHILD SIZE (38mm) Standard Goggles

38mm Sizes have PD approx. 55

Allows patients to train without instruments while reading, working, etc. Prisms are round, with flat or corrected curved plastic, and can be rotated to any position with finger screws. Power and base are marked on prisms. Adult size (48mm) goggles now come with an axis scale marked on front to more accurately set the angle of prism. We've made additional improvements on both frame sizes, including widening the cord and adding Velcro for easy on/off! This improvements are in addition to improved size screws which make for easier tightening of the lenses on both frames. Both now come with a 360 axis scale! 27" ELASTIC STRAP FOR PRISM GOGGLES WITH MINI LOCK (RP-1148EL


Flat Prisms: Allow training with maximum prismatic distortions of the patient's field.

Standard Flat Prism Goggles			Deluxe Flat Prism G	OGGLES WITH	Nose Pads
PART NUMBER	DIOPTERS	PRICE	PART NUMBER	DIOPTERS	PRICE
BC1138+ OR BC1148+	0.25, 0.50, 0.75, 1-7D		BC1138+NP OR BC1148+NP	0.25, 0.50, 0.75, 1-7D	
BC1138+ OR BC1148+	8 or 10D		BC1138+NP OR BC1148+NP	8 or 10D	
BC1138+ OR BC1148+	12, 15, or 20D		BC1138+NP OR BC1148+NP	12, 15, or 20D	
BC1138+ OR BC1148+	25D or 30D		BC1138+NP OR BC1148+NP	25D or 30D	
BC1138FR	38mm Empty Frame		BC1138NP	38mm Empty Frame	
BC1148FR	48mm Empty Frame		BC1148NP	48mm Empty Frame	

Corrected Curved Prisms: Ground lenses correct for optical and chromatic errors. These lenses are similar to prescription glasses to allow more accurate training.

STANDARD CORRECTED CURVED PRISM GOGGLES			DELUXE CORRECTED CURVED PRISM GOGGLES WITH NOSEPADS			
	PART NUMBER	DIOPTERS	PRICE	PART NUMBER	DIOPTERS	PRICE
	BC1138+C OR BC1148+C	0.25, 0.50, 0.75 1-7D		BC1138+CNP OR BC1148+CNP	0.25, 0.50, 0.75 1-7D	
	BC1138+C OR BC1148+C	8 or 10D		BC1138+CNP OR BC1148+CNP	8 or 10D	
	BC1138+C OR BC1148+C	12D		BC1138+CNP OR BC1148+CNP	12D	

ROUND FLAT & CORRECTED CURVED PRISMS



These bevel-edged lenses can be used in trial frames (38mm), flippers (38mm), or prism training glasses (both 38mm and 48mm). Base can be turned to any desired position. Powers and base are conveniently

marked on the prisms.

137



Round Flat Prisms

Some doctors prefer these to better accentuate prism distortion.

P				
PART #	DIOPTERS	PRICE		
AR38+ OR AR48+	0.25, 0.50, 0.75 1-7D			
AR38+ OR AR48+	8 or 10D			
AR38+ OR AR48+	12, 15, or 20D			
AR38+ OR AR48+	25 or 30D			
To Order: AR38+power or AR48+power		POWER		

CORRECTED CURVED PRISMS

Train without distortions and	aberrations for just a	couple more dollars
— "		_

Part #	DIOPTERS	PRICE	
CR38+ OR CR48+	0.25, 0.50, 0.75 1-7D		
CR38+ OR CR48+	8 or 10D		
CR38+ OR CR48+	12D		
To Order: CR38+POWER OR CR48+POWER			

Half power prisms available in corrected curved lenses as special orders.

Contact us for pricing. Sorry, we can only supply CC prisms up to 12PD.

VISION THERAPY

PRISM TRAINING



Deluxe Prism Goggles Kits and Accessories

Lightweight, premium frame designed to fit comfortably on both children and adults. Frames work with both our flat and corrected curve 48mm prisms. Prisms snap securely in place, yet can be easily rotated by hand to the desired location. No more thumb screws!

The frame has a lens well on both the front and back of the frame. It also features an internal groove for adding red/green filters to the goggles. The frame has white tick marks every 30 degrees. They come with a 1/2" wide flat elastic band and elastic adjustment toggle for a secure, yet comfortable fit.

VTP EDITION PRISM GOGGLES WITH				
FLAT PRISMS	DIOPTERS	PRICE		
	0.25, 0.5, 1-7PD			
VTP-1148+	8 or 10D			
VII 1140	12, 15, or 20D			
	25D or 30D			
VTP-1148FR	Empty Frame			
 Power and base are marked on prism. 360 axis scale with 30 degree markings. Includes microfiber drawstring pouch for frame storage. 				
VTP EDITION PRISM GOGGLES WITH				
CORRECTED CURVE PRISMS	DIOPTERS	PRICE		
	0.25, 0.5, 1-7PE)		
VTP-1148+C	8 or 10D			
	12			

" wide ble fit. Shown with flat prisms. Allows patients to train without instruments while reading, working, etc. Flat lenses and corrected curved prisms are round, CR39 plastic, and can be rotated to any position by hand. Power and base are marked directly on the prism. Goggles

FEATURES:

Prisms snap into place!

No more thumb screws!

· Designed for both children and adults.

Add 38mm trial

VTP Edition Prism

lenses to the

· Easily rotate prisms by hand to desired location.

· Works with both our flat and corrected curve 48mm prisms.

TRIAL LENS INSERTS

VTP-1148TI

Prism

Kits

Goggles

Goggles! Fits into existing 48mm goggle wells.

Sold as a pair of plastic inserts.

RED/GREEN FILTERS FOR VTP EDITION PRISM GOGGLE FRAME 1 Red and 1 Green 48mm round filters. 5 Filters snap into center grooves of VTP Edition Prism Goggle Frame.

come with an axis scale marked on the front and back to more accurately set the angle of the prism.

VTP-1148RG

VTP Edition Deluxe Prism Goggle Complete Kit - 28 prisms

VTP-1148KIT28

Powers Include: 0.25PD, 0.50PD, 1PD, 2PD, 3PD, 4PD, 5PD, 6PD, 8PD, 10PD, 12PD, 15PD, 20PD, 25PD


VTP Complete Kit has (28) prisms, (6) frames, 6 pairs of 38mm trial lens inserts (12 total), 6 pairs of red/green lens inserts (12 total) and a lens cloth color may vary. All items are housed in a custom plastic case that is lined with pre-cut foam to keep prisms safe & clean.

VTP Edition Deluxe Prism Goggle Mini Kit - 16 prisms

VTP-1148KIT16

Powers Include: 1PD, 2PD, 3PD, 4PD 5PD, 8PD, 10PD, 15PD

VTP Mini Kit has (16) prisms, (4) frames, 4 pairs of 38mm trial lens inserts (8 total), 4 pairs of red/green lens inserts (8 total) and a lens cloth color may vary. All items are housed in a custom plastic case that is lined with pre-cut foam to keep prisms safe & clean.



LOW VISION

REFERENCE


PRISM TRAINING & EQUIPMENT

E/RYPA+



VTE YokeD Rotating Prisms Adult and Child Size Available in PD: 0.25, 0.50, 1, 2, 3 E/RYPA+ Available in PD: 4, 5, 6, 8, 10 E/RYPA+ Available in PD: 4, 5, 6, 8, 10 E/RYPA+ Available in PD: 12, 15, 20, 25

E/FRRYPA



Free Space Septum Training Mirror

E/FRRYPC


A septum with two mirrored sides which is used for binocular training, antisuppression and fusional activities as well as fighting neglect in stroke/brain injury rehabilitation. This handheld two sided mirror has a custom handle for comfort, as well as to prevent finger prints from damaging the mirror surface. Includes activity guide for training.





Children's Size

Materials: Acrylic mirror with plastic handle Mirror Size: 8" x 5", 3/16" thickness. Overall Length: 12" Long



The Hemianopsia buttons help in clinic or hospital settings to let others know on which side the patient has troubles. Buttons state "I have trouble on my RIGHT (or LEFT) side." The correct side is blackened. Sold individually or both.

AFR225R (Right Side) AFR225L (Left Side) Buy 2+ for \$2.45 each

Rummel Hemianopsia Guides™



The Hemianopsia Reading Guide helps increase reading speed and comprehension. Without the guide, those with field loss or neglect will often leave off part of a line & miss or change the meaning of the sentence. With the guide, one reads until coming to stop or returns back until the "go" barrier.

Reading guides are made of durable vinyl.

ER049


ASSESSMENTS


EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY





Printed on Rigid White Vinyl. LVCDS LVCD2 Low Vision products are returnable with a restocking fee of 15%. Unfortunately, some doctors were ordering 5 mags for a patient to test and then trying to return 4! Order Now by Phone at (574) 259-2070 or Go To www.bernell.com


Very easy to use!

a font to make it easier for the low vision patient. Additionally, a black and white print

on the back can be used for patients who struggle with color vision or color confusion.

140

ELECTRONIC MAGNIFIERS





3.5" Portable Digital Magnifier

Easy to use, small portable magnifier has a 3.5" LCD screen and is very easy to use and take about.

Magnification ranges from 2X to 24X, for quick viewing. Other features include freeze frame for 1 image, and 15 different color options. Output to a TV available using an included AV output cable. 6 hour battery life.

PDM

VISUM ECARE 3.5" - ELECTRONIC MAGNIFIER

Ecare 3.5" is a full-function yet lightweight video magnifier that's easy to take with you anywhere you go. It offers a magnification range of 2.5x - 17x. Display images in your choice of color, B/W, reverse or one of the other false colors on the 3.5-in. widescreen LCD, and use the Freeze Frame function to capture an image onscreen and examine it more closely.

Ecare 3.5" is a handy tool to help those with low vision read food labels, tags in stores or restaurant menus, and is also versatile enough for longer reading tasks. 2.5 - 17x Magnification





Key Features:

Low Vision products are returnable with a restocking fee of 15%. Unfortunately, some doctors were ordering 5 mags for a patient to test and then trying to return 4!

141 Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com



EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION



LOW VISION

HANDHELD MAGNIFIERS LED LIGHTED 2.5x Power



Low Vision products are returnable with a restocking fee of 15%. Unfortunately, some doctors were ordering 5 mags for a patient to test and then trying to return 4! Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com 143

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

NOISIN MO

STAND MAGNIFIERS


EXAM ROOM NEEDS


ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION





ULTRADOME

distance, in perfect focus!! Size: 4X Magnifier - 2" round

SE402 (4x)

UltraDome® glides effortlessly across any surface for ease of use. View magnified objects from any

Includes a FREE Polishing Pouch for storage and cleaning.



Yellow guide allows patient to keep better track when reading. Bar is 6" long and magnifies 1.5x.

LED ILLUMINATED FOLDING MAGNIFIER 3X

- · 3x Power Crystal Clear Acrylic Lens (2.25 inch lens diameter)
- · 3 Super Bright LED lights Illuminate Your Viewing Area
- · Extendable neck for easy hands free usage
- Fully Collapsible Design for Safe Storage and Easy Portability Uses 3 AAA Batteries (not included)

LM10LFM

(RED & BLUE) (C585012R (Red) The VTM's have been specifically designed to aid patients suffering from Dyslexia or similar reading difficulties. The VTMs help patients to concentrate their vision on individual words or sentences

Lens Diameter: 50.5mm

Image Vergence: -4.0D

and

COIL VISUAL TRACKING DOME MAGNIFIER 1.7X



remove the confusion caused by nearby text or pictures. Magnification: 1.7x

Effective Focal Length: 52.90mm Field of View on Page: 143 x 104mm

Low Vision products are returnable with a restocking fee of 15%. Unfortunately, some doctors were ordering 5 mags for a patient to test and then trying to return 4! Order Now by Phone at (574) 259-2070 or Go To www.bernell.com 144 REFERENCE

LOW VISION

SPECTACLE MAGNIFIERS



OPTIVISOR® 1.75x, 14" FOCUS 2x, 10" FOCUS 2.5x, 8" FOCUS 2.75x, 6" FOCUS 3.5x, 4" FOCUS DDA3 DDA4 DDA5 DDA7 DDA10 A. LIGHT ATTACHMENT - Light with AA battery pack both attach to visor. B. OPTI-LOUPE ADD - Monocular attachable lens. Adds 2.5x to all Optivisor models C. REPLACEMENT LENSES DDALENS (2,3, 4 & 5) DDALENS (7 & 10)

Ideal for work or hobbies that require magnified viewing. Precision binocular magnifier permits unrestricted 3D vision while reducing eye strain & freeing the hands. Optical glass prismatic lenses, ground & polished to precision optical standards. Visor can be worn over prescription glasses. Pivots to allow visor to be tilted out of the way for distance vision. Headband conforms to fit all head sizes, and is comfortable even when worn for long periods.

Deluxe Binocular Headband Magnifier with LED Light Includes Five Interchangeable Lens

This lightweight, binocular magnification device allows for clear 3D viewing while limiting eye strain and allowing for full use of hands. The LED light and interchangeable lenses create additional ease of use. This device can be worn over prescription glasses. The adjustable headband allows for the user to fine-tune the fit to their desired comfort. Lenses are constructed of high quality ABS molded acrylic, providing clear and undistorted texts and images. One added benefit of the design is that it allows

the wearer to set two lenses at a time for easy flipping or switching from one magnification power to another, or doubling the magnification.

I HMAG1



Single-plate lens magnifications: 1.0x, 1.5x, 2.0x, 2.5x, 3.5x Bi-plate lens magnifications: 3.0x, 4.0x, 4.5x, 5.0x, 5.5x, 6.0x Battery: 3 x AAA battery (not included)





Low Vision products are returnable with a restocking fee of 15%. Unfortunately, some doctors were ordering 5 mags for a patient to test and then trying to return 4!

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com

OPTICAL

VISION

≥ 0

REFERENCE

145

EYE MODELS



SMALL EYE MODEL

This is a basic, economical eye model. It includes an information card, whole eye cut in half, cornea, lens, and realistic vasculature. Made of plastic, measures about 4 inches in diameter, and is mounted on a slightly larger plastic base.



LARGE EYE MODEL

This model features a whole globe which can be taken apart to reveal internal structures. In addition to the standard anatomy this model has added cords to represent myopia and hyperopia. Made of plastic. Measures 6 inches in diameter, and is mounted on a black plastic base.



IOL EYE MODEL

With this model, you can demonstrate different types of IOL implants and procedures. Comes with six IOL lenses. It is made of plastic, measures 3.75 inches in diameter, and is mounted on a larger white, opaque plastic base.

AMIOLEM



CATARACT EYE MODEL

The Cataract Eye Model is perfect for educating patients prior to surgery. Shown with optional lens stabilizer. Made of hard, durable plastic and measures 3.75 inches in diameter mounted on a slightly larger plastic base. Included are five interchangeable lenses that show various types of conditions including sub-capsular, capsular, mature, cortical, and nuclear.

AMCATEM



VISION THREATENING **DISEASE (VTD) EYE MODEL**

This model demonstrates the three most common eye diseases: Glaucoma, AMD, and Diabetic Retinopathy. The model is realistic and brightly colored. Perfect for patient education. Made of plastic, measures 4 inches in diameter, and is mounted on a slightly larger white plastic base.

AMVTD



POSTERIOR VITREOUS DETACHMENT (PVD) EYE MODEL

Use this model to demonstrate what happens in a PVD. The model shows a detached vitreous including the Vogt Ring. Also shown are a cross section of the eye in three layers; the retina, choroid, the sclera with veins and arteries. Made of plastic, measures 4 inches in diameter, and is mounted on a slightly larger white plastic base. The cornea, optic nerve, and insertion of muscle are shown.






ASTIGMATISM EYE MODEL

This is the perfect model to educate patients in astigmatism, myopia, and hyperopia. Use the sliding beads to demonstrate the focal length of the eye in different meridians. 6" in diameter on mounted on plastic white base.

AMAEM



LASIK EYE MODEL

Model helps you to explain the different aspects of Lasik and PRK surgery & saves the doctor time. Realistic flap in front of the model simulates the incision. Model is 14.5cm x 14.5cm in size. Made of hard durable plastic. Very portable to move from exam room to exam room.

G15120



TIMBERLAKE EYE MODEL

This model was developed to improve the teaching of fundoscopy skills. It can also be used to teach indirect ophthalmoscopy, direct ophthalmoscopy and the use of fundus lenses with a slit lamp. This inexpensive model features changeable pupil sizes (fully dilated 8mm and completely un-dilated 2mm), a calibrating scale that allows the eye to be set from -8 diopters to + 8 diopters to hypermetropia), easily inserted images of multiple retinal pathologies and durable plastic construction. тмте



Many patients find current explanations about corneas difficult to follow and may be apprehensive about procedures or conditions they don't understand. This unique model of the human eye allow the doctor or technician to quickly and completely familiarize the patient with their condition. Made of hard, durable plastic and measures 3.75 inches in

CORNEA EYE MODEL

AMCOREM

EXAM ROOM NEEDS

NEURO REHABILITATION

LOW VISION

EYE MODELS & CHARTS

VITREOUS FLOATER EYE MODEL



This Vitreous Model shows how floaters are caused and demonstrates the Weiss Ring. How shadows form on the retina can be demonstrated with a penlight. This model easily explains a difficult concept that happens to all patients and worries many. The model is 3 inches in diameter and is mounted on a white opaque plastic base.

RIVFEM

SYSTEMIC EFFECTS OF DIABETES MODEL



LAMINATED CHARTS

DISEASES OF THE EYE

AC1231

This innovative model demonstrates the effect of unsuccessfully treated diabetes on four areas: kidneys, nerves, arteries, and of course, retinas. One of the purposes of the model is the hope that by demonstrating the potential catastrophic effects of untreated diabetes on many body parts and organs patients will comply better with their prescribed medications, treatments, and diets - helping to save their eyesight.

AMSEDM

Clearly labeled eye charts help you explain problems in eye health. Make your office more impressive to your patients. Full color.

 $\rangle\rangle$



SIEPSER VITREOUS HEMORRHAGE DEMONSTRATOR EYE MODEL

This realistic model demonstrates the effects of vitreous hemorrhage and head position to patients. Contains simulated, suspended "blood particles" that float in the vitreous humor. The globe is glass and the base is plastic. Measures 2.25 inches in diameter.



DRY EYE MODEL

Demonstrate fluorescein staining on the cornea, lissamine green staining on the sclera, meibomian gland dysfunction, lacrimal gland dysfunction, and the lacrimal punctum. Excellent way to explain Dry Eye Disease to patients!

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

This laminated poster shows crossed eyes, chalazion, sty, basal cell, blepharitis, refractive error, blow out fracture, retinoblastoma, hypertensive retinopathy, SMD, optic neuritis, dry eye, lacrimal sac constriction, lacrimal gland inflammation, PVT, vitreous heme, conjunctivitis, foreign body, cataract, glaucoma, melanoma, ulcer, Herpes and other problems of vision.




GLAUCOMA POSTER

This popular poster readily explains the anatomy of the different glaucomas. Great for doctor/patient discussions. Defines Glaucoma and shows chronic, narrow angle, congenital and secondary. Poster is 20" x 26" flexible plastic lamination with metallic corners on the top for easy hanging. You can mark on it with write on/wipe off markers.

Size: 20" x 26"

ses of the Ey


AC1010




Guiden Ophthalmics' Dry Eye Mode

The Eye

AMDEM








EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY



Assessing Reading: Multiple Measures

This collection of tests assesses the skills related to reading. Written for the teacher or parent. Identifies the student with problems and determines the next step in instruction necessary to remediate that difficulty. Easy-to-follow instructions. Reproducible recording forms. Spanish assessments are included for the early grades.

ATP8120


The Optometry Team:

DEVELOPMENT MANUAL FOR SUPPORT STAFF IN OPTOMETRY PRACTICES - 4TH EDITION

BY: TONY HANKS, OD

"The Optometry Team" is the solution for training support staff. Contains over 400 easy-to-understand pages written by optometrist Tony Hanks. Staff will learn about eyecare, lenses, frames, contact lenses, ophthalmics, office procedures, etc. Includes hundreds of photos & illustrations to make learning easier.

A great staff training tool for any practice!



AUMAN

THE LOW VISION HANDBOOK FOR EYECARE PROFESSIONALS 2ND EDITION

BY: BARBARA BROWN, CO, MED



Clearly organized and simply presented, The Low Vision Handbook for Eyecare Professionals, Second Edition offers an introduction to all aspects of low vision, including a short history of low vision and the basic optics of magnifiers. Provides practical material on assessing low vision patients, the psychology of visual loss, and ways to alleviate patients' common fears. Additionally, the book explores current optical, non-optical, and electronic devices and their appropriate uses for various patients.

LVHBK

Visual Development, Diagnosis, and Treatment of the Pediatric Patient 2nd Edition



Robert H Duckman, Pam Schnell, Marc B Taub

Combines basic concepts of vision development with clinical diagnosis and treatment of vision disorders in infants, toddlers, children, and adolescents. Heavily updated, with new sections on timely issues and topics, the book is ideal for anyone who needs to know the practical aspects of evaluation and care of pediatric patients.

Paperback, 450 pages

VDDTPP

HIGH PERFORMANCE VISION


How To Improve Your Visual Acuity, Hone Your Motor Skills & Up Your Game

BY DONALD S. TEIG OD

HPVBK

During his work with professional athletes over the past forty years, Dr. Teig developed a series of visual and visualmotor performance tests to determine the strengths and weaknesses of their eyesight. After establishing an initial baseline of test results, the athletes were given specific exercises designed to improve their visual skills. He then tested them again and measured the results against the baseline. With each succession of exercises, their sight and motor coordination improved, as did their performance on the playing field. In High Performance Vision, Dr. Teig details his unique approach and offers his highly effective exercise regimen for improving your own vision.


Interested in adding Sports Vision to your practice? Visit our website to view the Teig Visual Calisthenics Sports Vision Training Kits!



SEE TO PLAY: THE EYES OF ELITE ATHLETES

BY MICHAEL A. PETERS OD





In See to Play, eminent optometrist Michael Peters, who works with many professional athletes and sports teams, addresses every aspect of this vital component of elite athletics. Included are in-depth discussions of all facets of vision accompanied by myriad exercises to help athletes hone and condition peripheral vision, acuity, speed of focus, eye- hand-body coordination, and more.



NEURO REHABILITATION

Books

EXAM ROOM NEEDS

ASSESSMENTS

OPTICAL

VISION THERAPY

NEURO REHABILITATION

LOW VISION

REFERENCE




VISUAL SECRETS FOR SCHOOL SUCCESS BY DR. BRENDA MONTECALVO

Handwriting, spelling, composition, math and reading can be



vision therapy. Once they are developed students often need to be coached on how to use the new visual skills when learning. Visual Secrets for School Success contains proven methods that help students read faster with good comprehension, improve handwriting and composition skills, spell better, and complete difficult math problems in less time. VSFSSHB



Do You See What I See? BY: DR. SANDY JOHAL, OD, FCOVD, FAAO

> them! Included in the book are testimonials by parents and patients praising the impact of vision therapy on their lives. This book is a tremendous resource for any optometry office to keep in their waiting room, library, or donate to community

> > Hardback: 32 Pages


When Your Child STRUGGLES THE MYTHS OF 20/20 VISION (WHAT EVERY PARENT NEEDS TO KNOW) When Written by Your David L. Cook. **O.D.** to explain Child the role of vision Struggles in learning to The Mythy of 20/20 Visio parents. Dr. David Cool 173 pages.

B200

Fax Orders to Save Time & Money at (574) 259-2102 or 2103 or Go To www.bernell.com 149

EXAM ROOM NEEDS

ASSESSMENTS



Applied Concepts in Vision Therapy 2.0

WRITTEN BY: LEONARD J. PRESS, OD, MARC B. TAUB, OD, MS, AND PAMELA H. SCHNELL, OD

This book is a valuable resource that will enhance your practice!

Twenty years after the ground-breaking original comes the 2.0 version that will serve as a valuable resource that will enhance your practice!

Covers a comprehensive array of topics on vision therapy, all in one place. Written and edited by over 50 experts with extensive experience incorporating vision therapy concepts successfully into their practices. Incorporated are new chapters on Optometric Phototherapy, Primitive Reflexes, Special Needs, Multisensory Care, Collaborating with Other Professions, Technology, and Visual Disability-Related Issues.

This book offers a unique blend of behavioral optometric care and classical optometric theory. The clinical approach helps you learn how to apply theoretical and didactic concepts to practice.

PRESSVT4

GUIDE TO VISION THERAPY™

The complete problem solving guide to contemporary vision therapy. Book covers the full scope of current techniques and is an excellent reference/source manual for either beginner or experienced professionals or assistants. Features individual training techniques with take-home instruction sheets provided. New format is a 3 ring binder to make copies of instruction sheets for patients to take home.



BCVT (ENGLISH)



The Vision Therapist's Toolkit Book

BY THOMAS HEADLINE, COVT, IRENE WAHLMEIER, COVT AND VICKI BEDES, COVT

A comprehensive text and optometric vision therapist training manual to acquaint the vision therapist to the tools of the profession: patches, lenses, prisms, anaglyphic filters and polarization. Authored by three vision therapists, it contains chapters on Anatomy; Visual Hygiene; Activities and Procedures for Gross Motor Skills, Accommodative Skills, Binocular Vision Skills, and Visual Perception/Visualization Skills. Spiral Bound, 176 pages. Has also been translated into Spanish.

HDVTSS (SPANISH)

CLINICAL MANAGEMENT OF BINOCULAR VISION 5TH EDITION

HDVTS (ENGLISH)

B0300 Paperback: 760 pages

HETEROPHORIC, ACCOMMODATIVE, AND EYE MOVEMENT DISORDERS BY: MITCHEL SCHEIMAN, OD AND BRUCE WICK, OD

The first new edition of Clinical Management of Binocular Vision in five years has been updated—across all chapters—with new guidelines and protocols based on the latest research in the field. You'll learn how to perform current, clinically accurate assessments, diagnoses, and therapies for a wide range of conditions. Each disorder is tackled in detail, touching on common symptoms, instrumentation, available treatment options, and more.


Features:

- New for this edition: 49 videos demonstrating proper diagnostic techniques and therapeutic procedures.
- · Now with new information on concussion-related vision disorders.
- · Enhanced guidelines for management of developing refractive error, including detailed description of preventative treatments for myopia.
- Formatted for easy reading, with succinct and bulleted text, numerous figures and illustrations, and step-by-step instructions.
- True-to-life case studies help you apply what you've read in a real-world clinical setting.
- Ideal for faculty, students who are new to the field, and veteran optometrists who need a refresh or easily accessible quick-reference.

INTRODUCTION TO BEHAVIORAL OPTOMETRY BY OEP

This series of paperback books are helpful additions to the professional library. Now offered individually by Bernell for non-members of OEP.

(INTRO TO CLINICAL OPTOMETRY	ІТСО)(Optometric Vision Therapy	OPVT	1	no litera Marian	
(THE FULL SCOPE OF RETINOSCOPY	TFSOR)(Working with the Brain-Inw red	WWTBI		and and	3
(STRABISMUS AND AMBLYOPIA	SAA)(VISION SCREENINGS	VISCR			
(VT: THE JOY OF OPTOMETRY	VTJOY)(PEDIATRIC OPTOMETRY: A PEDIATRICIAN'S PERSPECTIVE	POAPP		VISION	Ĩ
(DEVELOPMENTAL OPTOMETRY	DEVO)(PURSUITS & SACCADES: TESTING AND TRAINING	PSTT		THERAPS	
(Optometric Case Analysis	OCA)				-	




THE

	KLALKA	MALSA MALSA
ASION BELLAN	9399 	

Books

ASTRONAUT TRAINING: A SOUND ACTIVATED VESTIBULAR-VISUAL PROTOCOL FOR MOVING, LOOKING & LISTENING BY: MARY J. KAWAR, MS, OTR, & SHEILA M. FRICK, OTR



As a therapist, you will find this practical, spiral-bound handbook a valuable resource for expanding your treatment options. It presents clear, precise instructions and photographs of a vestibular-ocular protocol. It includes a companion CD that contains music and space sounds which support ocular control and provide energizing rhythms and inviting melodies to fully engage clients in the activities.

ATB


Developing Ocular Motor and Visual Perceptual Skills By: Development Skills Skills Development Skills Skills Development Skills Skills Development Skills Skills Development Skills

DEVELOPING OCULAR MOTOR & VISUAL PERCEPTUAL SKILLS BY: KENNETH A. LANE OD, FCOVD

Developing Ocular Motor and Visual Perceptual Skills contains daily lesson plans and practical tips on how to successfully start an activities program. Other helpful features include a glossary of terms and a reference list of individuals and organizations that work with learning disabled children to develop these skills.

A CLINICIAN'S GUIDE TO BALANCE AND DIZZINESS: EVALUATION AND TREATMENT

This book is an easy-to-use reference perfect for professionals who assess and treat balance impairments and dizziness. While it is an instructional text for physical therapy students and clinicians, it is also a great reference for established physicians, vestibular and balance therapy specialists, occupational therapists, nurse practitioners, physician assistants, audiologists, and athletic trainers.



CGBDET

DOMV

NEURO-OPTOMETRIC REHABILITATION

BY: WILLIAM PADULA, O.D.




Explore Neuro-Visual Processing Rehabilitation to learn new ways to think about vision. This book delves into the visual processing relationships of child development, motor and sensory interactions, and postural organization, led by vision, as the basis for understanding vision. The implications of understanding neuro-visual processing will change your thinking about vision as well as provide possibilities for helping millions of children and adults who have a neurological condition.

200 pages

NVPR

UNDERSTANDING & MANAGING VISION DEFICITS: A Guide for Occupational Therapists Third Edition



This book is a unique collaboration from occupational therapists, optometrists, and low vision rehabilitation specialists. The Third Edition of Understanding and Managing Vision Deficits is the go-to resource that will enable occupational therapists to develop a comprehensive understanding of vision, appreciate the various effects vision problems can have on the practice of occupational therapy, and to more effectively manage patients with vision disorders.


Hardback: 416 pages

UMVDOT

VISION REHABILITATION: MULTIDISCIPLINARY CARE OF THE PATIENT FOLLOWING BRAIN INJURY


BY: PENELOPE S. SUTER, OD & LISA H. HARVEY, OD

Providing the information required to understand, advocate for, and supply post-acute vision rehabilitative care following brain injury, Vision Rehabilitation: Multidisciplinary Care of the Patient Following Brain Injury bridges the gap between theory and practice. It presents clinical information and scientific literature supporting the diagnostic and therapeutic strategies applied in a comprehensive overview of current diagnostic and treatment strategies in adult post-brain injury vision rehabilitation.



Phantoms in the Brain

BY: V.S. RAMACHANDRAN, M.D., PH.D., & SANDRA BLAKESLEE



If you ever see patients with macular degeneration, or who have had a stroke, you need to know what is in this book! The author uses unusual clinical cases to explain how our vision really works. He explains blind sight, visual neglect, and other syndromes (as well as phantom limbs). This book covers it well. Dr. Ramachdran has lectured to N.O.R.A. about brain injuries & visual rehabilitation. He is considered one of the top two Neuroscientists of our time. "The Sherlock Holmes of Neuroscience". The other, Oliver Sacks MD, wrote the forward and says this book is "One of the most original and accessible neurology books of our generation." Well written and often guoted!

328 pages.

"The best book I've ever read about visual perception"- Craig Andrews, OD

OPTICAL

VISION THERAPY

TERMS AND CONDITIONS OF SALES INSIDE THE U.S.

TERMS: Standard payment terms are Net 30, subject to credit approval. Invoices may be subject to a 2% finance charge per month, if not paid within terms. Invoice payments are to be made via check or ACH. For any accounts approved for credit terms, invoices going past 60 days of due date will result in the account losing credit terms, unless written approval is granted by Bernell Management. All invoices going past 90 days will be turned over to our partner collections agency.

CREDIT CARDS: Bernell accepts MasterCard, Visa, Discover, and American Express.

PRICING: Current catalog prices supercede previous prices. Prices and item availability are subject to change without notice. Orders received are subject to approval. Bernell reserves the right to refuse orders and correct pricing and/or content errors in any of our sales and marketing materials.

ORDER PLACEMENT: Orders can be placed online at www.bernell.com, by phone at (574)259-2070, or by fax at (574)259-2102 or (574)259-2103. If placing a fax order, please include your name, company name, address, phone, fax number, Bernell account number, and country, if foreign. Web orders require registration to gain full access to professional pricing.

TELECOMMUNICATIONS SERVICES: Our phones are answered 8:00am - 5:30pm EST. After hours, the calls are recorded and orders are entered the next day. Internet orders received before 2pm EST will be processed that day, otherwise orders will process the next business day. Our customer service representatives may be reached at toll-free at (800) 348-2225, or +1 (574) 259-2070 during normal business hours. Live agents are also available for online chat during regular business hours.

SHIPPING: Unless specified when placing an order, at our discretion we will ship US Postal, UPS or FedEx. Shipments to known residential address will be charged a \$3.00 residential delivery fee, if requested to ship FedEx or UPS. For shipments with final destinations to Alaska, Hawaii, Puerto Rico, and the U.S. Territories, FedEx or UPS air shipping rates will apply and will be sent with a Pro Forma for your approval. Shipments that total 50lbs or more are subject to additional freight charges. The chart (shown on right) displays the ground shipping and handling charges that apply in the lower 48 states of the U.S. only.

THIRD-PARTY SHIPMENTS: Bernell will accommodate requests for the use of a third-party shipping company to pick up orders at our location. However, a \$2.50 packing and handling fee will be charged to the order. For international customers, there will be a \$5.00 per box handling fee assessed when using a third-party shipping company to ship outside the USA.

Shipping & Handling Charges

SPECIAL HANDLING: Bernell makes every effort to ship non-custom, in-stock items the same day if ordered by 2pm EST. However, orders requiring express shipment will be quoted carrier pricing, plus a priority handling fee of \$20 for express fulfillment, as these go outside our normal fulfillment process. All orders must be received and confirmed by 1pm EST in order to meet required time frames needed for the express fulfillment and carrier pick-ups. Refrigerated pharmaceuticals will require express shipping, but will not be charged the priority handling fee. These order will ship separately from your other items.

WARRANTY: All products are subject to manufacturer's instructions and warranties. Bernell does not provide any additional warranties.

RETURNS/DAMAGE or Loss Merchandise: It is the customers responsibility to examine merchandise carefully upon receipt. Claims for shortage or damaged goods must be made within 10 days of the invoice date. Merchandise can only be returned during the first 30 days from the invoice date and only after receiving a return merchandise authorization (RMA) number. This number may be obtained by calling our Customer Service Department at (574) 259-2070. All returns must come back in their original product packaging. Special Note: if claiming a damaged goods, then you must keep goods in the original shipping materials until items are inspected by the carrier. Bernell is not responsible for returns that we do not receive, or that do not include proper RMA identification. It is strongly recommended that you use a carrier which traces your shipment or that requires a signature. Products from the low vision pages may be returned with a RMA during the first 30 days, but there is a 15% restocking fee to do so. This is an unfortunate necessity because many offices loaned different magnifiers to their patients and then began returning any that the patient did not like. It is our belief that the low vision doctor should maintain a stock of demonstration magnifiers and we offer a Low Vision Demonstration Starter Kit (LVDEMO) at a discount for this purpose. We will only allow the return of unused products, unless a product is found to be defective upon first use. These must be reported in the first 30 days from the invoice as is required of all returns. Special ordered products are marked and are not returnable. They are items we ordered special for you and do not normally keep in inventory. All software is non-returnable with no exceptions. Some software items are offered in demos and we recommend this path if unsure about the software. Vision Therapy starter kits and primary care screening kits are returnable, but will be charged a 15% restocking fee.

Bernell wants to be a convenient and economical source for most of your office optical supplies. Bernell supports Vision Therapy research with gifts and product donations to several foundations, universities and through support of optometric groups such as COVD, OEPF and NORA, as well as our investments research & development at Bernell. Bernell welcomes your feedback regarding all areas of our business and how we can better serve your practice. If you want to send us your feedback, please use email us at info@bernell.com. Please provide only specific feedback on Bernell's existing products or marketing strategies. Bernell does support and assist with the development of new product ideas and bringing new products to market. In these instances, we are open to signing a non-disclosure agreement before listening to any product ideas. However, if a non-disclosure agreement is not in place, then any idea submissions will fall under the "Unsolicited Ideas Policy" of Bernell. To learn more about this policy, please visit our website at www.bernell.com and click on "Terms and Conditions" found under "Company Information." Bernell does not assume liability for practitioner or user compliance with state law or product instructions for use.

EMAIL CONTACTS

Customer Service: info@bernell.com	VP of Business Development: Joe Villari joev@bernell.com	Creative Director: Erik Warning	
President: Chris Andrews	VP of Operations: Al Martin	erik@bernell.com	

chris@bernell.com

al@bernell.com

FOR CUSTOMERS OUTSIDE OF THE USA YOU HAVE THE FOLLOWING OPTIONS:

- 1) You may contact any of our distributors listed on our website at www.bernell.com at the bottom of the home page.
- 2) You may order directly off our website. We will send an email quote including product costs, shipping and export handling charges and will require verification from you before we ship or charge your order.
- 3) You may fax to (574) 259-2102 or 2103
- 4) You may send an E-mail to info@bernell.com

Once your order has been received, we will provide you with a Pro Forma invoice. The proforma invoice will include the product costs, Export Handling Fees and actual shipping charges. All international orders require advanced payment.

WE ACCEPT THE FOLLOWING PAYMENT METHODS (US DOLLAR \$):

- MasterCard, VISA, Discover, or American Express credit cards card number and expiration date must be provided, as well as the card verification code (located on the back of MasterCard, Visa, or Discover Cards and on the front of American Express cards).
- Bank draft in the U.S. dollars drawn upon a U.S. bank.
- Wire transfer (\$25 fee will be assessed) or letter of credit (customer responsible for charges) are handled through: Vision Training Products, Inc., Mishawaka, IN, U.S.A. account.

Shipment will be made upon receipt of the advance payment. Orders are shipped: US Postal, UPS, or FedEx. You may also provide arrange pick-up via your own carrier or freight forwarder. No C.O.D.

GET SOCIAL WITH US:

Follow us on social media to join the **#conversationswithbernell**. We want to hear from you, the experts. Send and keep up with **#bernelltuesdays** tips! Join the **#bernellcommunity** and keep up with our special promotions, contests and education series! We have so much in store for you during 2022! Don't miss out! A DIREON OF VIEW THANKIN PRODUCTS. INC. 4016 N. HOME STREET MISHAWAKA, IN 46545

- f @bernellcorp
- (O) @bernellcorp
 - in @company/bernellcorp

@bernellcorp

Øbernellcorp